CareFusion

MicroLab

Operating Manual – English

Operating Manual – CareFusion USA

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Bruksanvisning – Svenska

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Contents

I.	Introduction – English.....	7
	Overview.....	9
	Getting started	10
	Calibration Check	15
	Customisation.....	17
	Administration Mode.....	18
	Paper Loading.....	19
	External printing	20
	PC connection using SPCS	20
	Cleaning the Transducer.....	22
	Servicing.....	23
	Consumables/ supporting products.....	24
	Specification of the MicroLab	25
II.	Introduction – CareFusion USA	28
	Overview.....	30
	Getting started	31
	Calibration Check	36
	Customisation.....	38
	Administration Mode.....	39
	Paper Loading.....	40
	External printing	41
	PC connection using SPCS	42
	Cleaning the Transducer.....	43
	Servicing.....	44
	Consumables/ supporting products.....	45
	Specification of the MicroLab	46

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Operating manual

I. Introduction - English

The MicroLab is mains/battery operated portable spirometer with the unique combination of ease of use and sophistication. Ease of use is assured through the use of context sensitive help screens, accessed at a touch of a button, that explain every MicroLab feature.

The MicroLab uses the CareFusion Digital Volume Transducer, an extremely stable form of volume transducer, which measures expired air directly at B.T.P.S (Body Temperature and Pressure with Saturated water vapour) thus avoiding the inaccuracies of temperature corrections. The transducer is insensitive to the effects of condensation and temperature and avoids the need for individual calibration prior to performing a test.

Test results may be uploaded to a PC using the optional Spirometry PC Software (SPCS) and patient details may be downloaded to the MicroLab.

Stored data may be printed to the integral thermal printer or to an external HP printer using the USB cable supplied.

Customise - System		
Region		
United Kingdom		▼
United Kingdom		
France		
Deutschland		
Nederland		
Espana		
USA		
ddmmyyyy		▼
Date Separator		
00/00/00		▼
Colour Printing		
Yes		▼
Cancel	Next	Finish
	Page 1 of 2	 15:30

When the unit is turned on for the first time the system customise menu is displayed. Touch the name of the required region with the stylus stored on the left hand side of the unit and then touch 'Finish'. This will set the regional defaults including language, if available.

If this screen is not displayed then touch the 'Customise' icon on the main menu followed by the 'System' icon.

Overview


The MicroLab uses a touch screen with icons representing each function available. A stylus, housed in the left hand side of the unit, is provided for icon screen activation and an optional mouse is available.

Touch the displayed time to adjust time and date.

Unused icons may be disabled by touching the blue background and selecting from the list displayed.

Touch the toolbox icon to adjust volume and brightness.

Four levels of battery charge are indicated by the segmented battery icon.

When this icon turns red the battery is nearly exhausted and the batteries must be charged – see Charging Procedure. The complete functionality is described on the help screen.

This is obtained by pressing the help button (?).

Help text exists for every screen viewed during the operation of the MicroLab.

You are recommended to make full use of the extensive Help screens provided.

Getting Started

When performing a spirometry test the recommended workflow is to enter the patient's details, or retrieve them from memory, perform the required test and then print and save the results.

Please ensure that the turbine transducer is plugged in to either of the first two sockets on the right hand side of the instrument.

Patient Select

Search	Database Usage
ID:	1%
Name:	

ID	Name
123	Smith, Adam
234HI889	Lawson, Christopher
54tjn666889	Denton, James

Cancel
Add
Today

10:51

Touch the 'Patients' icon to enter the patient database. The required patient may be selected from the stored patient list. If the patient details have not been previously stored then touch Add to enter the new patient's details. The patient details may also be downloaded from the optional SPCS software.

Once selected, the patient's name will appear at the bottom of the screen.

Use the help button to obtain further information.

New Patient

Patient Details

ID	<input style="width: 80%;" type="text"/>
Last Name	<input style="width: 80%;" type="text"/>
First Name	<input style="width: 80%;" type="text"/>
Sex	Male ▼
Origin	Caucasian ▼
Height (cm)	

1	2	3	4	5	6	7	8	9	0
q	w	e	r	t	y	u	i	o	p
↑	a	s	d	f	g	h	j	k	l
	z	x	c	v	b	n	m		
âü	.	,				←	←		

Cancel
Finish


10:16

To add a patient to the database use the on screen keyboard to type a unique patient ID and then touch the enter key. You will then be prompted for Last Name, First Name, Sex, Ethnic Origin, Height, Weight, date of Birth and Factor. A factor can be applied when testing individuals of other ethnic origins who would not normally be tested against the countries set of predicted values eg ECCS predicted values are used but Non Caucasian individuals tested. The factor alters the predicted value set on volume indices by the percentage applied.


The following factors are recommended when using ECCS normal values:

Hong Kong Chinese	100%
Japanese American	89%
Polynesians	90%
North Indians and Pakistanis	90%
South Indians and those of African descent	87%

Ref: Lung Volumes and Forced Ventilatory Flows. P.H.Quanjer et al.
Eur Respir J, 1993, 6, Suppl. 16p5-40


Once all the patient details are added the patient is added to the database and the main menu is displayed with the patient name displayed at the bottom of the screen. From the main menu select the required test, by touching the icon with the stylus. If the displayed patient is not required for testing then touch the patients name and options to change or remove the current patient will become available


If Relaxed Spirometry is selected then a volume/time graph will be displayed. Note that the unit may be customised to perform a relaxed Vital Capacity with tidal breathing or from a single expiration or single inspiration.

When a manoeuvre has been obtained touch 'Results' to view the indices, 'Again' to repeat the manoeuvre, 'Reject' to delete the manoeuvre or 'Done' to end the test.

All the active indices are displayed for any of the manoeuvres selected together with an option to review the volume/ time curves. The active indices listed can be changed by using the customisation option.


Select 'Done' to proceed to the Spirometry Main Menu.

From this menu the results of the test may be viewed, saved, or printed and notes may be added.

It is also possible to proceed to a forced baseline spirometry test, or a post medication relaxed spirometry test.

Select 'Exit' when all the required functions have been used.


If forced spirometry is selected the default graph will be displayed. This may be changed by touching the arrows at the top of the screen. Flow/Volume, Volume/time or child incentive default displays may be selected using the customise option from the main menu.

When the spirometry manoeuvre has been completed options to repeat the test, reject the test, and view results will be available.

At the end of the test options to view results, save results, print results, and to add notes will be available from the spirometry main menu.


Calibration Check

Connect a 3 litre syringe to the transducer with the adapter supplied and empty by pushing the handle fully in. Select 'Calibration Check' from the main menu and then select 'check Calibration'.

Fill the syringe by pulling the handle at a constant rate until the end stop is reached and then immediately empty the syringe completely. Try to maintain a flow rate that keeps the trace within the grey bands on the display.

Press 'Reject' to retry the calibration check at the required flow rate.

Press 'Again' to repeat the calibration check at a low flow rate.

Press 'Again' to repeat the calibration check at a high flow rate.

When a calibration check at all three flow rates has been completed press 'Done' to view the calibration check report screen.

Calibration Check Report

Expiratory calibration: -0.7% (mid flow 4.0l/s)
Inspiratory calibration: 1.2% (mid flow 5.5l/s)

Expiratory calibration: -0.3% (mid flow 0.7l/s)
Inspiratory calibration: -0.5% (mid flow 1.0l/s)

Expiratory calibration: 1.5% (mid flow 10.8l/s)
Inspiratory calibration: -0.7% (mid flow 11.4l/s)

Calibration check status: Pass

Done

  16:17

The calibration error for expiration and inspiration at each flow rate are displayed. The calibration error should be less than 3%. If a greater error is shown, repeat the procedure ensuring that the syringe is emptied and filled in a smooth manner without jerking the handle. If an error greater than 3% is still shown, inspect the turbine transducer and clean if necessary.

Customisation

The 'Customise' option from the main menu may be used to configure many of the features of your MicroLab and are divided into system and spirometry options.

System options allow you to configure the following:

Language.

Height and weight units.

Date format.

Date separator.

Colour or monochrome printing (on external printer).

Personalised printout heading.

Spirometry options allow you to configure the following:

Relaxed spirometry mode (with or without tidal breathing).

Predicted value sets.

Predicted area or line display.

Display default.

Incentive display type.

Printed graphs.

Best test criteria

Interpretation and Lung Age indication.

Dyspnoea score and smoking status.


Daily calibration reminder.

Manual temperature adjustment.

Indices selection.

Note that when the language is selected, the height and weight units, date format, and date separator will be automatically changed. However this automatic selection may be overridden manually.

Administration Mode

Administration mode allows the administrator to restrict the availability of functions to the user by disabling icons on the main menu. For example, after the unit has been configured to the administrator's requirements, disabling of the 'Customise' icon will prevent any further adjustment by the user. Similarly, disabling of the 'Database Management' icon will prevent the user from deleting any patient details or test results.


To enter administration mode turn the unit on whilst holding down the help key. The default access code is 0000. Type this number in using the on-screen keyboard. A number of functions are now available.

Touch the 'Change Access Code' icon in order to enter your personal access code and thereby restrict the use of administration to authorised personnel.


Press the help button to obtain a full description of the functions.

Paper Loading

To load a new roll of thermal paper lift the paper cover using the side levers, place the paper into the compartment as shown and close the cover firmly. It is recommended that only CareFusion thermal printer paper (Cat No. 36-PSA1600) is used with the MicroLab to avoid damage to the thermal printer head.

To tear off the paper pull the paper towards you and to the right as shown below:

External Printing

When printing, the MicroLab will default to an external USB printer, if connected.

Using the USB cable provided connect the mini A/B socket on the right hand side of the MicroLab to the input socket on the printer. For a list of compatible printers refer to the web site www.carefusion.com/micromedical or contact your local dealer.

It is recommended that during printing the batteries are on charge with the unit connected to the mains adapter.

NOTE: Keep the printer out of reach of the patient at all times.

Charging Procedure

The MicroLab's internal batteries are discharged when shipped from the factory and should be fully charged on first use. Plug the AC adapter into the mains supply and plug the adapter output plug into the power input socket on the right hand side of the instrument. The orange charging light next to the power input socket will flash to indicate charging and will turn on constantly to indicate full charge.

The batteries will take approximately 4 hours to become fully charged.

Note: Use only the AC adapter supplied. Use of any other type may cause permanent damage to the MicroLab and cause a fire or electrical hazard. Do not plug in and remove the power lead from the AC adapter repeatedly.

PC connection using SPCS

SPCS is an easy to use PC based windows application that interfaces to the MicroLab via a USB or serial port. It incorporates a database

into which patient details can be entered and downloaded to the MicroLab or test results may be uploaded from the MicroLab to the PC.

Using SPCS and the MicroLab, live blows can be performed with the PC directly controlling the operation of the MicroLab.

The results and graphs produced are displayed directly on the PC screen.

The spirometer is connected from the serial or USB port on the PC, to the port on the right hand side of the instrument using the USB or serial cable provided with SPCS.

Note: The MicroLab should only be connected to a computer that is manufactured in accordance with EN60950 1992/1993 – ‘Safety of Information Technology Equipment including Electrical Business Equipment’.

Keep the PC out of reach of the patient at all times.

It is recommended that whilst the unit is connected to a computer the mains adapter is used.

Looking after your Spirometer


Please observe the following precautions:

- Do not touch the screen with fingers. Use only the stylus provided.
- Use only a damp, lint free, cloth to clean the screen.
- Do not keep the spirometer in a damp place or expose it to extremes of temperature.
- Do not direct the transducer holder towards a strong light source whilst operating the spirometer.
- Check the AC charger for compatibility with local power rating.

Cleaning the Transducer

The transducer requires no routine maintenance or servicing. However, if you wish to disinfect or clean the transducer it may be removed by means of the following procedure:

1. Rotating the turbine transducer anti-clockwise until the locating pip lines up with the small rectangular cut-out in the housing as shown below.
2. Gently pull the transducer away from the housing.
3. The transducer may now be immersed in warm soapy water for routine cleaning or immersed in cold disinfecting solutions e.g. Perasafe for a maximum of 10 minutes (Alcohol and chloride solutions should be avoided). After cleaning/disinfecting, the transducer should be rinsed briefly in distilled water and dried. Perasafe is available from CareFusion in convenient 81g containers, Cat No. 36-SSC5000A.
4. Re-assemble the mouthpiece holder.

Servicing

There is no routine maintenance required for the MicroLab and there are no user serviceable parts in this instrument. Please return the unit to CareFusion or an authorised agent if servicing is required.

Troubleshooting

In the unlikely event that the display 'freezes' and the unit does not respond to any key presses, hold the on/off button down for 10 seconds until the unit switches off and then turn on again.

Symbols

Type B device

In accordance with Directive 93/42/EEC

0086

Disposal in compliance with WEEE

Environment

This instrument complies with directive EN60601-1-2 electromagnetic compatibility but can be affected by cellular phones and by electromagnetic interference exceeding levels specified in EN 50082-1:1992

This instrument is not for use with flammable anaesthetic gases or in an oxygen enriched atmosphere

Electrical classification

Class I equipment.

Consumables/supporting products

Cat. No.	Description
36-PSA1600	Thermal Printer Paper (pack of 5 rolls)
36-PSA1000	Adult Disposable Mouthpieces (500 per box)
36-SST1000	One-way Safety Mouthpieces (500 per box)
36-SST1250	One-way Safety Mouthpieces (250 per box)
36-PSA1200	Paediatric Disposable Mouthpieces (250 per box)
36-PSA1100	Paediatric Adaptor
36-MGF1025	MicroGard Pulmonary Filters (25 per box)
36-MGF1100	MicroGard Pulmonary Filters (100 per box)
36-SSC5000A	PeraSafe Sterilising Powder 81g (to make up 5 litres of solution)
36-VOL2104	Nose Clips (pack of 5)
36-SM2125	3 Litre Calibration Syringe
36-SPC1000	Spirometry PC Software

Caution

- Pulmonary filters are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may increase air resistance and lead to an incorrect measurement.”
- Mouthpieces are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may degrade materials and lead to an incorrect measurement.”

Training

CareFusion run spirometry training courses throughout the UK, sessions include:

- Anatomy and Physiology of Respiratory System
- Lung Volume Terminology
- Objectives of Spirometry
- Setting up the test & performing a good manoeuvre
- Errors of technique
- Interpretation of results
- Using Spirometry as an aid to diagnosing airways disease
- Choosing a Spirometer
- Case Studies

Specification of the MicroLab

General:

Storage:	>2000 tests including Flow/Volume loops and Volume/Time curves
Printer Output:	Compatible Hewlett Packard USB printers (Specified models only).
Display:	Colour 1/4VGA LCD.
Power supply:	Input 100 to 240V, 50 to 60Hz. Output 12V 2.5A (Class 1)
Battery Pack:	Rechargeable NiMH 8.4V 1A-hours.
Dimensions:	350x255x120 mm. Transducer 50x60x90mm.
Weight:	630g.
Operating Temperature:	0 to +40 °C
Operating Humidity:	30% to 90% RH
Transport and Storage Temperature:	-20 to + 70 °C
Transport and Storage Humidity:	10% to 90% RH

Spirometry Measurements:

Relaxed Expiratory Vital Capacity (VC)

Forced Expired Volume in 0.75 seconds (FEV.75)

Forced Expired Volume in 1 second (FEV1)

Forced Expired Volume in 3 second (FEV3)

Forced Expired Volume in 6 seconds (FEV6)

Forced Vital Capacity (FVC)

Peak Expiratory Flow Rate (PEF)

FEV_{0.75} as a percentage of VC (FEV.75/VC)

FEV_{0.75} as a percentage of FVC (FEV.75/FVC)

FEV₁ as a percentage of VC (FEV1/VC)

FEV₁ as a percentage of FVC (FEV1/FVC)

FEV₃ as a percentage of VC (FEV3/VC)

FEV₃ as a percentage of FVC (FEV3/FVC)

FEV_{0.75} as a percentage of FEV₆ (FEV.75/FEV6)

FEV₁ as a percentage of FEV₆ (FEV1/FEV6)

Maximum Expired Flow at 75% of FVC remaining (MEF75)

Maximum Expired Flow at 50% of FVC remaining (MEF50)

Maximum Expired Flow at 25% of FVC remaining (MEF25)

Mean Mid-Expiratory Flow Rate (MMEF)

Forced expiratory flow at 50% of volume as a percentage of VC (FEF50/VC)

Forced expiratory flow at 50% of volume as a percentage of FVC (FEF50/FVC)

Maximal voluntary ventilation indicated (MVV (ind))

Forced inspired volume in 1 second (FIV1)

Forced inspiratory Vital Capacity (FIVC)

Peak Inspiratory Flow Rate (PIF)

FIV₁ as a percentage of FIVC (FIV1/FIVC)

Forced inspiratory flow at 25% of inhaled volume (FIF25)

Forced inspiratory flow at 50% of inhaled volume (FIF50)

Forced inspiratory flow at 75% of inhaled volume (FIF75)

Forced expiratory flow at 50% of volume as a percentage of FIF50 (FEF50/FIF50)

The time taken between 25% and 75% of the forced expired volume (MET2575)

Forced Expiratory Time (FET)

Tidal Volume (TV)

Expiratory reserve volume (ERV)

Inspiratory reserve volume (IRV)

Inspiratory capacity (IC)

Expiratory Relaxed vital capacity (EVC)

Inspiratory vital capacity (IVC)

Breathing frequency rate (FR)

Inspiratory time (Ti)

Expiratory time (Te)

Ti as a % of total breath time (Ti/Ttot)

Tidal volume as a % of Ti (TV/Ti)

Tests per subject:	5 VC manoeuvre 8 FVC manoeuvres
Predicted Values:	Various – depends upon national preference
Transducer:	CareFusion Bi-Directional Digital Volume.
Resolution:	10ml volume 0.03l/s flow
Accuracy:	+/-3%. To ATS recommendations – Standardisation of spirometry 1994 update for flows and volumes.

CareFusion

MicroLab

Operating manual

Federal (USA) law restricts this device to sale by or on the order of a physician or licensed practitioner.

CareFusion,
22745 Savi Ranch Parkway,
Yorba Linda
CA 92887-4668
USA

II. Introduction – CareFusion USA

The MicroLab is mains/battery operated portable spirometer with the unique combination of ease of use and sophistication.

Ease of use is assured through the use of context sensitive help screens, accessed at a touch of a button, that explain every MicroLab feature.

The MicroLab uses a Digital Volume Transducer, an extremely stable form of volume transducer, which measures expired air directly at B.T.P.S (Body Temperature and Pressure with Saturated water vapour) thus avoiding the inaccuracies of temperature corrections. The transducer is insensitive to the effects of condensation and temperature and avoids the need for individual calibration prior to performing a test.

Test results may be uploaded to a PC using the optional Spirometry PC Software (SPCS) and patient details may be downloaded to the MicroLab.

Stored data may be printed to the integral thermal printer or to an external USB printer.

The MicroLab utilizes a **single patient use** disposable mouthpiece that must be disposed of after use.

The MicroLab provides a suggested interpretation **that must be supported by clinical judgement.**

Overview


The MicroLab uses a touch screen with icons representing each function available. A stylus, housed in the left hand side of the unit, is provided for icon screen activation.

Touch the displayed time to adjust time and date.

Unused icons may be disabled by touching the blue background and selecting from the list displayed.

Touch the toolbox icon to adjust volume and brightness.

Four levels of battery charge are indicated by the segmented battery icon.

When this icon turns red the battery is nearly exhausted and the batteries must be charged – see Charging Procedure. The complete functionality is described on the help screen.

This is obtained by pressing the help button.

Help text exists for every screen viewed during the operation of the MicroLab.

You are recommended to make full use of the extensive Help screens provided.

Getting Started

When performing a spirometry test the recommended workflow is to enter the patient's details, or retrieve them from memory, perform the required test and then print and save the results.

Please ensure that the turbine transducer is plugged in to either of the first two sockets on the right hand side of the instrument.

Patient Select

Search	Database Usage
ID:	1%
Name:	

ID	Name
123	Smith, Adam
234HI889	Lawson, Christopher
54tjn666889	Denton, James

Cancel
Add
Today

10:51

Touch the 'Patients' icon to enter the patient database. The required patient is selected from the stored patient list. If the patient details have not been previously stored then touch Add to enter the new details. The patient details may also be downloaded from the optional SPCS software.

Once selected patient's name will appear at the bottom of the screen.

Use the help button to obtain further information.

New Patient

Patient Details

ID	<input style="width: 95%;" type="text"/>
Last Name	
First Name	
Sex	Male ▼
Origin	Caucasian ▼
Height (cm)	

1	2	3	4	5	6	7	8	9	0
q	w	e	r	t	y	u	i	o	p
↑	a	s	d	f	g	h	j	k	l
	z	x	c	v	b	n	m		
âü	.	,				←	←		


Cancel
Finish

10:16

To add a patient to the database use the on screen keyboard to type a unique patient ID and then touch the enter key. You will then be prompted for Last Name, First Name, Sex, Ethnic Origin, Height, Weight and date of Birth.

Once all the patient details are added the patient is added to the database and the main menu is displayed with the patient name displayed at the bottom of the screen. From the main menu select the required test, by touching the icon with the stylus. If the displayed patient is not required for testing then touch the patients name and options to change or remove the current patient will become available.

If Relaxed Spirometry is selected then a volume/time graph will be displayed. Note that the unit may be customised to perform a relaxed vital capacity with tidal breathing or from a single expiration or single inspiration.

When a satisfactory manoeuvre has been obtained touch 'Results' to view the indices and 'Done' to end the test.

Relaxed Result

Base

1*


Indice	Value	%Pred	[Min	Pred]
EVC	4.70				
IC	3.61				
TV	0.36				
ERV	1.09				
IRV	3.25				
FR	27.0				
Ti	1.16				
Te	0.97				
Ti/Ttot	54				

Graph Done

Christopher Lawson 13:28

All the active indices are displayed for any of the manoeuvres selected together with an option to review the volume/ time curves.


Select 'Done' to proceed to the Spirometry Main Menu.

From this menu the results of the test may be viewed, saved, or printed and notes may be added.

It is also possible to proceed to a forced baseline spirometry test, or a post medication relaxed spirometry test.

Select 'Exit' when all the required functions have been used.


If forced spirometry is selected the default graph will be displayed. This may be changed by touching the arrows at the top of the screen. Flow/volume, volume/time or child incentive default displays may be selected using the customise option from the main menu.

When the spirometry manoeuvre has been completed options to repeat the test, reject the test, and view results will be available.

At the end of the test options to view results, save results, print results, and to add notes will be available from the spirometry main menu.

Calibration Check


Connect a 3 litre syringe to the transducer with the adapter supplied and empty by pushing the handle fully in.

Select 'Calibration Check' from the main menu and then select 'check Calibration'.

Fill the syringe by pulling the handle at a constant rate until the end stop is reached and then immediately empty the syringe completely.


Press 'Done' to view the calibration check results.

Press 'Again' to repeat the calibration check at a low flow rate.

Press 'Again' to repeat the calibration check at a high flow rate.

When a calibration check at all three flow rates has been completed press 'Done' to view the calibration check report screen.

The calibration error for expiration and inspiration are displayed. The calibration error should be less than 3%. If a greater error is shown, repeat the procedure ensuring that the syringe is emptied and filled in a smooth manner without jerking the handle. If an error greater than 3% is still shown, inspect the turbine transducer and clean if necessary.

Customisation

The 'Customise' option from the main menu may be used to configure many of the features of your MicroLab and are divided into system and spirometry options.

System options allow you to configure the following:

Language.

Height and weight units.

Date format.

Date separator.

Colour or monochrome printing (on external printer).

Personalised printout heading.

Spirometry options allow you to configure the following:

Relaxed spirometry mode (with or without tidal breathing).

Predicted value sets.

Predicted area or line display.

Display default.

Incentive display type.

Printed graphs.

Best test criteria

Interpretation and Lung Age indication.

Dyspnoea score and smoking status.

Calibration reminder and manual temperature adjustment.

Indices selection.

Note that when the language is selected, the height and weight units, date format, and date separator will be automatically changed. However this automatic selection may be overridden manually.

Administration Mode


Administration mode allows the administrator to restrict the availability of functions to the user by disabling icons on the main menu. For example, after the unit has been configured to the administrator's requirements, disabling of the 'Customise' icon will prevent any further adjustment by the user. Similarly, disabling of the 'Database Management' icon will prevent the user from deleting any patient details or test results.


To enter administration mode turn the unit on whilst holding down the help key. The default access code is 0000. Type this number in using the on-screen keyboard. A number of functions are now available. Press the help button to obtain a full description of the functions.

Paper Loading

To load a new roll of thermal paper lift the paper cover using the side levers, place the paper into the compartment as shown and close the cover firmly.

To tear off the paper pull the paper towards you and to the right as shown below:

External Printing

When printing the MicroLab will default to an external USB printer if connected.

Using the USB cable provided connect the mini A/B socket on the right hand side of the MicroLab to the input socket on the printer. For a list of compatible printers contact your local dealer

It is recommended that during printing the batteries are on charge with the unit connected to the mains adapter.

NOTE: Keep the printer out of reach of the patient at all times.

Charging Procedure

The MicroLab's internal batteries are discharged when shipped from the factory and should be fully charged on first use. Plug the AC adapter into the mains supply and plug the adapter output plug into the power input socket on the right hand side of the instrument. The orange charging light next to the power input socket will flash to indicate charging and will turn on constantly to indicate full charge.

The batteries will take about 4 hours to become fully charged.

Note: Use only the AC adapter supplied. Use of any other type may cause permanent damage to the MicroLab and cause a fire or electrical hazard. Do not plug in and remove the power lead from the AC adapter repeatedly.

PC connection using SPCS

SPCS is an easy to use PC based windows application that interfaces to the MicroLab via a USB or serial port. It incorporates a database into which patient details can be entered and downloaded to the MicroLab or test results may be uploaded from the MicroLab to the PC.

Using SPCS and the MicroLab, live blows can be performed with the PC directly controlling the operation of the MicroLab.

The results and graphs produced are displayed directly on the PC screen.

The spirometer is connected from the serial or USB port on the PC, to the port on the right hand side of the instrument using the USB or serial cable provided with SPCS.

Note: the MicroLab should only be connected to a computer that is manufactured in accordance with EN60950 1992/1993 – ‘Safety of Information Technology Equipment including Electrical Business Equipment’.

Keep the PC out of reach of the patient at all times.

It is recommended that whilst the unit is connected to a computer the mains adapter is used.

Looking after your Spirometer


Please observe the following precautions:

- Do not touch the screen with fingers. Use only the stylus provided.
- Use only a damp, lint free, cloth to clean the screen.
- Do not keep the spirometer in a damp place or expose it to extremes of temperature.
- Do not direct the transducer holder towards a strong light source whilst operating the spirometer.
- Check the AC charger for compatibility with local power rating.

Cleaning the Transducer

The transducer requires no routine maintenance or servicing. However, if you wish to sterilise or clean the transducer it may be removed by means of the following procedure:

1. Rotating the turbine transducer anti-clockwise until the locating pip lines up with the small rectangular cut-out in the housing as shown below.
2. Gently pull the transducer away from the housing.
3. The transducer may now be immersed in warm soapy water for routine cleaning or immersed in cold sterilising solutions e.g. Perasafe for a maximum of 10 minutes (Alcohol and chloride solutions should be avoided). After cleaning/sterilising, the transducer should be rinsed briefly in distilled water and dried. Perasafe is available from CareFusion in convenient 81g containers.
4. Re-assemble the mouthpiece holder.

Servicing

There is no routine maintenance required for the MicroLab and there are no user serviceable parts in this instrument. Please return the unit to CareFusion or an authorised agent if servicing is required.

Troubleshooting

In the unlikely event that the display 'freezes' and the unit does not respond to any key presses, hold the on/off button down for 10 seconds until the unit switches off and then turn on again.

Symbols

Type B device

Environment

This instrument complies with directive EN60601-1-2 electromagnetic compatibility but can be affected by cellular phones and by electromagnetic interference exceeding levels specified in EN 50082-1:1992

This instrument is not for use with flammable anaesthetic gases or in an oxygen enriched atmosphere

Electrical classification

Class I equipment.

Consumables / Supporting Products

Description

Thermal Printer Paper (pack of 5 rolls)

Archive Printer Paper ~ 25 Year Anti-Fade (pack of 5 rolls)

Adult Disposable Mouthpieces (500 per box)

One-way Safety Mouthpieces (500 per box)

One-way Safety Mouthpieces (250 per box)

Paediatric Disposable Mouthpieces (250 per box)

Paediatric Adaptor

MicroGard Pulmonary filters

MicroGard Pulmonary filters

PeraSafe Sterilising Powder 81g (to make up 5 litres of solution)

Nose Clips (pack of 5)

3 Litre Calibration Syringe

SPCS Software

For further information or to place an order for Disposables/
Supporting products please contact CAREFUSION, your local
distributor or view our website www.carefusion.com/micromedical

For US Customer Care: Toll Free 1-800-231-2466, phone 1-714-
283-2228

Caution

- Pulmonary filters are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may increase air resistance and lead to an incorrect measurement."
- Mouthpieces are single patient use. If used on more than one patient there is a risk of cross-infection. Repeat use may degrade materials and lead to an incorrect measurement."

Specification of the MicroLab

General:

Storage :	>2000 tests including Flow/Volume loops and Volume/Time curves
Printer Output:	Compatible Hewlett Packard USB printers (Specified models only).
Display:	Colour 1/4VGA LCD.
Power supply:	Input 100 to 240V, 50 to 60Hz. Output 12V 2.5A (Class 1)
Battery Pack:	Rechargeable NiMH 8.4V 1A-hours.
Dimensions:	350x255x120 mm. Transducer 50x60x90mm.
Weight:	630g.
Operating Temperature:	0 to +40°C
Operating Humidity:	30% to 90% RH
Transport and Storage Temperature:	-20 to + 70°C
Transport and Storage Humidity:	10% to 90% RH

Spirometry Measurements:

Relaxed Expiratory Vital Capacity (VC)
Forced Expired Volume in 0.75 seconds (FEV.75)
Forced Expired Volume in 1 second (FEV1)
Forced Expired Volume in 3 second (FEV3)
Forced Expired Volume in 6 seconds (FEV6)
Forced Vital Capacity (FVC)
Peak Expiratory Flow Rate (PEF)
FEV_{0.75} as a percentage of VC (FEV.75/VC)
FEV_{0.75} as a percentage of FVC (FEV.75/FVC)
FEV₁ as a percentage of VC (FEV1/VC)
FEV₁ as a percentage of FVC (FEV1/FVC)

FEV₃ as a percentage of VC (FEV3/VC)
FEV₃ as a percentage of FVC (FEV3/FVC)
FEV_{0.75} as a percentage of FEV₆ (FEV.75/FEV6)
FEV₁ as a percentage of FEV₆ (FEV1/FEV6)
Maximum Expired Flow at 75% of FVC remaining (MEF75)
Maximum Expired Flow at 50% of FVC remaining (MEF50)
Maximum Expired Flow at 25% of FVC remaining (MEF25)
Mean Mid-Expiratory Flow Rate (MMEF)
Forced expiratory flow at 50% of volume as a percentage of VC (FEF50/VC)
Forced expiratory flow at 50% of volume as a percentage of FVC (FEF50/FVC)
Maximal voluntary ventilation indicated (MVV(ind))
Forced inspired volume in 1 second (FIV1)
Forced inspiratory Vital Capacity (FIVC)
Peak Inspiratory Flow Rate (PIF)
FIV₁ as a percentage of FIVC (FIV1/FIVC)
Forced inspiratory flow at 25% of inhaled volume (FIF25)
Forced inspiratory flow at 50% of inhaled volume (FIF50)
Forced inspiratory flow at 75% of inhaled volume (FIF75)
Forced expiratory flow at 50% of volume as a percentage of FIF50 (FEF50/FIF50)
The time taken between 25% and 75% of the forced expired volume (MET2575)
Forced Expiratory Time (FET)
Tidal Volume (TV)
Expiratory reserve volume (ERV)
Inspiratory reserve volume (IRV)
Inspiratory capacity (IC)

Tests per subject:	5 VC manoeuvre 8 FVC manoeuvres
Predicted Values:	Various – depends upon national preference
Transducer:	CareFusion Bi-Directional Digital Volume.
Resolution:	10ml volume 0.03l/s flow
Accuracy:	+/-3%. To ATS recommendations – Standardisation of spirometry 1994 update for flows and volumes.

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Relaxed Expiratory Vital Capacity (VC)
Forced Expired Volume in 0.75 seconds (FEV.75)
Forced Expired Volume in 1 second (FEV1)
Forced Expired Volume in 3 second (FEV3)
Forced Expired Volume in 6 seconds (FEV6)
Forced Vital Capacity (FVC)
Peak Expiratory Flow Rate (PEF)
FEV_{0.75} as a percentage of VC (FEV.75/VC)
FEV_{0.75} as a percentage of FVC (FEV.75/FVC)

FEV₁ as a percentage of VC (FEV1/VC)
FEV₁ as a percentage of FVC (FEV1/FVC)
FEV₃ as a percentage of VC (FEV3/VC)
FEV₃ as a percentage of FVC (FEV3/FVC)
FEV_{0.75} as a percentage of FEV₆ (FEV.75/FEV6)
FEV₁ as a percentage of FEV₆ (FEV1/FEV6)
Maximum Expired Flow at 75% of FVC remaining (MEF75)
Maximum Expired Flow at 50% of FVC remaining (MEF50)
Maximum Expired Flow at 25% of FVC remaining (MEF25)
Mean Mid-Expiratory Flow Rate (MMEF)
Forced expiratory flow at 50% of volume as a percentage of VC (FEF50/VC)
Forced expiratory flow at 50% of volume as a percentage of FVC (FEF50/FVC)
Maximal voluntary ventilation indicated (MVV(ind))
Forced inspired volume in 1 second (FIV1)
Forced inspiratory Vital Capacity (FIVC)
Peak Inspiratory Flow Rate (PIF)
FIV₁ as a percentage of FIVC (FIV1/FIVC)
Forced inspiratory flow at 25% of inhaled volume (FIF25)
Forced inspiratory flow at 50% of inhaled volume (FIF50)
Forced inspiratory flow at 75% of inhaled volume (FIF75)
Forced expiratory flow at 50% of volume as a percentage of FIF50 (FEF50/FIF50)
The time taken between 25% and 75% of the forced expired volume (MET2575)
Forced Expiratory Time (FET)
Tidal Volume (TV)
Expiratory reserve volume (ERV)
Inspiratory reserve volume (IRV)
Inspiratory capacity (IC)

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Customer contact information

UK Customers only

For all Sales Order processing for products, training and spare parts, Service and Technical Support enquiries, please contact the following:

CareFusion UK 232 Ltd
UK Customer Service & Support
The Crescent
Jays Close
Basingstoke
RG22 4BS

Customer Service Sales Enquiries:

Telephone: 01256 388550

Email: micro.uksales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 01256 388552

Email: micro.ukservice@carefusion.com

Technical Support Enquiries:

Telephone: 01256 388551

Email: support.rt.eu@carefusion.com

International customers only

For all Sales Order processing for products and Spare parts, Service and Technical Support enquiries, please contact the following:

Carefusion Germany 234 GmbH
Customer Service & Support International
Leibnizstrasse 7
97204 Hoechberg
Germany

Customer Service Sales Enquiries:

Telephone: 0049 931 4972 670

Email: micro.international-sales@carefusion.com

Factory Repair and Administration Enquiries:

Telephone: 0049 931 4972 867

Email: support.admin.eu@carefusion.com

Technical Support Enquiries:

Telephone: 0049 931 4972 127

Email: support.rt.eu@carefusion.com

CareFusion UK 232 Ltd.,
Quayside
Chatham Maritime
Kent ME4 4QY
U.K.

CE
0086

JDE no. 36-MAN1273
Drg no. 085-40WW
Issue 1.0
December 2009
© CareFusion 2009