




The Sanibel Supply Catalogue

Quality Supplies for your Clinic



Sanibel
Supply

Contents



Pages 5 – 7

Disposable ear tips



Pages 9 – 10

Probe tips



Pages 11 – 12

Supplies for insert earphones



Page 13

Sanibel™ Infant EarCup™



Page 14

REM probe tubes for HA fitting



Pages 15 – 17

Cerumen management accessories



Pages 18 – 19

Impression materials and supplies



Page 20

Hearing aid fitting



Pages 21–22

Hearing aid maintenance and care



Pages 23 – 25

Disposable electrodes



Page 27

Electrode cables



Pages 28 – 29

Patient preparation



Pages 30 – 31

Infection control



Pages 32 – 35

Parts for audiometry



Page 36

Balance supplies



Page 37

Paper



Page 38

Otoscopy supplies



Page 39

Batteries and various items

List of abbreviations



Auditory Brainstem Response



Audiometry



Balance



Hearing Instrument Fitting



Impedance



Otoacoustic Emissions



Single patient use only!



Product is not made with natural rubber latex



Why use ear tips for single use only?

Sanibel Supply® recommends using disposable ear tips for single use only to stop ear infections from spreading

At Sanibel Supply® we think that most audiologists would agree with us that it is crucial to avoid inadvertent transmission of disease within the audiology clinic. The good news is that it is quite easy to ensure that ear infections are not spread from one patient to another; use a new ear tip for each patient, and you will significantly reduce or eliminate the incidences of cross contamination.

Cerumen may be an infectious agent

Cerumen is not considered an infectious agent unless contaminated with blood or mucus, and due to its color and viscosity, visual detection of blood or ear discharge contaminants may be difficult. Therefore, you should treat cerumen as a potentially infectious agent. It is especially important to be mindful of using new tips when testing children and those with weaker immune systems, such as the frail and elderly as they are more susceptible to infections.

Are disposable ear tips expensive?

Sanibel Supply® offers volume discounts on ear tips. This will often offset the time and cost involved in the ear tip cleaning process.

Be more cost effective in your clinic

By using a new ear tip pr. every patient, you reach a seal much more reliable than with re-used ear tips. How often have you not been fumbling for minutes before you get a seal with a patient? Using the right ear tip and a fresh one pr. patient ensures you a reliable fit with every patient, meaning you can calculate less minutes for testing a patient and in most cases allow you to have an extra patient pr. day without having stressed out personnel. The cost of using one ear tip pr. patient is far less than what you earn on your next patient, plus you can promote your clinic as taking patient safety and hygiene seriously.

6 strong reasons why you should use a new ear tip after each patient:

1

You reduce the risk of cross contamination between your patients

Unfortunately, it is not possible to be 100% certain that complete decontamination has taken place when cleaning ear tips.

2

You minimize the risk of non-sealing ear tips

Ear tips degrade over time, which can make it more difficult to attain a good seal, hence leading to less accurate results.

3

Your patients avoid allergic reactions from cleaning agent residue

Residue from cleaning agents may cause allergic reactions in some patients.

4

Your patients will appreciate it

Patients feel safer and are more comfortable knowing that new ear tips are being used.

5

It's so much easier

Cleaning ear tips takes time, so why not give your staff the time to do what really matters?

















6

It's so much faster

A fresh ear tip creates a seal quicker and more reliable, so you do not waste time performing a standard procedure.



Selecting the correct ear tip

Age	Ear tips inserted against ear canal	Ear tips inserted against entrance of ear canal			Ear tips for Echo-Screen®	Ear tips for AccuScreen®
0.5 y/o		 Green, pediatric			 4.0 mm, ANA	 3.7 mm, AOT
					 4.5 mm, ANA	 4.0 mm, AOT
		 3-5 mm, flanged	 3-6 mm, green	 6 mm, green	 5.0 mm, ANA	 4.5 mm, AOT
	 7.5 mm, umbrella	 4-7 mm, flanged	 6 mm, yellow		 infant tree, ANA	 5.0 mm, AOT
	 9 mm, umbrella	 5-8 mm, flanged			 mushroom, ANA	 4-7 mm, AOT
	 10 mm, umbrella	 7 mm, mushroom	 7 mm, red	 7 mm, yellow		
	 11 mm, umbrella	 8 mm, mushroom	 8 mm, blue	 8 mm, red		
	 13 mm, umbrella	 9 mm, mushroom	 9 mm, green	 9 mm, blue		
		 10 mm, mushroom	 10 mm, yellow	 10 mm, green		
	 15 mm, umbrella	 11 mm, mushroom	 11 mm, red	 11 mm, yellow		
10 y/o		 12 mm, mushroom	 12 mm, blue	 12 mm, red		
	 19 mm, umbrella	 13 mm, mushroom	 13 mm, green	 13 mm, blue		
		 14 mm, mushroom	 14 mm, yellow	 14 mm, green		
		 15 mm, mushroom	 16 mm, red	 16 mm, yellow		
	 22 mm, umbrella	 19 mm, mushroom	 18.5 mm, blue			
Adults						

Selecting the correct size of ear tip takes practice

In general, flanged or mushroom-shaped ear tips will be used for most testing (Impedance, DPOAE, TEOAE and Automated ABR).

Sometimes, an umbrella-shaped ear tip will be used for performing a quick impedance screening.

The table can be used as a guide to help you select the appropriate ear tip size.

Disposable Sanibel™ ADI ear tips

For use with the following acoustic immittance instruments and otoacoustic emission systems with CLEAR probe tip: **Amplivox**: all Otowave versions. **Interacoustics**: Titan, AT235, AA220, Eclipse, MT10, Sera, OtoRead II. **MAICO**: easyTymp, touchTymp MI 22, MI 23, MI 24, MI 26, MI 34, MI 44, easyScreen, EroScan Plus,



ADI Assortment Box

Part no.	Description
8511445	Various sizes for tympanometry and OAE
8501572	Screening ear tip kit, various sizes for newborn screening



Flanged infant ear tips

Part no.	Description	Unit
8502837	Pediatric (green)	Bag/25 pcs
8502838	Pediatric (green)	Bag/100 pcs
8012963	3-5 mm (red)	Bag/25 pcs
8012964	3-5 mm (red)	Bag/100 pcs
8012965	4-7 mm (blue)	Bag/25 pcs
8012966	4-7 mm (blue)	Bag/100 pcs
8012967	5-8 mm (yellow)	Bag/25 pcs
8012968	5-8 mm (yellow)	Bag/100 pcs

VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags



Medical grade silicone rubber
Biocompatible



Umbrella ear tips

Part no.	Description	Unit
8012985	7.5 mm (green)	Bag/25 pcs
8012986	7.5 mm (green)	Bag/100 pcs
8012987	9 mm (red)	Bag/25 pcs
8012988	9 mm (red)	Bag/100 pcs
8012989	10 mm (blue)	Bag/25 pcs
8012990	10 mm (blue)	Bag/100 pcs
8012991	11 mm (yellow)	Bag/25 pcs
8012992	11 mm (yellow)	Bag/100 pcs
8012993	13 mm (green)	Bag/25 pcs
8012994	13 mm (green)	Bag/100 pcs
8012995	15 mm (red)	Bag/25 pcs
8012996	15 mm (red)	Bag/100 pcs
8012997	19 mm (blue)	Bag/25 pcs
8012998	19 mm (blue)	Bag/100 pcs
8012999	22 mm (yellow)	Bag/25 pcs
8013000	22 mm (yellow)	Bag/100 pcs



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags



Mushroom ear tips

Part no.	Description	Unit
8013001	7 mm (blue)	Bag/25 pcs
8013002	7 mm (blue)	Bag/100 pcs
8013003	8 mm (yellow)	Bag/25 pcs
8013004	8 mm (yellow)	Bag/100 pcs
8012969	9 mm (green)	Bag/25 pcs
8012970	9 mm (green)	Bag/100 pcs
8012971	10 mm (red)	Bag/25 pcs
8012972	10 mm (red)	Bag/100 pcs
8012973	11 mm (blue)	Bag/25 pcs
8012974	11 mm (blue)	Bag/100 pcs
8012975	12 mm (yellow)	Bag/25 pcs
8012976	12 mm (yellow)	Bag/100 pcs
8012977	13 mm (green)	Bag/25 pcs
8012978	13 mm (green)	Bag/100 pcs
8012979	14 mm (red)	Bag/25 pcs
8012980	14 mm (red)	Bag/100 pcs
8012981	15 mm (blue)	Bag/25 pcs
8012982	15 mm (blue)	Bag/100 pcs
8012983	19 mm (yellow)	Bag/25 pcs
8012984	19 mm (yellow)	Bag/100 pcs



Single use ear tips for immittance instruments & OAE screeners

For use with many acoustic immittance instruments and otoacoustic emission systems, including the **MAICO** EroScan, **Interacoustics** OtoRead and **Grason-Stadler** TympStar.



VOLUME DISCOUNT AVAILABLE for minimum orders of 10 bags

Infant ear tips

Part no.	Description	Unit
8013005	3 mm (yellow)	Bag/25 pcs
8013006	3 mm (yellow)	Bag/100 pcs
8002036	4 mm (red)	Bag/25 pcs
8002035	4 mm flanged (red)	Bag/100 pcs



IA series ear tips

For use with many acoustic immittance instruments and otoacoustic emission systems, including: **GSI**: TymStar Pro, Allegro, Novus and Corti. **Interacoustics** (clear probe): AT235, AA222 and Titan and Otoread II, MT10. **MAICO**: touchTym Racecar, touchTym MI24, MI26, MI34, MI36, easyTym. **Madsen/Otometrics**: Zodiac (new), Accuscreen. **Path Medical**: all Sentiero products.



Mushroom ear tips

Part no.	Description	Unit
8515491	Sample Kit (3-18.5mm), 3, 6-10 count each	
8520348	3-6 mm (green)	Bag/100 pcs
8520349	6 mm (yellow)	Bag/100 pcs
8520350	7 mm (red)	Bag/100 pcs
8520351	8 mm (blue)	Bag/100 pcs
8520352	9 mm (green)	Bag/100 pcs
8520353	10 mm (yellow)	Bag/100 pcs
8520354	11 mm (red)	Bag/100 pcs
8520355	12 mm (blue)	Bag/100 pcs
8520356	13 mm (green)	Bag/100 pcs
8520357	14 mm (yellow)	Bag/100 pcs
8520358	16 mm (red)	Bag/100 pcs
8520359	18.5 mm (blue)	Bag/100 pcs



IMPEDANCE EAR TIPS – fading product

For use with the following systems with BLACK probe tip:
Interacoustics: Old MT10, AT235, AA22x, MI 22, MI 23, MI 24, MI 26, MI 34, MI 44. **FADING OUT** product: consider AKR instead.



VOLUME DISCOUNT AVAILABLE
for minimum orders of 10 bags

Umbrella ear tips

Part no.	Description	Unit
8107617	10 mm (green)	Bag/100 pcs
8107618	13 mm (green)	Bag/100 pcs
8107619	16 mm (green)	Bag/100 pcs
8107620	19 mm (blue)	Bag/100 pcs
8107621	22 mm (blue)	Bag/100 pcs



Save time pr. patient, earn more money and help more patients!

Use a fresh ear tip pr. patient and be sure to reach a seal right away, when inserting the ear tip. A fresh ear tip contains the optimal conditions for reaching a seal when inserting it into the patients' ear, whereas a multi-used one often will cause you trouble and delay you in your process. By allowing yourself to plan less time pr. patient for testing, you can see more patients pr. day. That is worth much more than cleaning an ear tip, both for your clinic and your patients.

AKR mushroom ear tips

Use with the following acoustic immittance instruments and otoacoustic emission systems with BLACK probe tip from **Interacoustics**: AT235, AA220, AA220, Eclipse, MT10, Otoread (I), **MAICO**: touchTym RaceCar MI 22, MI 23, MI 24, MI 26, MI 34, MI 44. **Amplivox**: Otowave (all with white probe). **Madsen/Otometrics**: Capella 2, Zodiac 901, Biologic AuDX Pro (+Flex). **Neurosoft**: Neuro-Audio, **Oscilla**: TSM 300/400 + T800 series. All **GSI** tympanometers from GSI 20 – 39 + TymStar, Allegro.



Mushroom ear tips

Part no.	Description	Unit
8518409	AKR Assortment Box, various sizes	
8500570	3-6 mm (yellow)	Bag/100 pcs
8107510	4 mm (red)	Bag/100 pcs
8107511	5 mm (blue)	Bag/100 pcs
8500560	6 mm (green)	Bag/100 pcs
8500550	7 mm (yellow)	Bag/100 pcs
8500540	8 mm (red)	Bag/100 pcs
8500530	9 mm (blue)	Bag/100 pcs
8500520	10 mm (green)	Bag/100 pcs
8500510	11 mm (yellow)	Bag/100 pcs
8500500	12 mm (red)	Bag/100 pcs
8500490	13 mm (blue)	Bag/100 pcs
8107520	14 mm (green)	Bag/100 pcs
8500480	16 mm (yellow)	Bag/100 pcs
8505657	Mushroom ear tips, 5mm (blue)	Bag/25 pcs
8505658	Mushroom ear tips, 6mm (green)	Bag/25 pcs
8505659	Mushroom ear tips, 7mm (yellow)	Bag/25 pcs
8505660	Mushroom ear tips, 8mm (red)	Bag/25 pcs
8505661	Mushroom ear tips, 9mm (blue)	Bag/25 pcs
8505662	Mushroom ear tips, 10mm (green)	Bag/25 pcs
8505663	Mushroom ear tips, 11mm (yellow)	Bag/25 pcs
8505664	Mushroom ear tips, 12mm (red)	Bag/25 pcs
8505665	Mushroom ear tips, 13mm (blue)	Bag/25 pcs
8505666	Mushroom ear tips, 14mm (green)	Bag/25 pcs
8505667	Mushroom ear tips, 16mm (yellow)	Bag/25 pcs



DISPOSABLE SANIBEL™ AOT and ANA ear tips



VOLUME DISCOUNT AVAILABLE
for minimum orders of 10 bags

ANA ear tips for Echo-Screen®

For use with **Natus**: Echo-Screen®. **Madsen/Otometrics**: AccuScreen® Classic. **Interacoustics**: Eclipse with pink probe tip

Part no.	Description	Unit
8106363	4.0 mm (blue)	Bag/50 pcs
8106362	4.5 mm (red)	Bag/50 pcs
8106364	5.0 mm (green)	Bag/50 pcs
8106366	Infant tree (yellow)	Bag/50 pcs
8106365	Mushroom (green)	Bag/50 pcs



DISPOSABLE SANIBEL™ AOT and ANA ear tips



**VOLUME
DISCOUNT
AVAILABLE**
for minimum
orders of
10 bags

AOT ear tips

For use with **Madsen/Otometrics** AccuScreen®, Biologic Bio-logic® AuDX®, **Path Medical**: all Sentiero products

Part no.	Description	Unit
8106367	3.7 mm (yellow)	Bag/50 pcs
8106368	4.0 mm (light green)	Bag/50 pcs
8106369	4.5 mm (light pink)	Bag/50 pcs
8106370	5.0 mm (light blue)	Bag/50 pcs
8106371	4-7 mm (transparent)	Bag/50 pcs



Grason & Associates OT & MO select series ear tips



OT select ear tips w/ triangle holes

Ear tips with triangle holes for use with **Otodynamics** UGD, SND, SGD probes

Part no.	Description	Unit
8102644	4.8 mm (blue)	Bag/100 pcs
8102645	6 mm (green)	Bag/100 pcs
8102646	8 mm (red)	Bag/100 pcs
8102647	10 mm (blue)	Bag/100 pcs
8102648	12 mm (green)	Bag/100 pcs
8102649	14 mm (red)	Bag/100 pcs
8102650	16 mm (blue)	Bag/100 pcs



OT select ear tips w/two probe holes

Ear tips with two probe holes for use with **Otodynamics** ECP, UGS, SNS, SGS probes.

Part no.	Description	Unit
8102632	3.5 mm (blue)	Bag/100 pcs
8102633	4.5 mm (blue)	Bag/100 pcs
8102634	5.5 mm (blue)	Bag/100 pcs
8102635	6.5 mm (blue)	Bag/100 pcs
8102636	4.8 mm (blue)	Bag/100 pcs
8102637	6 mm (green)	Bag/100 pcs
8102638	8 mm (red)	Bag/100 pcs
8102639	10 mm (blue)	Bag/100 pcs
8102640	12 mm (green)	Bag/100 pcs
8102641	14 mm (red)	Bag/100 pcs
8102642	16 mm (blue)	Bag/100 pcs



MO select ear tips for Madsen OTOflex

For **MADSEN®** OTOflex™.

Part no.	Description	Unit
8102653	3-6 mm (yellow)	Bag/100 pcs
8102654	6 mm (red)	Bag/100 pcs
8102655	7 mm (blue)	Bag/100 pcs
8102656	8 mm (green)	Bag/100 pcs
8102657	9 mm (yellow)	Bag/100 pcs
8102658	10 mm (red)	Bag/100 pcs
8102659	11 mm (blue)	Bag/100 pcs
8102660	12 mm (green)	Bag/100 pcs
8102661	13 mm (yellow)	Bag/100 pcs
8102662	14 mm (red)	Bag/100 pcs
8102663	16 mm (blue)	Bag/100 pcs
8102664	18 mm (green)	Bag/100 pcs



Ear tips for OAE screeners

For use with **Biologic AuDX**, **Welch Allyn** OAE screener and **Natus** Echo-Screen® III.



OAE ear tips, AuDX

Part no.	Description	Unit
8107885	AuDX Preemie 3 mm	Bag/100 pcs
8515505	AuDX Preemie 3 mm	Bag/20 pcs
8107889	AuDX Tree 3.5-5mm	Bag/100 pcs
8515506	AuDX Tree 3.5-5mm	Bag/20 pcs
8107887	AuDX Infant 6 mm	Bag/100 pcs
8515507	AuDX Infant 6 mm	Bag/20 pcs
8102681	Adult OAE foam	Bag/100 pcs
8515436	Adult OAE foam	Bag/20 pcs
8102682	Pediatric OAE foam	Bag/100 pcs
8515435	Pediatric OAE foam	Bag/20 pcs
8103123	Jumbo OAE foam	Bag/100 pcs
8515438	Jumbo OAE foam	Bag/20 pcs



Snap probe accessories

Accessories specifically designed to work with SnapPROBE®. For use with **Interacoustics** Sera, **GrasonStadler** Novus, **Maico** EasyScreen.

Part no.	Description	Unit
8522758	Ear tip Infant A (Conical)	Bag/100 pcs
8522759	Ear tip Infant A (Conical)	Bag/25 pcs
8522760	Ear tip Infant C (Narrow)	Bag/100 pcs
8522761	Ear tip Infant C (Narrow)	Bag/25 pcs
8516436	Probe Waxguard	1 pc
8517889	Probe Waxguard	Set 10 pcs
8522461	Assortment of accessories including Infant tip A and C and waxguard	





We know that, first and foremost, great quality is at the top of your list, and you'll find it in each and every purchase you make from Sanibel Supply.

Bringing you the products you need is what drives us every day.

If we can ever be of assistance, don't hesitate to let us know.

Jakob Skovgaard
Business Unit Manager



Probe tips



IMP

Probe tip for Sanibel™ ADI™ ear tips

Fits all Otowaves, new MT10 and GSI Allegro.

Part no.	Unit
8002592	1 pc



ABR

IMP

OAE

Probe tip for Sanibel™ ADI™ ear tips

Fits Titan™ and easyTymp™.

Part no.	Unit
8001701	1 pc



IMP

OAE

Probe tip kit

Fits AT235™ or newer, AA222™, Sera™, Tympstar Pro™ Novus™, easyScreen™, and touchTymp™. Includes probe tip and gasket.

Part no.	Unit
8107762	1 pc



OAE

Internal probe tip replacements for ERO•SCAN™

Fits ERO•SCAN™.

Part no.	Unit
8102339	Bag/4 pcs



OAE

ERI0D probe tips for ERO•SCAN™ Pro and OtoRead™

It is recommended to clean the probe tip after every patient and to change the probe tip after 20 tests (40 ears). No re-calibration necessary.

Part no.	Unit
8012836	Bag/4 pcs



OAE

Internal probe tip replacements for ERO•SCAN™ Pro

Fits ERO•Scan™ Pro.

Part no.	Unit
8110364	Bag/4 pcs



OAE

Probe tip kit for OAE Ø3.3 mm ear tips

Fits all Eclipse™ – use with ADI™ ear tips. Requires no re-calibration.

Part no.	Unit
8012960	Bag/8 pcs



IMP

Sanibel™ probe tip adaptor kit

Converts AT235™, AT235h™, AA220™, AA222™ (blue/new), MI24™, MI26™, MI34 old series, MI44™ for use with ADI™ ear tips. Includes probe tip, washer, damper. No re-calibration required.

Part no.	Unit
8500430	Set/3 pcs



OAE

Remote probe tip replacements for ERO•SCAN™

Fits ERO•Scan™ Classic.

Part no.	Unit
8102340	Bag/4 pcs





Probe tips




100 pcs

OAE

Sanibel™ probe tubes for ADI™ ear tips

Fits Corti™, ERO•SCAN™ Plus and the new OtoRead™.

Part no.	Unit
8104159	Box/100 pcs




OAE

Probe Tube Removal Tool

Fits Corti™, ERO•SCAN™ Plus and the new OtoRead™.

Part no.	Unit
8513131	Bag/1 pc



1 pc +
100 pcs

OAE

Probe Tube Removal Tool kit

Fits Corti™, ERO•SCAN™ Plus and the new OtoRead™.

Part no.	Unit
8513132	Bag/100 pcs + 1 pc

Sanibel™ Probe Tube Removal Tool


- A = Probe
- B = Probe Tube Installer
- C = Attachment to probe cable
- D = Tube Remover

To Remove a Probe Tube

1. Place the front hole of the Probe Tube Tool over the end of the probe, as shown. The probe should be seated against the face of the tool and snap in place.
2. Squeeze the tool closed and hold it.
3. Twist the tool a couple of times while holding the tool closed and pulling away from the probe.
4. The probe tube will pull out from the probe. Discard the probe tube.

To Install a New Probe Tube

5. Place a new tube into the hole on the top of the tool.
6. Align the probe with the tube.
7. Push the probe onto the probe tube until the probe contacts the face of the probe tube tool and a snapping sound is heard.
8. Slide the tool off, leaving the tube in place in the probe.





**VOLUME
DISCOUNT
AVAILABLE**
for minimum
orders of
10 bags

ABR AUD

3A E-A-RLINK™ foam ear tips

Complies to all insert phones.

Part no.	Description	Unit
8500090	Adult/Medium size 13.7 mm	Bag/50 pcs
8106226	Adult/Medium size 13.7 mm	Bag/100 pcs




ABR

E-A-RLINK™ infant ear tips

For use with infants or in very small ear canals.

Part no.	Description	Unit
8500300	3.5 mm	Bag/20 pcs
8500310	4.0 mm	Bag/20 pcs



ABR

Sanibel™ IP30™ tube connectors

For use with IP30™ and E-A-Rtone® 3A.

Part no.	Unit
8107391	Bag/10 pcs



**VOLUME
DISCOUNT
AVAILABLE**
for minimum
orders of
10 bags

ABR AUD

3B E-A-RLINK™ foam ear tips

Complies to all insert phones.

Part no.	Description	Unit
8500110	Pediatric/Small size 9.7 mm	Bag/50 pcs
8106227	Pediatric/Small size 9.7 mm	Bag/100 pcs




ABR AUD

Sanibel™ ADI ear tip adaptor

Makes Sanibel™ ADI ear tips compatible with 3A insert phones.

Part no.	Unit
8500580	Bag/6 pcs



ABR AUD

Tube nipples for E-A-Rtone® 3A

For use with sound tubes 3A.

Part no.	Unit
8500140	Bag/10 pcs




ABR AUD

3C E-A-RLINK™ foam ear tips

Complies to all insert phones.

Part no.	Description	Unit
8500100	Jumbo/Large size 17.8 mm	Bag/24 pcs




ABR AUD

Sound tubes for E-A-Rtone® 3A ABR 50 ohm

Includes adaptor nipples, 1 pc blue and 1 pc red 50 ohm.

Part no.	Unit
8500720	Set/2 pcs



ABR AUD

Sanibel™ sound tubes and Sanibel™ ADI™ ear tip adaptor kit

Fits E-A-Rtone® 3A inserts and IP30™ insert phones. Includes ADI™ ear tip adaptors, 1 pc blue and 1 pc red 50 ohm.

Part no.	Unit
8100590	Set/2 pcs




ABR

Sound tubes for E-A-Rtone® 3A 10 ohm

Includes adaptor nipples.

Part no.	Unit
8107396	Bag/4 pcs




ABR

Sound tubes for E-A-Rtone® 3A 10 ohm

Includes adaptor nipples.

Part no.	Unit
8103605	Bag/4 pcs



ABR

Sanibel™ IP30 and E-A-Rtone® sound tubes and nipples

Fit all IP30™ and E-A-Rtone® versions. Includes tube connectors, 1 pc blue and 1 pc red 50 ohm.

Part no.	Unit
8501072	Set/2 pcs



ABR

E-A-Rtone® replacement fastener clip

Clip for fastening insert phone on clothes.

Part no.	Unit
8000744	1 pc



ABR

Sanibel™ fastener kit for IP30™

Includes velcro neckstrap and clip.

Part no.	Unit
8507410	Set



ABR OAE

Sanibel™ ear tip adaptor cleaning tool

Small brush for cleaning the Sanibel™ probe tip adaptors.

Part no.	Unit
8500830	Bag/10 pcs

100 tests

ABR

ABR consumable kit

Includes electrodes, insert phone tips, skin scrub and gaze.

Part no.	Description	Unit
8109670	Adult 3A	Kit/100 tests
8109671	Pediatric 3B	Kit/100 tests



2

LATEX

ABR

Sanibel™ Infant EarCup™ with snap or tab electrodes

Reduce test time! Test on both ears in one screening. Compatible with Titan™, Eclipse™, AccuScreen® (1st generation), MB II Classic and Echo-Screen®. Use only with equipment calibrated for over-ear style disposables. Disposable over-the-ear earphones are high quality, low cost, and ideal for bilateral newborn hearing screenings. 1 set consists of 2 Infant EarCups™, 3 snap or tab electrodes and 1 alcohol pad.

For use with a range of different providers of newborn hearing screening equipment:

Grason Stadler: Novus. **Interacoustics:** Titan, Sera. **Maico:** easyScreen, MB II Classic.

Madsen/Otometrics: Accuscreen. **Natus:** Biologic Abear, Echo-Screen.

Part no.	Description	Unit
8049600	Tab	Box/20 sets
8049700	Snap	Box/20 sets
8508848	w/o electrodes	Box/20 sets

Sanibel Infant EarCup™ w/ Hydrogel

New Hydrogel version to ensure a more comfortable solution for the baby while testing.

Choose any electrode type you wish to add to meet your preference.

Order electrodes separately.

Part no.	Unit
8516095	Box/20 sets

Exchange the electrodes with your preferred version!
CONTACT US FOR SUPPORT




ABR

Sanibel™ Infant EarCup™ adaptor set

Connects insert phones with Infant EarCups™. NB: Do not contain filters.

Part no.	Unit
8500450	Set/2 pcs



ABR

Sanibel™ sound tubes and Sanibel™ Infant EarCup™ adaptor kit

Includes adaptors, 1 blue and 1 red tube.

Part no.	Unit
8100577	Set/2 pcs

20 sets



The “REM effect”
– better fitting & fewer follow-ups!

Be sure of an
individual and
ideal fitting

Assurance of your
fitting gives you
confidence in your work

Confidence increases
your credibility and
client satisfaction

Client satisfaction means
**fewer follow-up
appointments**

Use Sanibel™ probe tubes for any equipment

We often hear from the industry that only the manufacturers' specific probe tubes should be used with their specific equipment. However, a recent comparative study by Sanibel Supply® disproves this myth.

In this study, our in-house audiologists and engineers compared the effects of different probe tubes on an Affinity^{2D} system. They compared AudioScan tubes (red), Otometrics tubes (blue), Sanibel™ tubes (orange) and Sanibel™ tubes (grey).

This involved running each tube via the 'tube calibration' process and then comparing the following REUG curve when the tube had been left in the calibration position.

The investigation showed no significant difference in the REUG curve and all probe tubes achieved the same result when the calibration process had been run correctly.




HIF

Sanibel™ REM probe tubes, type A

Universal fit, length 76 mm, beveled 45° angle tip and black marker ring.

Part no.	Unit
8104116	Tray/50 pcs



HIF

Sanibel™ REM probe tubes, type B

Fit all brands, **Interacoustics**, **Natus Aurical**, **MedRx** and **AudioScan**.

Part no.	Unit
8500360	Bag/40 pcs



50 pcs

HIF


Silicone tubes for SPL60™

For use with Affinity™ REM probe.

Part no.	Unit
8500700	Bag/50 pcs




Cerumen Management accessories



Bionix® Safe Ear Curettes

Bionix® Safe Ear Curettes allow healthcare providers to safely and gently remove cerumen with styles to match different types of wax and personal preference. Bionix® Safe Ear Curettes are flexible, to reduce injury to the ear canal, making them safer than stainless steel curettes. Seven color coded tip styles to meet specific needs, giving you “the right tool for any job” for every cerumen situation. Also available as a variety pack.

Part no.	Description	Unit
8515585	Green MicroLoop®	50 pcs/box
8515586	Orange ControLoop®	50 pcs/box
8515587	Red AngleLoop®	50 pcs/box
8519359	Purple VersaLoop®	50 pcs/box
8519360	Blue InfantScoop®	50 pcs/box
8519372	Yellow CeraSpoon®	50 pcs/box
8519361	White FlexLoop®	50 pcs/box
8519362	Variety Kit, (10 each style & 15 white)	



Bionix® Lighted Ear Curettes

Bionix® Lighted Ear Curettes combine illumination and magnification to help healthcare providers take less time to see and safely remove cerumen from the ear canal. The Lighted Ear Curette with magnification brilliantly illuminates the ear canal, eliminating the practice of blindly curettaging the sensitive area. Complete visualization shortens procedure time and increases efficacy. Box of 50 pcs.

All versions include 1 Light Source and 1 magnifying lens.

Part no.	Description
8519364	Lighted AngleLoop®, 4mm diameter
8519365	Lighted MicroLoop®, 3mm diameter
8519366	Lighted WaveCurette™, 2.5mm wide tip
8519367	Lighted FlexLoop®, 4mm wide tip
8519368	Lighted CeraSpoon®, 4mm wide tip
8519369	Lighted InfantScoop®, 2mm wide tip
8519373	Lighted VersaLoop®, 3mm wide tip
8519370	Lighted CeraPik™, 2.5mm concave tip. Great for hard, impacted cerumen and foreign bodies.



Lighted Ear Curette Variety Pack

Part no.	Description
8519363	With 1 light source and magnifying lens includes all tip styles except CeraPik™.



Lighted Forceps

Bionix Lighted Forceps combine illumination and magnification to help healthcare providers take less time to see and safely remove ear canal and nasal passage obstructions.

Part no.	Description
8519392	10 forceps, 1 light and 1 magnifying lens



OtoClear® Spray Wash Kit

OtoClear Ear Irrigation system gently directs water to the ear canal wall with three divergent streams, instead of directly at the Tympanic Membrane.

Part no.	Description
8519371	kit of 20x irrigator tips, 1x 470ml spray wash bottle, 1x ear basin
8519390	Otoclear Ear Basin, Box/3 pcs

Otoclear Irrigation Tips

Bionix OtoClear Ear Irrigation helps provide safe and comfortable ear irrigation without the backsplash, at half the cost of a self-constructed device.

Part no.	Unit
8519375	Box/40 pcs





Micro Suction Products


Otopront VAC 18 is a transportable suction unit, specific designed for ear wax removal or other temporary suction purposes within hospitals and clinics. Available in selceted coutries only, please contact us for details.

Part no.	Description
8515194	Otopront VAC 18 suction unit

Accessories for Otopront Suction units

To optimize the use of the Otopront VAC 18. Please use the appliances below per demand.

Part no.	Description
8517966	Suction container system (to be added to the VAC 18)
8515200	Refill bags for disposables suction system, box/60 pcs, Value box of refill bags
8515372	Refill bags for disposables suction system, 10 pcs
8520345	Silicon suction hose medical grade, 50 m
8515195	Air filter VAC 18
8515197	Internal filter VAC 18 / VAC 28
8515199	Lid for suction container VAC 18 / VAC 28
8522023	Turmakin Ear Speculum 5,5 mm/1 pc
8522024	Turmakin Ear Speculum 7,0 mm/1 pc
8522025	Suction cannula, 2mm x 80mm/1 pc






Training
kit
available
upon
request

2


EarWay Pro – Wax Removal Helix

Very easy and safe way to remove ear wax. The helix grabs around the ear wax lump and pulls it out in a comfortable way. Available in different languages, please contact us for details.


Part no.	Description	Unit
8518025	Pro M, 5,4 mm, purple	25 units
8518024	Pro S, 4,2 mm, blue	25 units



Before insertion



Engulfing of Cerumen



Extraction





Impression materials and supplies



HIF

Impression pads, standard, pink

Diameter approx. 1.1 cm, length 1,2 cm.

Part no.	Unit
8515559	Box/100 pcs



HIF

Impression pads CIC, Yellow

Diameter approx. 1.1 cm, length 0,7cm.

Part no.	Unit
8522182	Box/100 pcs




HIF

Foam Impression Pads, conical, light pink

Diameter approx. 1.25 cm.

Part no.	Unit
8515560	Bag/100 pcs




HIF

Foam Impression Pads Slim, small, grey

Diameter approx. 0.95 cm.

Part no.	Description	Unit
8515429	Grey	Bag/50 pcs




HIF

Foam Impression Pads Slim, medium, white

Diameter approx. 1.1 cm.

Part no.	Description	Unit
8515430	White	Bag /50 pcs



HIF

Foam Impression Pads Slim, large, blue

Diameter approx. 1.25 cm.

Part no.	Description	Unit
8515423	Blue	Bag/50 pcs



HIF

Cotton Impression Pads, small with blue safety thread

Diameter approx. 0.65 cm.

Part no.	Description	Unit
8515571	Blue	Bag/50 pcs



HIF

Cotton Impression Pads, medium with white safety thread

Diameter approx. 0.75 cm.

Part no.	Description	Unit
8515572	White	Bag/50 pcs



HIF

Cotton Impression Pads, large with red safety thread

Diameter approx. 0.85 cm.

Part no.	Description	Unit
8515573	Red	Bag /50 pcs

Different varieties available, please ask us for details



Otoform® A flex

Soft, flexible, easy processing.
Recommended for CIC impressions, narrow auditory canals or people sensitive to pain.
Final hardness approx. 30Shore A.

Note: requires an impression gun and mixing tip.

Part no.	Description	Unit
8521806	48 ml double cartridges	1 pc
8522433	50 ml double cartridges	1 pc



Manual Impression Gun

Part no.	Description	Unit
8522108	Injection for 48 ml double cartridges	1 pc
8522434	Injection for 50 ml double cartridges	1 pc



Mixing tray with dosing scale

Part no.	Unit
8522105	1 pc



Mixing Tips for 48 ml double cartridge

Use with 8521806.

Part no.	Description	Unit
8522430	5.4mm	Bag/100pcs



HIF

Impression material A-ZOFT

Premium quality and smooth initial consistency. Extremely soft, easy to remove and allows highest impression precision.
Final hardness approx. 30Shore A.

Part no.	Unit
8515424	2x800 g jars
8515426	40 foil bags (sets 2x8 g)



Mixing pad

Part no.	Unit
8522106	Bag/100 sheets



Mixing Tips for 50ml double cartridge

Use with 8522433.

Part no.	Description	Unit
8522432	5,4mm	Bag/100pcs



Large Mixing Spatula

Part no.	Unit
8522104	1 pc



Impression Syringe, grey

Easy to fill, inject and clean. Provides a stable and easy injection. 4 mm.

Part no.	Unit
8515445	1 pc





Earlight with Scaled Top

Part no.	Description
8519357	Luxamed Ear light with acrylic scaled top, black
8519358	Luxamed Acrylic Top with Scale



Stethoclip

Hearing Instrument listening stethoscope for reassuring the fitting process

Part no.	Description	Unit
8515600	Stethoclip, black	1 unit
8520346	Ear funnels	1 unit
8520347	ITE Adaptor	1 unit



Surgical Blades and Handles

Suitable for cutting HA tubing.

Part no.	Description	Unit
8515427	Surgical blade handle No.3	1 pc
8515428	Surgical blades No11	100 pcs



Libby Horn

Special tube designed to increase gain over conventional tubing. Ideal when using an analog hearing instrument and fitting a high frequency hearing loss. Also designed to allow for a friction fit into earmold without the need to damage soundbore when removing it. Colour.

Part no.	Description	Unit
8515663	4 mm,	Bag/10pcs
8515657	3 mm,	Bag/10pcs



Mini Cleaning Brush

Clean hearing instruments in an easy way.

Part no.	Unit
8519100	1 pc



Loctite 406 Instant adhesive

An instant adhesive designed especially for the fast bonding of plastics and elastomeric materials where very fast fixturing is required. This product offers exceptionally low viscosity.

Part no.	Unit
8005067	1 unit



Hearing aid maintenance and hearing care

Take care of your hearing aids

Take good care of your hearing aid. Your Hearing instrument is a high tech device needing regular care in order to provide you the best performance every day.



Cleansing tablets

Part no.	Unit
8520148	Pk/4 pcs



Dehydrating capsules

Part no.	Unit
8520149	4 pcs



Drying cup

Part no.	Unit
8520150	1 pc



Cleaning cup

Part no.	Unit
8520151	1 pc



Otofloss

Part no.	Unit
8520152	1 pc



Perfect Dry LUX box for all types of hearing aids

Part no.	Unit
8520154	1 pc



Starter kit

Part no.	Unit
8520155	1 kit



Drying spray

Part no.	Unit
8520153	1 pc



Earmold Air Blowers

Part no.	Unit
8515562	1 pc





Hearing aid maintenance and care



Cleaning spray

Part no.	Unit
8520147	1 pc



Cleaning wipes

Part no.	Unit
8520146	Pk/30 pcs



Protect your ears

Part no.	Description	Unit
8520052	Natural lubricant	1 pc



Protect your ears

Part no.	Description	Unit
8520053	Dermoprotective Gel	1 pc



Protect your ears

Part no.	Description	Unit
8520063	Earspray	1 pc



Available in different languages, ask us for details

Sanohra Earplugs - Fly

Sanohra Fly earplugs, suitable for children and adults. Reduces ear pressure experiences when flying.

Part no.	Description	Unit
8519351	EarPlugs Fly Kids, EN	1 set
8519352	EarPlugs Fly Adult, EN	1 set



Available in different languages, ask us for details

Sanohra Earplugs - Swim

Sanohra Swim earplugs, suitable for children and adults. Prevent water from getting into the ear canal and protect ears when exposed to water.

Part no.	Description	Unit
8519353	EarPlugs Swim Kids, EN	1 set
8520143	EarPlugs Swim Adult, EN	1 set

SANOHR fly ear plugs

Beyond "noise protection" it can be important to protect the ear drum (tympanic membrane) from changes in external pressure that are too rapid, e.g., during take off and landing phases in airplanes. For airplane passengers experiencing a problem with the pressure equalisation, the use of a special ear protection product means freedom from an unpleasant feeling or piercing pain during flight.

SANOHR swim ear plugs

In case of injury to the ear drum or insertion of a vent tube, it is necessary to prevent water from reaching the ear drum through the outer auditory canal. Swimming, bathing or showering might cause inflammation of the middle ear. Staying in chlorinated water for an extended period of time can cause the skin in the outer auditory canal to swell and cause inflammation (otitis externa).



Disposable electrodes

Large, soft
electrode size!

60 pcs



ABR

Sanibel™ snap electrodes

Pre-gelled with a very gentle hydrogel, perfect for ABR testing on infants. Size 2.2 x 3 cm. Item no. 8107137 is packed in 20 smaller pouches of 3 electrodes. Intended use areas: Newborn, diagnostic follow-up, ABR, C-Vemp.

Part no.	Unit
8107137	Box/60 pcs (20x3 pcs)



Large, soft
electrode size!

60 pcs



ABR

Sanibel™ tab electrodes

Pre-gelled cloth electrode with a very gentle hydrogel, perfect for ABR testing on infants. Size 2.2 x 3 cm. Packed in 20 smaller pouches of 3 electrodes. Intended use areas: Newborn, diagnostic follow-up, ABR, C-Vemp.

Part no.	Unit
8107410	Box/60 pcs (20x3 pcs)



100 pcs



ABR

Sanibel™ vinyl tab electrode

For alligator clip cable. Solid gel. Size 2.5 x 2.2 cm. Packaged in convenient re-sealable pouches.

Part no.	Unit
8107426	Bag/100 pcs (10x10 pcs)



20 pcs

ABR

E-A-RLINK™ gold foil tip trodes

For ECoG. Use with cable 8011343.

Part no.	Description	Unit
8500370	Adult 13 mm	Bag/20 pcs
8500440	Pediatric 10 mm	Bag/20 pcs



10 pcs




ABR

Sanibel™ TM Electrode

For ECoG. Use with cable 8500650.

Part no.	Unit
8011489	Box/10 pcs



Sanibel Supply® offers various kinds of surface electrodes for ECG as well as EEG and ABR. Surface electrodes vary in size, adhesion, connectivity and frequency (test category).

The frequency plot to the above shows the electrical stimuli sent directly through the system and compared to

Choose the correct
electrode for optimal
performance

the surface electrode (here ABR). The two lines are similar and demonstrate that the surface electrode can be used for ABR testing.

At Sanibel Supply®, we always test and validate our surface electrodes to ensure optimal performance for your testing!

50 pcs



ABR

Vermed Economy line 3-pack solid gel snap Electrode

Solid gel type snap electrode suitable for newborn screening and diagnostic follow-up. Latex free. Intended use areas: Newborn, diagnostic follow-up all ages.

Part no.	Unit
8516252	Bag/50 pcs





Disposable electrodes



75 pcs

ABR

Vermed Wetgel Snap Electrode

Pre-gelled hydrogel electrode with Ag/AgCl sensor. Very suitable for newborn screening or diagnostic follow up. Packaged 3 pcs. per pouch. Intended use areas: Newborn, diagnostic follow-up, ABR, C-Vemp.

Part no.	Unit
8516221	Box/75 pcs



60 pcs

ABR

Vermed Hydrogel Snap Electrode

Smaller than the rectangular version, but otherwise similar. Good, soft electrode size making it very suitable for newborn screening and ABR measurements. Cloth material makes it convenient to wear. Intended use areas: Newborn, ABR.

Part no.	Unit
8102022	Bag/60 pcs (3x20 pcs)



300 pcs

ABR

Vermed Prewired neonatal eeg

Electrode for use in pediatric and neonatal settings. 3 electrodes per pouch. Square shape 22.2 x 22.2 mm. Intended use areas: Newborn, ABR.

Part no.	Unit
8107425	Bag/300 pcs (3x100 set)



100 pcs

ABR

Vermed Tab Electrode

Rectangular tab electrode for newborn screening. 12.7 x 22.9 mm. Conductive Adhesive Hydrogel. 10x10 pcs. per pouch. Intended use areas: Newborn, O-Vemp.

Part no.	Unit
8516225	Bag/100 pcs



100 pcs

ABR

Vermed Tab Electrode

Rectangular tab electrode for newborn screening. 22.9 x 34 mm. Conductive Adhesive Hydrogel. 10x10 pcs. per pouch. Intended use areas: Newborn, O-Vemp.

Part no.	Unit
8516227	Bag/100 pcs



60 pcs

ABR

Vermed Hydrogel Snap Electrode

Large soft electrode size making it very suitable for newborn screening and ABR measurements. Cloth material makes it convenient to wear. Size 3.4 x 2.2cm. Intended use areas: Newborn, ABR.

Part no.	Unit
8516907	Bag/60 pcs (3x20 pcs)



60 pcs

ABR

Vermed Wet Gel foam

Electrode with an Ag/AgCl sensor. 10x6 pcs per pouch. Diameter 38.1 mm

Part no.	Unit
8108563	Bag/60 pcs




50 pcs

ABR

Vermed Economy Line Snap Wet Gel Electrode

Non-woven cloth electrode with conductive wet gel. Passed ISO 10993 and is Latex free. Intended use areas: Newborn, diagnostic follow-up all ages.

Part no.	Unit
8516251	Bag/50 pcs



5 pcs

ABR

Vermed Allround ECG electrode

Provides very good impedance and it well suited for both newborns and adults. Can be used both as snap and tab electrode with alligator cables. Pouches come in 3 pcs. Value packs available in 600 pcs. per order. Intended use areas: Newborn, diagnostic follow-up all ages.

Part no.	Unit
8516906	Bag/5 pcs





100 pcs

ABR

Ambu® 0415M tab electrode

Medium size tab electrode with solid gel for 12-lead ECG. Presented conveniently with 10 pieces to a card. Fish tail tab design reduces electrode lift-offs. Suitable for both adults and children. Intended use areas: ABR diagnostics of all kinds with Tab electrodes.

Part no.	Unit
8516138	Bag/100 pcs (10x10 pcs)



1000 pcs

ABR

Ambu® 2300 Blue Sensor Tab electrode

The medium size Ambu BlueSensor 2300 electrode features a skin-friendly solid gel that provides optimal signal quality recordings during 12-lead diagnostic testing ECG application. Very suitable for ABR diagnostics of all kinds. Intended use areas: ABR diagnostics of all kinds.

Part no.	Unit
8516139	Box/1000 pcs (10x100 pcs)




25 pcs

ABR

Ambu® Blue Sensor™ N-00-S

Foam ECG electrode for use with button cables. Wetgel. Intended use areas: ABR diagnostics of all kinds.

Part no.	Unit
8013029	Bag/25 pcs



ABR

Ambu® Ground electrode 71410-M/1

Part no.	Unit
8500690	1 pc



25 pcs

ABR

Ambu® Neuroline 720

Ambu Neuroline 720 electrodes are pre-gelled and self adhesive. They are ideal for use in sleep registration where the strong, yet skin-friendly, adhesion ensures that the electrode stays in place for the entire examination, even if the patient sweats. Intended use areas: ABR diagnostics of all kinds.

Part no.	Unit
8500060	Bag/25 pcs




50 pcs

ABR

Ambu® Blue Sensor™ P-00-S

High quality Ag/AgCl sensor – ensures excellent trace quality. Intended use areas: ABR diagnostics of all kinds.

Part no.	Unit
8504782	Bag/50 pcs



ABR

Ambu® Neuroline 720 – with 100 cm cable and k-connector

Ambu Neuroline 720 electrodes are pre-gelled and self adhesive. They are ideal for use in sleep registration where the strong, yet skin-friendly, adhesion ensures that the electrode stays in place for the entire examination, even if the patient sweats. Intended use areas: ABR diagnostics of all kinds.

Part no.	Unit
8503244	Bag/10 pcs



25 pcs

ABR

Ambu® Blue Sensor™ N-10-A

Foam ECG electrode with 10 cm wire and 4 mm connector, for use with universal electrode cable. Wetgel.

Part no.	Unit
8500670	Bag/25 pcs





Surface electrodes tailored to your need

At Sanibel Supply® we believe that it is important to choose the correct surface electrode for each specific recording and patient type.

Wetgel or hydrogel electrodes?

A wetgel electrode consists of water, thickener, ionic salts, surfactants and bactericides or fungicides. Wetgel migrates easily and rapidly onto the skin surface, reducing impedance levels, which in turn will generally result in lower noise in the recording.

A hydrogel surface electrode is less adhesive when compared to wetgel, making it particularly suitable for newborns.

Chloride salt concentration

The higher the chloride salt content, the more conductive the electrode. Higher salt content (10% and above), pre-gelled surface electrodes are useful for making fast, high quality measurements of biopotentials. High conductivity electrodes generally have reduced artifact, due to the lower impedance generated between electrode and skin surface. As the chloride salt content drops, there is a corresponding reduction in electrode conductivity.

If the chloride content is below 10%, the electrode can then be employed for longer recording times, with a reduced risk of skin irritation. As a rule of thumb, wetgel has a chloride content of 10% and hydrogel 4%.

It is important to use electrodes with similar chloride salt content when testing.

The impedance of the electrode/skin junction is highly dependent on the electrolyte type and the chloride salt concentration. For example, a hydrogel electrode with 4% chloride concentration will have approximately 10x higher impedance than a wet liquid gel electrode with 10% chloride concentration, after first application to the skin.

Adhesion

An advantage of a conductive adhesive is that the active surface can be as large as the electrode itself because no additional adhesive strip is needed. The Sanibel disposable electrodes come with a range of adhesive qualities. All electrodes are designed to adhere well to skin surfaces, but

the adhesion characteristics can be grouped according to the application:

- Strong adhesive
- Moderate adhesive

Strong adhesion is not recommended for newborns.

Skin preparation

For optimum electrode to skin conductivity, the skin should be lightly brushed with a gentle abrasive gel/wipe. **An alcohol wipe is not recommended to improve conductivity as this will only dry out the skin surface.**





ABR

Sanibel™ ECoG starter kit

For Eclipse™ and Bio-logic® AEP system. Includes cable, gel and 2 pcs TM electrodes.

Part no.	Unit
8013209	Set/4 pcs



ABR

Sanibel™ ECoG cable

For Eclipse™ and Bio-logic® AEP system. 1.5 mm/.060" female connector. For use with Sanibel™ TM electrodes.

Part no.	Unit
8500650	1 pc

After use ...

After use, cables should be wiped with a soft cloth that has been dipped in a mild detergent solution or alcohol.

Contact between the cleaning solution and the connectors should be avoided.




ABR

Sanibel™ EP electrode jumper cable

Allows two of the preamplifier inputs to be connected, e.g. for linking mastoid electrodes for cortical testing.

Part no.	Unit
8107872	1 pc




ABR

Disc electrodes

For use with Ten20™ electrode paste. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50 cm cable.

Part no.	Unit
8011498	Set/3 pcs RWB
8011345	Set/5 pcs RWBYK




ABR

Gold-plated cup electrodes

Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). Ø 10 mm gold-plated. 50 cm cable.

Part no.	Unit
8107418	Set/3 pcs RWB
8107145	Set/5 pcs RWBYK




ABR

Alligator clip cable

For use with tab electrodes. 1.5mm/.060" female connector. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K).

Part no.	Unit
8516581	Set/3 pcs RWB
8513776	Set/3 pcs WYK
8516584	Set/5 pcs RWBYK



ABR

Pinch style cable

For use with snap electrodes. 1.5mm/.060" female connector. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K).

Part no.	Unit
8500390	Set/3 pcs RWB
8500400	Set/3 pcs WYK
8500410	Set/5 pcs RWBYK




ABR

Tip trode cables

For use with TTE25 gold foil tip trodes. 1.5 mm with touch proof connectors and plated clamp clips. 50 cm long.

Part no.	Unit
8011343	Set/2 pcs



ABR

Universal electrode cable

Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50 cm cable. Non-tangling.

Part no.	Unit
8011353	Set/3 pcs RWB
8011344	Set/5 pcs RWBYK



ABR

Standard button cable

For use with snap electrodes. Available in Red (R), White (W), Blue (B), Yellow (Y), Black (K). 50 cm cable. Non-tangling.

Part no.	Unit
8011354	Set/3 pcs RWB
8011346	Set/5 pcs RWBYK



ABR

NuPrep™ skin prep gel, 114 g

Efficient, yet gentle skin abrasive. Chloride and acetone free. Pack of 1x4 oz/114 g tube or 3x4 oz/114 g tube.

Part no.	Unit
8500010	1 tube
8500320	Box/3 tubes



ABR

Lectron® II electrode gel, 250 ml

Highly conductive, chloride-free electrode gel. Recommended for use with Sanibel™ electrodes for ECoG.

Part no.	Unit
8102626	Tube/250 ml



ABR

Spectra 360® electrode gel, 250 g

Salt- & chloride-free conductive electrode gel suitable for long-term applications.

Part no.	Unit
8500740	Tube/250 g



ABR

Lemon Prep™ abrasive lotion, 114 g

Pack of 3x4 oz/114 g tubes.

Part no.	Unit
8500770	Box/3 tubes



ABR

SignaGel®, 250 g

Bacteriostatic conductive saline electrode gel.

Part no.	Unit
8102628	Tube/250 g



ABR

Ten20™ conductive EEG paste, 114 g/228 g

Adhesive gel for disc electrodes. Available in pack of 3x4 oz/114 g tubes, 3x8 oz/228 g jar or 1x8 oz/228 g jar.

Part no.	Unit
8012924	Box/3 tubes
8500070	Box/1 jar
8500330	Box/3 jars



ABR

Everi™ abrasive paste, 160 g

Reduces the skin impedance during EEG, EP, ECG and Holter Stress test examination.

Part no.	Unit
8012952	Bottle/160 g



ABR

Sonavelle® electrode gel, 250 ml

Formaldehyde-free and glycerol-free.

Part no.	Unit
8500380	Bottle/250 ml



ABR

3M™ Red Dot™ trace prep, 5 m

One-step skin preparation. 18 mm x 5 m / 0.70" x 196".

Part no.	Unit
8500030	Dispenser/5 m



12 rolls

ABR

3M™ Transpore™ Surgical tape, 9.14 m

Transparent, latex-free, water-resistant and porous for breathability. 25 mm x 9.14 m / 1" x 10 LY.

Part no.	Unit
8500750	Box/12 rolls



12 rolls

ABR

3M™ Micropore™ Surgical tape, 9.14 m

Latex-free and hypoallergenic. Holds well on damp skin. 25 mm x 9.14 m / 1" x 10 LY.

Part no.	Unit
8500760	Box/12 rolls



10 packs


ABR

OAE

Q-tips

Cotton tipped applicators for patient preparation. Non-sterile. Each pack contains of 100 pcs.

Part no.	Unit
8500810	Box/10 packs



ABR

OAE

B Braun® alcohol pads

70% isopropyl alcohol.

Part no.	Unit
8500050	Box/100 pcs



ABR

OAE

Gauze swabs

Non-sterile, 7.5 x 7.5 cm and 8 layered.

Part no.	Unit
8500040	Pack/100 pcs





2

LATEX

AUD



500 pcs

3" Disposable ear cushion cover

Replaces disinfectants. Headphone covers for classic audiometric headphones.

Part no.	Unit
8107419	Box/500 pcs



ABR IMP OAE

Probe cleaning tool

For **Interacoustics** Titan™/AT235™/AA222™/Sera™/MT10™/**Amplivox** Otowaves/**MAICO** touchTym™ MI22, MI23, MI24, MI26, MI34, MI36, MI44 + touchTym™ and **GSI** TymStar™ Pro, Novus and Allegro.

Part no.	Unit
8011351	1 pc

2

LATEX

AUD



100 pcs

4.5" Disposable ear cushion cover

Replaces disinfectants. Headphone covers for noise cancelling audiometric headphones. For HDA200/DD450/DD65v2/Peltor headsets.

Part no.	Unit
8513404	Bag/100 pcs



Disinfectant detergent, 1000 ml

Alcohol-free, ready to use solution for a quick and extensive disinfection of non-invasive medical devices.

Part no.	Unit
8501706	Bottle/1000 ml



Keyboard protective foil

Place a protective foil over your keyboard to minimize contamination with grime, bacteria and germs, or liquid.

Part no.	Unit
8500460	1 pc



ABR IMP OAE

Bridge & Implant floss

Used to clean probe tips.

Part no.	Unit
8003566	Box/30 pcs

2

How to clean probe tips!



ABR IMP OAE

3-in-1 floss

Used to clean probe tips.

Part no.	Unit
8003565	Box/100 pcs



Disposable head covers

Perform high hygiene practices in front of your patients or clients.

Part no.	Unit
8519350	Box/100 pcs

Disposable face masks

Perform high hygiene practices in front of your patients or clients.

Part No.	Description	Unit
8515555	Safety mask with loops	Box/50pcs

Disposable apron

Keep a professional appearance also after servicing a hearing instrument.

Part no.	Unit
8519348	Box/100 pcs

Disposable shoe covers

Perform high hygiene practices in front of your patients or clients.

Part no.	Unit
8519349	Box/100 pcs



Safety goggles

Part no.	Unit
8522417	1 pc



Nitrile examination gloves

For use at medical sites and at research institutes etc. Do not interact with petroleum oils, grease, solvents. Do not contain latex.



Part no.	Description	Unit
8515551	Small	Box/00 pcs
8515552	Medium	Box/100 pcs
8515553	Large	Box/100 pcs
8515554	Extra large	Box/100 pcs

Disinfection gel and wipes

Easy disinfection of skin and surfaces

Part no.	Description	Unit
8522407	Dax Alcogel 85	600 ml
8520146	Wipes canister	30 pcs

Infection control

A clean clinic is also part of your professional image

Infection control is so much more than cleaning probe tips, protecting your personal computer's keyboard, and maintaining high standards within hygiene using e.g. disinfectant detergent for a quick and extensive disinfection of non-invasive medical devices.

When it comes to infection control, preventing germs, dirt and dust from spreading via other sources, is also of critical importance – both to the safety of

your patients and yourself! Gloves, shoe covers, facemasks, and aprons are just some of the "accessories", which can help you succeed in this.

Maintain high, professional standards and make your clinic an ever-safe environment for your patients, making the wearing of gloves, facemasks and other body protection wear a natural part of the "patient interaction" process.






Parts for audiometry



AUD

USB cables A to B

Part no.	Description	Unit
8011241	2 m/6.56 feet	1 pc
8011243	5 m/16.40 feet	1 pc



AUD

TDH39™/DD45™ replacement cable, symmetric

For TDH39™ and DD45™. 30° jack plug.

Part no.	Description	Unit
8010822	2 m	1 pc
8507204	5 m	1 pc



AUD

HDA200 replacement cable, symmetric

For HDA200 and DD65 headphones. 2 m cable and 6.3 mm/ 30° jack plug.

Part no.	Unit
8004172	1 pc



AUD IMP

TDH39™/DD45™ contra headset replacement cable

For TDH39™ and DD45™. Single mono. 30° jack plug.

Part no.	Unit
8010826	1 pc



AUD

B71W Bone conductor replacement cable

2 m cable and 6.3 mm/ 30° jack plug.

Part no.	Unit
8011094	1 pc



AUD

RadioEar IP30 Insert earphone cable

Cable length: 2 m. Plug type: 30°, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8004432	1 pc



AUD

IP30 Contra Insert earphone cable

Cable length: 2 m. Plug type: 30°, 1x6.3 mm gold-plated mono jack.

Part no.	Unit
8107358	1 pc



AUD

RadioEar Foldable audiometric headband

Black, foldable.

Part no.	Unit
8010882	1 set



AUD

RadioEar RE-7 Headband

Black.

Part no.	Unit
8517379	1 set



AUD

Peltor® headset/earmuff replacement hygiene kit

For DD450, HDA200 and DD65.

Part no.	Unit
8500270	1 set



AUD

RadioEar DD45 Audiometric headset

Headband: HB7. Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8010954	1 set



AUD

RadioEar DD45 Adult audiometric headset

Headband: HB3045. Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8105268	1 set

RadioEar DD45 Adult audiometric headset

Headband: HBA, foldable. Transducer: DD45. Cable length: 2 m. Plug type: 30°, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8106340	1 set



AUD

RadioEar DD65 v2 Passive attenuation audiometric headset

Cable length: 2 m. Plug type: Straight, 2x6.3 mm gold-plated mono jack.

Part no.	Unit
8510117	1 set



AUD

RadioEar DD45 transducer

Audiometric earphones engineered to match the performance and function of the most popular earphones available for audiological testing. 10 ohm receiver. Easy installation. Reliable. RoHS compatible.

Part no.	Unit
8507268	1 pc



AUD

Sanibel™ TDH39/DD45 earphone cushion red & blue, model 51

For TDH39, 49, 50 & DD45 earphones.

Part no.	Unit
8107611	1 Set



AUD

DD 51 cushion

For DD45™ and TDH39™ headsets.

Part no.	Unit
8010894	1 pc



AUD

RadioEar DD45 Contra headband

Tension: 4,5 N ±0.5N

Part no.	Unit
8507211	1 set



AUD

Foam pad for contra headset

Compatible with e.g. TDH39c™ and DD45c™.

Part no.	Unit
8507292	Bag/10 pcs



AUD HIF ABR

RadioEar replacement pads (3288-1)

For RadioEar Spring Steel headband.

Part no.	Unit
8504777	Bag/10 pcs



AUD

Ear cushion HZP 09 pair

Ear cushion foam (1 pair) suitable for: PC 20, PC 25, PC 30, PC 31, PC 35, PC 36, PC 130, PC 131, PC 135, PC 136, PC 140, PC 141, PC 145, PC 146.

Part no.	Unit
8104416	1 pair



AUD

Headband cushion

For RadioEar HBA foldable audiometric headband.

Part no.	Unit
8010886	1 pc



AUD

RadioEar B71 Bone transducer headset

10 ohm. Cable length: 79"/200 cm molded cord. Plug type: 30°, 1x6.3 mm gold-plated mono jack. Headband: P3333

Part no.	Unit
8510370	1 pc



AUD HIF ABR

RadioEar Spring Steel headband (P3333)

Part no.	Unit
8507213	1 pc



AUD

Headband cushion

For Sennheiser HDA200, RadioEar DD450, DD65v2 and Peltor audiometric headbands.

Part no.	Unit
8000770	1 pc



AUD HIF ABR

RadioEar AMBAND bone conductor headband (P1188)

For adult, children and infants.

Part no.	Unit
8507214	Bag/10 pcs



Replacement cable for bone conductor B71/B81 (with safety fixation point)

Part no.	Unit
8515713	1 pc

New



24 pcs

BAL

Disposable VNG goggle foam pads

For **Interacoustics**® VNG mask.

Part no.	Unit
8011527	Box/24 pcs



New

Disposable specula for Difra Airstar

Part no.	Description	Unit
8011457	2.5 mm spec. (green)	Bag/100 pcs
8011456	4 mm spec. (blue)	Bag/100 pcs



New

Paper Kidney tray

Part no.	Description
8520139	1 pc




34 pcs

BAL

Disposable ear specula

For AirFx™.

Part no.	Description	Unit
8013013	2.75 mm	Bag/34 pcs
8013014	4.25 mm	Bag/34 pcs



50 pcs

BAL

Water irrigator tips

For use with **Interacoustics**® Aqua Stim™.

Part no.	Unit
8013018	Bag/50 pcs





**VOLUME
DISCOUNT
AVAILABLE**
for minimum
orders of
10 rolls

OAE

Thermal paper for ERO•SCAN™/ OtoRead™, 25 m

Width: 57 mm
Inner diameter: 12 mm
Outer diameter: 45 mm

Part no.	Unit
8500200	1 roll/25 m



**VOLUME
DISCOUNT
AVAILABLE**
for minimum
orders of
10 rolls

IMP

TPR10™ thermal paper for MTP10™, 25 m

Width: 79 mm
Inner diameter: 12 mm
Outer diameter: 45 mm

Part no.	Unit
8500220	1 roll/25 m



IMP

Thermal paper LPQ80, 80 m

For **GSI** TympStar™ Pro.
Width: 80 mm
Inner diameter: 13 mm
Outer diameter: 80 mm

Part no.	Unit
8509893	1 roll/80 m



**VOLUME
DISCOUNT
AVAILABLE**
for minimum
orders of
10 rolls

ABR IMP OAE

AP1300 thermal paper for Titan™, 10 m

Width: 57 mm
Inner diameter: 12 mm
Outer diameter: 35 mm

Part no.	Unit
8500230	1 roll/10 m



ABR IMP OAE

Sanibel™ MPT-III thermal paper, 20 m

For Sanibel™ MPT-III printer which is supplied
with **Interacoustics®** AT235™.

Width: 79 mm
Inner diameter: 12 mm
Outer diameter: 45 mm

Part no.	Unit
8500610	1 roll/20 m



IMP

Adhesive thermal paper LPQ80, 195 labels

For **GSI** TympStar™ Pro.
Width: 82 mm
Inner diameter: 25 mm
Outer diameter: 88 mm

Part no.	Unit
8509894	1 roll/195 labels



**VOLUME
DISCOUNT
AVAILABLE**
for minimum
orders of
10 rolls

IMP OAE

Thermal paper for Sanibel™ MPT-II printer, 20 m

For the Sanibel™ MPT-II printer which
is supplied with **MAICO®** ERO•SCAN™,
Grason-Stadler® Corti™ and
Interacoustics® OtoRead™.

Width: 57 mm
Inner diameter: 9 mm
Outer diameter: 40 mm

Part no.	Unit
8029305	1 roll/20 m



AUD IMP

TPR26™ thermal paper for Interacoustics®, 25 m

Fits most popular audiometers and
tympanometers.

Width: 112 mm
Inner diameter: 12 mm
Outer diameter: 45 mm

Part no.	Unit
8500210	1 roll/25 m



For OtoRead™ PRT MPT-II printer, 10 m

Part no.	Unit
8013009	1 roll/10 m



Heine/
KaWe
100 pcs

Welch
Allyn
250 pcs



No sharp
edges

AUD IMP

Otoscope specula

Disposable ear speculums for Heine/KaWe and Welch Allyn otoscopes size 2.5 and 4 mm.

Part no.	Description	Unit
8507618	4 mm Speculum Economy for Heine/KaWe	100 pcs
8507594	2.5 mm Speculum Economy for Heine/KaWe	100 pcs
8514367	4 mm Speculum for Heine/KaWe	100 pcs
8514368	2.5 mm Speculum for Heine/KaWe	100 pcs
8514369	4 mm Speculum Genuine for Heine/KaWe	50 pcs
8514370	2.5 mm Speculum Genuine for Heine/KaWe	50 pcs
8107750	2.5 mm Speculum for Welch Allyn	250 pcs
8108252	4 mm Speculum for Welch Allyn	250 pcs

2



New

AUD IMP

Luxamed Otoscope Auris LED 2.5V

Otoscope from German Luxamed. LED light to ensure longevity and great light source, dimable. Fits Heine/KaWe compatible speculums.

Part no.	Unit
8519433	1 pc

Luxamed Otoscope Auris LED 3.7V USB/PC Charged

Otoscope from German Luxamed. LED light to ensure longevity and great light source. USB charging. Fits Heine/KaWe compatible speculums.

Part no.	Unit
8519434	1 pc



No sharp
edges

50 pcs

AUD BAL IMP

Disposable speculas for video otoscope

Fits the AURICAL™ OTocam 300, **Interacoustics**® VIOT and XION™ video otoscopes.

Part no.	Unit
8500470	Bag/50 pcs

2



No sharp
edges

100 pcs

AUD BAL IMP

Video Otoscopes specula

Fits the new **Interacoustics**® Viot and OriVision otoscope.

Part no.	Description	Unit
8515083	KaWe 2.5 mm speculum	Bag/100 pcs
8511728	KaWe 4 mm speculum	Bag/100 pcs

2



AUD IMP

HEINE® Mini3000® Otoscope

Classic pocket otoscope from Heine in Germany.

Part no.	Unit
8108833	1 pc



HIF

Ear coupler sealing wax

Use to create a seal between coupler and hearing instrument.

Part no.	Unit
8500080	1 pc



ABR AUD OAE

Alkaline battery 1.5V AA

Part no.	Unit
8011226	1 pc




5 pcs

EARNimals®

Cartoon faces that are placed on the otoscope. Make ear examinations more fun for kids.

Part no.	Unit
8011501	Bag/5 pcs (assorted)



HIF

Battery adaptor

For **Interacoustics** Affinity™, Callisto™, Equinox™, **MedRx** Avant Hit+ and **Siemens** UNITY 2.

Part no.	Description	Unit
8011198	BAA675	1 pc
8011199	BAA13	1 pc
8011200	BAA312	1 pc
8011201	BAA10	1 pc
8011202	BAA5	1 pc



IMP HIF

Lithium battery 3.7V

For **Interacoustics**® Callisto™ and Titan™. Restrictions apply to shipping, contact us for more information.

Part no.	Unit
8009978	1 pc



New

Accessories for Otoscopy

Stylish wall dispenser for ear funnels. Can contain 6 different sizes.

Part no.	Description	Unit
8515678	Specula wall dispenser	1 pc



HIF

Hearing Instrument batteries Rayovac

Fits all types of Hearing Instruments.

Part no.	Description	Unit
8520134	Rayovac Extra Zap 10	6 pcs
8520135	Rayovac Extra Zap 13	6 pcs
8520136	Rayovac Extra Zap 312	6 pcs
8520137	Rayovac Extra Zap 675	6 pcs



HIF

Battery accessories

Fits all types of Hearing Instruments.

Part no.	Description	Unit
8515452	Battery tester	1 pc
8520138	Battery "carrier"	1 pc



New

Accessories for Otoscopy


Replacement bulb for Heine mini3000 Otoscope.

Part no.	Description	Unit
8515448	Replacement bulb	1 pc



Sanibel Supply

To order, please contact your local
special instrument distributor via
e-mail: info@diatec-diagnostics.nl



© Sanibel Supply A/S - 8506431 - Mar-57 - Rev. 7 - 1/2020

Sanibel Supply A/S • info@sanibelsupply.com • www.sanibelsupply.eu

Disclaimer

This catalogue has been published for worldwide distribution. Some parts described may NOT BE AVAILABLE in your area. Your local distributor will be pleased to give you details. For advice on compatibility and instructions of use, please consult your operation manual or local distributor. Sanibel Supply A/S assumes no legal liability or responsibility for the accuracy, typing errors or mistakes, correctness or completeness of any information in this flyer. The information in this flyer was correct to the best of our knowledge at the time of printing. Sanibel Supply A/S reserves the right to make changes at any time in specification of products described in this flyer without notice or obligation. MADSEN AccuScreen® and MADSEN OTOflex® are registered trademarks of Otometrics A/S. Echo-Screen®, AuDx®, Navigator® Pro, and ABAer® are registered trademarks of Natus Medical Incorporated. Sanibel™ is a registered trademark of Sanibel Supply®. All other Instrument brands or model names are copyrighted or trademarked by their respective owners and are used only for reference.

