


PRODUCT INFORMATION OTICON CHILI SP9, SP7 AND SP5

Oticon | Chili

Oticon Chili is a discreet super power instrument designed to meet the special and diverse needs of people with greater hearing losses, ranging from moderately-severe to profound. The slim durable 13 battery based BTE style offers prompt control and visual volume indication.

Powered by the RISE 2 platform, Chili applies Speech Guard signal processing to optimally preserve vital speech dynamics for the individual. Chili offers second generation wireless ConnectLine support for phone, mobile and TV use as well as full conventional FM/DAI and T-coil connectivity.



The power of understanding

Speech Guard is designed to prioritise conversations and facilitate selecting and following a speaker and to minimise listening effort. It maintains signal fidelity and uniquely preserves the richness of speech dynamics as well as responding instantaneously to sudden changing sound levels helping maintain focus on the speaker.

Spatial Noise Management eases listening in challenging noise environments. Chili uses binaural analysis of the listening environment to effectively identify situations where the listener can benefit from active Spatial Noise Management. In conversations with dominating noise from one side, the ear with the clearest speech signal is automatically prioritised to reduce listening effort.

Natural interaction

ConnectLine complements Oticon Chili in forming a unique and intuitive hearing solution providing the possibility for extended interaction with less stress. ConnectLine offers an intuitive wireless interface to modern entertainment and communication appliances like phone, mobile and TV.

Connect[+] adds Power Bass and Music Widening to music streaming to provide users a richer and more open listening experience.

Family Features

- Speech Guard
- Spatial Noise Management
- Binaural Dynamic Feedback Cancellation 2 (DFC2)
- Binaural Processing
- Binaural Synchronisation
- Binaural Push Button Coordination
- My Voice
- Multi Band Adaptive Directionality
- AutoPhone program
- Analogue volume control with mute
- Program button with mute
- Memory (datalogging)
- Streamer & ConnectLine options
- Power Bass
- Music Widening
- DAI and FM option
- T-coil

Fitting

- DSEsp, DSL v5.0a m[i/o], NAL RP and NAL-NL1 rationales
- Overall loudness trimmer
- In-situ audiometry
- nEARcom Cordless enabled



People First

People First is our promise to empower people to communicate freely, interact naturally and participate actively

www.oticon.com



PRODUCT OVERVIEW

Oticon | Chili

BTE SUPER POWER

OTICON CHILI SP9
OTICON CHILI SP7
OTICON CHILI SP5

Oticon | Chili

| MODEL FEATURES | CHILI SP9 | CHILI SP7 | CHILI SP5 |
|---|--------------------|--------------------|----------------------|
| Speech Guard | Yes | Yes | No |
| Spatial Noise Management | Yes | No | No |
| Binaural Dynamic Feedback Cancellation 2 (DFC2) | Yes | Yes | Yes |
| Binaural Synchronisation (DIR, NR) | Yes | No | No |
| Binaural Push Button Coordination | Yes | Yes | Yes |
| My Voice | Yes | No | No |
| Adaptive Directionality | Multiband/Tri mode | Multiband/Tri mode | Singleband/Dual mode |
| Noise Management | Tri State | Tri State | Modulation |
| Music Program | Yes | No | No |
| Bandwidth | 6.5 kHz | 6.5 kHz | 6.5 kHz |
| Streamer & ConnectLine options | Yes | Yes | Yes |
| Power Bass | Yes | No | No |
| Music Widening | Yes | No | No |
| DSEsp, DSL v5.0a m[i/o], NAL RP, NAL-NL1 | Yes | Yes | Yes |
| Identities for DSEsp | 4 | 4 | 1 |
| Fitting bands | 9 | 8 | 6 |

FITTING

Oticon Chili instruments are programmed using the Genie 2010.2 Fitting Software or higher compatible with NOAH 3 or higher. They can be programmed using either programming cables #3 and Programming Shoe or cordlessly using nEARcom.


Cordless fitting - nEARcom

nEARcom provides a cordless link between NOAHlink and one or two wireless enabled hearing instruments. In addition nEARcom provides a pass-through connection to accommodate programming cables and replaces the existing NOAHlink neck loop.

OPTIONS AND ACCESSORIES

| | |
|---------------------------------|---|
| Tamper resistant battery drawer | Available in all shell colours |
| Sound Hook | Interchangeable standard and paediatric hook <i>Note: use the dedicated black interchange tool</i> |
| Damper | Damping element option |
| DAI Adaptor | AP 900 |
| Dedicated FM Receiver | Amigo R12 |
| FM Adaptor | FM 9 Compatible with Amigo R1, R2 and other universal receivers |

COLOUR SELECTION





Scale 1:1

Technical Information

Omnidirectional mode is used unless otherwise stated.

Warning to the instrument dispenser

The maximum output capability of the hearing instrument may exceed 132 dB SPL (IEC 711). Special care should be exercised in selecting and fitting the instrument as there may be risk of impairing the remaining hearing of the hearing instrument user.



| | | | |
|----------------------------------|---------------|----------------|----------------|
| OSPL90 | Peak | 143 dB SPL | 139 dB SPL |
| | 1600 Hz | 131 dB SPL | 122 dB SPL |
| | Average | 136 dB SPL | 126 dB SPL |
| Full-on gain | Peak | 82 dB | 78 dB |
| | 1600 Hz | 73 dB | 65 dB |
| | Average | 76 dB | 67 dB |
| Frequency range | | 100-6500 Hz | 100-6200 Hz |
| Telecoil output (1600 Hz) | 1 mA/m field | 105 dB SPL | - |
| | 10 mA/m field | 125 dB SPL | - |
| | SPLITS L/R | - | 109/109 dB SPL |
| Total harmonic distortion | 500 Hz | 2.5 % | 2.5 % |
| (Input 70 dB SPL) | 800 Hz | 1.0 % | 1.0 % |
| | 1600 Hz | 2.0 % | 2.5 % |
| Equivalent input noise level (A) | Omni | 20 dB SPL | 21 dB SPL |
| | Dir | 31 dB SPL | 32 dB SPL |
| Battery consumption | Quiescent | 1.1 mA | 1.1 mA |
| | Typical | 1.2 mA | 1.8 mA |
| Battery life* | Calculated | 217 hours | |
| (Size 13, IEC PR48) | | | |
| IRIL (IEC 60118-13) | GSM/DECT | -49/-28 dB SPL | |

*) Based on the standardised battery consumption measurement. The actual battery life depends on battery quality, use pattern, active feature set, hearing loss and sound environment