EN - Instructions for use, pages 2 to 9.  
FR - Instructions d'utilisation, pages 10 à 17.  
DE - Bedienungsanweisungen, Seiten 18 bis 25.  
IT - Istruzioni per l'uso, pagine da 26 a 33.  
ES - Instrucciones de funcionamiento, páginas 34 a 41.  
SK - Návod na použitie, strany 42 až 49.  
PL - Instrukcja użytkowania, strony od 50 do 57.  
PT - Instruções de utilização, páginas 58 a 65.  
RO - Instrucțiuni de utilizare, paginile 66 până la 73.  
NL - Gebruiksaanwijzingen, pagina's 74 tot 81.  
TR - Kullanım talimatları, sayfa 82 ila 89.  
ΕΛ - Οδηγίες χρήσης, σελίδες 90 έως 97.  
RU - Инструкции по применению, стр. 98 - 105.  
112 إلى 106 تعليمات الاستخدام، الصفحات - AR

# CONTENTS

<b>1. DESCRIPTION OF THE DEVICE</b> .....	3
Model name.....	3
Commercial identification .....	3
Device description.....	3
<b>2. PACKAGE CONTENTS</b> .....	3
<b>3. LIST OF COMPATIBLE ACCESSORIES</b> .....	3
<b>4. SYMBOLS AND MEANINGS</b> .....	3
<b>5. INDICATIONS</b> .....	3
<b>6. MEDICAL CONTRAINDICATIONS</b> .....	4
<b>7. UNDESIRABLE SIDE EFFECTS</b> .....	4
<b>8. WARNINGS FOR USE</b> .....	4
<b>9. CAUTION, INTERACTION WITH TREATMENT AND CLINICAL INVESTIGATIONS</b> .....	5
Caution .....	5
High-voltage electrical field.....	5
Electrotherapy .....	5
Electrosurgery .....	5
Electroshock therapy and defibrillation.....	5
Medical diathermy .....	5
Diagnostic tests or treatments using ultrasound .....	5
Ionizing radiation therapy .....	5
Non-ionizing electromagnetic radiation .....	5
MRI (Magnetic Resonance Imaging) .....	5
<b>10. INSTRUCTIONS FOR USE</b> .....	6
Electrosurgery .....	6
Preoperative management .....	6
Intraoperative management .....	6
Determining the location of the incision and the position of the receiver .....	6
Handling the implant.....	6
Orientation of the implant .....	6
Fixation of the implant .....	7
Inserting the electrode-array.....	7
Positioning the reference electrode .....	7
Confirming the implant is functioning.....	7
<b>11. EXPLANTATION</b> .....	7
<b>12. TECHNICAL SPECIFICATIONS OF THE DIGISONIC® SP EVO IMPLANT</b> .....	8
<b>13. PERFORMANCE CHARACTERISTICS</b> .....	8
<b>14. SPECIFICATIONS AND CHARACTERISTICS OF THE ELECTRODE-ARRAY</b> .....	8

## 1. DESCRIPTION OF THE DEVICE

### • Model name

Digisonic® SP EVO cochlear implant

### • Commercial identification

I-SP-SD-EVO

### • Device description

The Digisonic® SP EVO is a multi-channel cochlear implant designed to rehabilitate patients with severe (2nd degree) to profound bilateral perceptive hearing loss. This cochlear implant works in conjunction with a behind-the-ear external processor (DIGI SP/ SAPHYR® SP/ZEBRA/DIGI SP'K). Electromagnetic induction between an external antenna and an internal one transmits the acoustic signal processed by the processor to the Digisonic® SP EVO receiver implanted under the skin just behind the auricle of the ear.

The Digisonic® SP EVO includes an electrode array containing 20 electrodes. This array is inserted into the implanted cochlea. Each individual electrode stimulates different groups of auditory nerve fibers and each corresponds to a different frequency band. The implant and the external antenna contain magnets that hold them in the proper relative position for optimal performance. The Digisonic® SP EVO does not contain a power source and only receives power when the external antenna is in place. If the external antenna is not present, the Digisonic® SP EVO is not powered and does not work.

## 2. PACKAGE CONTENTS

1 Digisonic® SP EVO (I-SP-SD-EVO) implant, 1 implant template made of silicon, 1 fixation system, 1 identification card and accompanying documentation.

## 3. LIST OF COMPATIBLE ACCESSORIES

The implant is supplied with the following accessories:

- 1 silicone implant template used during the surgery to verify the correct positioning of the implant under the skin.
- 2 self-tapping screws used to fix the implant into position; these do not require any pre-drilling.

## 4. SYMBOLS AND MEANINGS

	Catalog reference		Important. Read the accompanying documentation
	Sterilization method: ethylene oxide Sterile product		Maintain humidity level between 15% and 75%.
	Single-use device, do not reuse		Date of manufacture
	STORAGE condition: store in a cool dry place. Storage temperature between -20° C and +50° C.		Batch code
	Fragile; handle with care		Expiration date
	Do not use if the package is damaged		Serial number
	Do not re-sterilize		Manufacturer
	Waste electrical and electronic equipment (WEEE)		Instructions for use

## 5. INDICATIONS

The multi-channel Digisonic® SP EVO cochlear implant is designed to rehabilitate patients who have 2nd degree severe to total bilateral hearing loss with a speech intelligibility level of less than 50% at 60 dB HL in open-set format with the use of a hearing aid (professional medical opinion required).

## 6. MEDICAL CONTRAINDICATIONS

The Digisonic® SP EVO cochlear implant is not indicated in patients with perceptive hearing loss accompanied by large lesions in the cochlea (major cochlear malformation, fracture of the petrous pyramid, significant ossification of the cochlea), the auditory nerve (axonal neuropathy, tumor near or on the auditory nerve such as an acoustic neuroma, complete destruction of both auditory nerves), a severe anomaly of the auditory pathways, acute or chronic middle ear conditions (including tympanic membrane perforation), is psychologically unstable or has a contact allergy to implant materials (silicone, platinum iridium, titanium). Other types of implants may be recommended.

## 7. UNDESIRABLE SIDE EFFECTS

For patients who meet the indications for implantation must accept typical risks associated with surgery (effects from general anesthesia, infections, etc.) which is independent of any product. However, there is also a risk that the patient's body may reject the implant or a part of the implant. This risk has been reduced by using biocompatible materials in the design.

Complications associated with the cochlear implant surgery (temporary or permanent facial paralysis, risk of meningitis, changes in taste, dizziness, tinnitus, etc.) is rare, but should be considered carefully. It is important to inform every prospective implant candidate about these potential risks. Specific information should be given to the patient regarding the symptoms and the initial signs of meningitis. According to current recommendations, pneumococcal vaccination is also strongly recommended.

Once the implant is in place, there remain certain risks that may result in explantation. Explantation require additional surgical intervention under general anesthesia. Explantation may occur in the following cases:

- Medical complication
- Implant malfunction
- Displacement of the device as a result of trauma
- Extrusion of the implant

These postential problems were evaluated during product design and the materials and design of the implant have been chosen to minimize these risks.

Finally, the long-term effects from trauma associated with the insertion of electrodes and chronic electrical stimulation are unknown at the present time. These effects may include ossification of the cochlea or degeneration of the nerve fibers, and may require replacement of the implant or lead to a reduced response to stimulation.

## 8. WARNINGS FOR USE

### **Information to provide to the Patient:**

- The patient should be informed of the benefits of a cochlear implant, but also of its possible undesirable side effects (see §7).
- The supplied identification card must be completed.
- Inform the patient that they must present the identification card prior to any medical examination or treatment.
- Advise the patient to carefully read the user instructions supplied with his/her external processor, in particular the section relating to the warnings for use.
- In case of failure or malfunction of the cochlear implant system, the patient should contact his/her implantation center.
- Contact sports (rugby, boxing, etc.) are strongly discouraged, since strong impacts to the area can damage the implant.

- Scuba diving: Recreational scuba diving is not recommended at depths below 20m. Excess pressure can damage the implant. Moreover, it is strongly discouraged to engage in professional deep-sea diving activities. The implant is not guaranteed against repeated high pressure stresses.

- Access to restricted areas: Please consult a physician before entering restricted areas (MRI exam room, walk-through metal detectors, 3D scanning booths, etc.).

## 9. CAUTION, INTERACTION WITH TREATMENT AND CLINICAL INVESTIGATIONS

- **Caution:** Implantable device parts should not be reused if they have been previously implanted in another patient.

### - High-voltage electrical field:

- **Electrotherapy:** Electrotherapy may send currents of varying strengths. The use of high-voltage electrotherapy techniques is prohibited due to the risk of damage to the implant system. However, low-voltage electrotherapy may be considered only if the electrodes are not placed in areas near the head or neck.

- **Electrosurgery:** Avoid the use of monopolar electrosurgical instruments. These instruments can create radio-frequency fields through electrical charges that could cause the tip of the instrument and the electrode-array to interact. Generated currents might damage the cochlear tissue or result in permanent damage to the implant. As soon as a cochlear implant is removed from its packaging in the operating room, any monopolar surgical systems should be turned off to avoid any damage to the implant. **However, bipolar electrosurgical systems** can be used as long as it is not near to and does not come into direct contact with the cochlear implant.

- **Electroshock therapy and defibrillation:** Transmitting electrical shocks of several thousand volts through the body is not advised in a patient wearing a cochlear implant. Electrical shocks can cause tissue damage in the cochlea or permanently damage the implant.

In a life-threatening situation, only the medical team is in a position to make a decision. If the situation is not life-threatening for the patient, please contact Neurelec

- **Medical diathermy:** Electromagnetic diathermy cannot be used on patients with implants containing metal. This may cause irreparable damage to tissue (burning of the cochlea) or the implant. **However, ultrasonic diathermy may be considered** on areas of the body that do not touch the head or neck.

- **Diagnostic tests or treatments using ultrasound:** The implant should not be exposed to therapeutic levels of ultrasonic energy. The device may inadvertently concentrate the ultrasonic field and be damaged.

- **Ionizing radiation therapy:** During radiation therapy sessions, we strongly advise against direct irradiation of the implant zone. Direct, massive exposure of the implant to the rays could lead to partial or total loss of implant functions. This damage is not necessarily immediately noticeable.

### - Non-ionizing electromagnetic radiation:

- **MRI (Magnetic Resonance Imaging) :** The Digisonic® SP EVO implant contains a permanent magnet. An MRI exam or a strong magnetic field applied near the implant may cause damage to the implant and/or the patient. It is however possible to perform an MRI exam at 1.5 Tesla if recommendations are followed. In this case, the radiologist should fill out an exam form, which is available at [www.neurelec.com](http://www.neurelec.com).

Refer to the processor's user manual for more information.

## 10. INSTRUCTIONS FOR USE

### Please note!

**Electrosurgery:** Once the implant is put into place, the rest of the surgical procedure should not include the use of any monopolar electrosurgical instruments. These instruments can create radio-frequency fields by generating electrical charges that could damage the implant. **Bipolar electrosurgical instruments** can be used as long as they are not used near nor do come into contact with the cochlear implant.

### Preoperative management

- Before implanting the Digisonic® SP EVO for the first time, the surgeon should become familiar with the particular technical specifications of the Digisonic® SP EVO cochlear implant and the associated surgical technique.
- Before implantation, patients should be informed of the benefits of a cochlear implant, but also of its potential risks (see §7).


### Intraoperative management

**- Determining the location of the incision and the position of the receiver:** Before making the incision for the skin flap, it is recommended to determine the optimal site of the implanted system. In doing so, the location of the incision should account for the fact that the implant should be placed far enough from the auricle to allow space for the behind the ear processor without interfering with the external antenna. In addition, the incision line should be far enough away from the implant to prevent the risk of extrusion or postoperative infection. It is therefore recommended to use a sterile skin marker to draw the retroauricular incision line about 1cm away from the retroauricular crease, and to draw the receiver's position about 2cm away from the auricle of the ear. These locations can be determined by using the implant template and the processor template (supplied separately upon request) placed over the ear.

After making the incision, use the implant template for preparing the site so that the receiver is properly positioned. The surface of the bone should be checked to make sure it is flat enough to hold the receiver in place with screws.

**- Handling the implant:** Remove the implant from the sales packaging only after completing the standard surgical procedures up until the cochleostomy. Remove the outer packaging and make sure the implant is not damaged. Carefully cut the suture threads holding the implant in place and remove it from the sterile implant holder. Avoid touching and/or bending the electrode array. Do not use sharp surgical instruments that could damage the electrode-array.

**- Orientation of the implant:** The engraved side of the Digisonic® SP EVO implant (NEURELEC, BOTTOM) should be faced toward the skull and should not be visible. This metal titanium plate contains important information that identifies the implant as follows:

- **Fixation of the implant:** No milling of the bone is needed for the implant bed. Simply slide the receiver under the temporal muscle into the prepared site using the silicone template. Make sure the implant is positioned properly and does not move when finger pressure is applied. Only access to the implant's mounting

tabs is required for system fixation. The implanted receiver must be fixed with the supplied screws to prevent any possible displacement, which could create stress and possibly damage the electrode wires. It is always recommended to secure it in place with the two self-tapping screws provided in the packaging, and pre-positioned in the implant holder. Follow the steps below to remove the screws from the holder:

- Insert the screwdriver (provided with the initial surgical kit) into the screw using firm axial pressure.
- Unscrew the screw while slowly withdrawing it from the holder.
- The screw is now attached to the screwdriver and can be used.

Gently position the first screw into one of the fixation system's titanium inserts. It is recommended to hold the screwdriver vertical to the implant axis for fixation. Tighten the screw all the way, then gradually withdraw the screwdriver making a slight circular motion. Check that it is secure, then repeat the same procedure for the second screw.

Note: Neither the surgeon nor any person not authorized by Neurelec may make changes to the implant (such as to the fixation system).

- **Inserting the electrode-array:** Insert the electrode-array so that it will follow the cochlear spiral when inserted. Guide the tip of the electrode-array towards the base of the scala tympani using the insertion claw or micro-jeweler forceps. Then gradually advance the electrode-array using minimal force. Finish insertion by using the silicone extracochlear push rings as reference. Once insertion is complete, the rings should block the cochleostomy. The electrode array may be secured to prevent the risk of migration. The fixation method and fixation points will depend on the surgical access and the surgeon's preferences.

- **Positioning the reference electrode:** Place the head of the extracochlear reference electrode against the bone under the temporal muscle.

- **Confirming the implant is functioning:** Measuring impedance, before or after the incision is closed, can ensure that the implanted device is functioning properly.

## 11. EXPLANTATION

If malfunction of the Digisonic® SP EVO implant is suspected, the system must be examined with the help of Neurelec clinical support. If malfunction of the implant is confirmed and the medical team decides to explant, it is important to contact Neurelec, who will inform you of the explantation procedure to follow. It is especially important to request an explantation kit from Neurelec so that the explanted system may be sent back, thereby making it possible to examine the device.

<b>12. TECHNICAL SPECIFICATIONS OF THE DIGISONIC® SP EVO IMPLANT</b>			
Primary function	Cochlear implant		
Stimulation mode	Balanced biphasic stimulation		
Rate of stimulation	24,000 maximum pps (pps: pulse per second)		
Other available functions	<ul style="list-style-type: none"> <li>- Impedance measurement</li> <li>- Measurement of the implant's power</li> <li>- Integrity test*</li> <li>- EABR (Evoked Auditory Brainstem Response)*</li> <li>- Stapedius reflex</li> <li>- Psychoacoustic tests (gap test, etc.)*</li> </ul> * with associated equipment		
Weight	10.5 g		
Dimensions	Diameter 30.2 mm – Thickness: ranging from 4.9 (edge) to 5.75 mm (center)		
Volume	4 cm <sup>3</sup>		
Material in direct contact with human tissue	<ul style="list-style-type: none"> <li>- LSR 40 shore A silicone</li> <li>- HCR 35 shore A &amp; HCR 50 shore A</li> <li>- Adhesive silicone</li> <li>- Platinum iridium 10%</li> <li>- Titanium grade 2</li> <li>- Titanium grade 5</li> </ul>		
<b>13. PERFORMANCE CHARACTERISTICS</b>			
Characteristics of the output signal (on resistance of 1 kΩ)	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 60%;">Max 2V -255µs</td> <td style="width: 40%;">           U: 0 to 4V            I: 10 µA to 2 mA            Δt: 10 µs to 255 µs         </td> </tr> </table>	Max 2V -255µs	U: 0 to 4V I: 10 µA to 2 mA Δt: 10 µs to 255 µs
Max 2V -255µs	U: 0 to 4V I: 10 µA to 2 mA Δt: 10 µs to 255 µs		
Impedance measurement	Yes: Normal values: 500Ω – 5kΩ		
MRI safety level	Compatible with 1.5 Tesla. Refer to the recommendations		
Methods recommended for determining the proper functioning of the system	Yes. Impedance measurement and integrity test (with collection equipment)		
<b>14. SPECIFICATIONS AND CHARACTERISTICS OF THE ELECTRODE-ARRAY</b>			
Monopolar electrical configurations	Common ground		
Number of independent active electrodes	20		
Reduced cochleostomy size	Diameter of 0.8mm		
General shape	Straight with shape memory Straight : distance between electrodes and silicon, inferior to 0.1mm		
Shape at apex	Rounded shape		
Shape at the base	2 1.5mm-diameter push rings		
Material components	<ul style="list-style-type: none"> <li>- Connecting wire: Platinum iridium 10%</li> <li>- Stimulation electrodes: Platinum iridium 10%</li> </ul>		
Insulation	PE (polyester)		
Dimensions	<ul style="list-style-type: none"> <li>- Length; 25mm</li> <li>- Diameter at the base and at the apex: 0.5mm and 0.4mm</li> <li>- Minimum and maximum surface area of the stimulation electrodes: 0.46mm<sup>2</sup> (0.59mm<sup>2</sup> theoretical cylinder) et 0.60mm<sup>2</sup> (0.74mm<sup>2</sup> theoretical cylinder)</li> <li>- Distance between electrodes: 0.7mm</li> <li>- Maximum distance between proximal and distal electrodes: 22.3mm</li> </ul>		
Flexibility	- The insertion force is less than 32% above insertion of 200°, compared to the Digisonic® SP implant.		

- In case of doubt about the specified performance, please contact Neurelec customer service for your safety.
- Should you have any comments or if the information provided is incomplete, please contact your manufacturer or your local distributor.

**Digisonic®SP EVO is a registered trademark of Neurelec-France**

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