EN - Instructions for use, pages 3 to 9
FR - Instructions d'utilisation, pages 11 à 17
DE - Bedienungsanweisungen, Seiten 19 bis 25
IT - Istruzioni per l'uso, pagine da 27 a 33
ES - Instrucciones de funcionamiento, páginas 35 a 41
SK - Návod na použitie, strany 43 až 49
PL - Instrukcja użytkowania, strony od 51 do 57
PT - Instruções de utilização, páginas 59 a 65
RO - Instrucțiuni de utilizare, paginile 67 până la 73
NL - Gebruiksaanwijzing, pagina's 75 tot 81
TR - Kullanım talimatları, sayfa 83 ila 89
EL - Οδηγίες χρήσης, σελίδες 91 έως 97
RU - Инструкции по применению, стр. 99 - 105
113 إلى 108 الصفحات، الاستخدام، تعليمات - AR

CONTENTS

1. DESCRIPTION	4
Model Name	4
Commercial Identification	4
Device Description	4
2. PACKAGE CONTENTS.....	4
3. LIST OF COMPATIBLE ACCESSORIES.....	4
4. SYMBOLS AND MEANINGS.....	4
5. INDICATIONS.....	4
6. MEDICAL CONTRAINDICATIONS.....	5
7. ADVERSE SIDE EFFECTS.....	5
8. PRECAUTIONS.....	5
9. CAUTION, INTERACTION WITH CLINICAL TREATMENTS AND INVESTIGATIONS	6
Caution	6
High-Energy Electrical Field	6
Electrotherapy	6
Electrosurgery	6
Electroshock Therapy and Defibrillation	6
Medical Diathermy.....	6
Diagnostic Tests and Treatment with Ultrasound	6
Radiation Therapy with Ionizing Radiation	6
Non-Ionizing Electromagnetic Radiation	6
MRI (Magnetic Resonance Imaging)	6
10. INSTRUCTIONS FOR USE.....	7
Electrosurgery	7
Preoperative Management	7
Intraoperative Management	7
Placing the Incision and Positioning the Receiver	7
Handling the Implant.....	7
Orientating the Implant.....	7
Attaching the Implant.....	8
Inserting the Electrode-Array	8
Positioning the Reference Electrode.....	8
Confirming Functioning of the Implant.....	8
11. EXPLANTATION	9
12. TECHNICAL SPECIFICATIONS OF THE DIGISONIC® SP IMPLANT	9
13. PERFORMANCE CHARACTERISTICS	9
14. SPECIFICATIONS AND CHARACTERISTICS OF THE ELECTRODE-ARRAY	9

1. DEVICE DESCRIPTION

• Model Name

Digisonic® SP Cochlear Implant

• Commercial Identification

I-SP-SD

• Device Description

The Digisonic® SP is a multi-channel cochlear implant intended to rehabilitate bilateral severe (second degree) to profound perceptive deafness. This transcutaneous cochlear implant functions with a behind the ear external processor (DIGI SP/ SAPHYR® SP) or with a micro-behind the ear processor with remote battery (DIGI SP'K). By means of electromagnetic coupling, an external antenna transmits the acoustic signal processed by the processor to the Digisonic® SP implant implanted under the skin behind the auricle of the ear.

Le Digisonic® SP includes an electrode-array with 20 electrodes that is inserted into the cochlea (inner ear). Each electrode stimulates a different set of auditory nerve fibers and is associated with a frequency band of the sound signal that is processed by the external processor. The implant and the external antenna contain a magnet; the external antenna is thus kept in place relative to the implant by magnetic attraction. The Digisonic® SP does not contain its own electric power source; it receives the energy necessary for functioning by electromagnetic coupling. Consequently, when the external antenna is not correctly positioned relative to the Digisonic® SP, the implant is passive.

2. PACKAGE CONTENTS


1 Digisonic® SP (I-SP-SD) implant, 1 silicone implant template, 1 attachment system, 1 wearer-ID card, and accompanying documents.

3. LIST OF COMPATIBLE ACCESSORIES

The implant is supplied with the following accessories:

- 1 silicone implant gauge used during the surgery to verify the correct positioning of the implant under the skin.
- 2 self-tapping screws used to fix the implant in position; these do not require any pre-drilling.

4. SYMBOLS AND MEANINGS

	Catalog reference		Important. Read the accompanying documentation
	Sterilization method: ethylene oxide Product delivered sterile		Maintain humidity level between 15% and 75%.
	Single-use device, do not reuse		Manufacturing date
	STORAGE: store in a cool dry place. Storage temperature between -20° C and +50° C.		Batch code
	Fragile; handle with care		Use by date
	Do not use if package is damaged		Serial number
	Do not re-sterilize		Manufacturer
	Waste electrical and electronic equipment (WEEE)		Instructions for use

5. INDICATIONS

The multi-channel Digisonic® SP cochlear implant is designed to rehabilitate patients who have 2nd degree severe to total bilateral hearing loss with a speech intelligibility level of less than 50% at 60 dB HL in open-set format with the use of a hearing aid (professional medical opinion required).

6. MEDICAL CONTRAINDICATIONS

The Digisonic® SP cochlear implant is not indicated in patients with perceptive hearing loss accompanied by significant lesions of the cochlea (major cochlear malformation, fracture of the petrous part of the temporal bone, significant ossification in the cochlea), auditory nerve (axonal neuropathy, tumor near or on the auditory nerve such as a neurinoma, complete destruction of both auditory nerves), a severe anomaly of the auditory pathways, acute or chronic middle ear conditions (including tympanic membrane perforation), is psychologically unstable or has a contact allergy to implant materials (silicone, platinum iridium, titanium). Other types of implants may be recommended.

7. ADVERSE SIDE EFFECTS

For patients who meet the indications, implantation has typical risks associated with surgery (effects from general anesthesia, infections, etc.) which is independent of the product itself. However, there is also a risk that the patient's body may reject the implant or a part of the implant; this risk has been reduced by using biocompatible materials in the design.

Complications associated with the cochlear implant operating technique (temporary or permanent facial paralysis, risk of meningeal irritation syndrome, changes in taste, dizziness, tinnitus, etc.) are rare, but should be considered carefully. It is important to inform every cochlear implant candidate about these potential risks. Specific information should be given to the patient regarding the symptoms and the initial signs of meningeal irritation syndrome. According to current recommendations, pneumococcal vaccination is also strongly recommended.

Once the implant is in place, there remain certain risks that may result in explantation. Explantation requires another surgical intervention under general anesthesia. Explantation may occur in the following cases:

- Medical complication
- Implant malfunction
- Displacement of the device as a result of trauma
- Extrusion of the implant

These potential problems were evaluated during product design and the materials and design of the implant have been chosen to minimize these risks.

Finally, the long-term effects from trauma associated with the insertion of electrodes and chronic electrical stimulation are unknown at the present time. These effects may include cochlear ossification or degeneration of the nerve fibers, and may require replacement of the implant or lead to a reduced response to stimulation.

8. PRECAUTIONS

Information to give the Patient

- Inform the patient of the benefits of a cochlear implant, and also the adverse side effects which could occur (see §7).
- The supplied identification card must be fully completed and given to the wearer.
- Inform the patient that they must present the identification card prior to any medical examination or treatment.
- Advise the patient to carefully read the user instructions supplied with his/her external processor, in particular the section relating to the warnings for use.
- Tell the patient to contact the implantation center in case of failure or malfunction of the cochlear implant system.
- Contact sports (rugby, boxing, etc.) are strongly discouraged, since strong impacts to the area of the implant can damage it.

- Scuba diving: Recreational scuba diving is not recommended at depths below 20m. Excess pressure can damage the implant. Moreover, engaging in professional deep-sea diving activities is strongly discouraged, since the implant is not guaranteed against repeated high pressure stresses.

- Access to restricted areas: Please consult a physician before entering restricted access areas (MRI exam room, walk-through metal detectors, 3D scanning booths, etc.).

9. CAUTION, INTERACTION WITH CLINICAL TREATMENTS AND INVESTIGATIONS

- **Caution:** Implantable device parts should not be reused if they have been previously implanted in another patient.

- High-energy electrical field:

- **Electrotherapy:** Electrotherapy may send currents of varying strengths. The use of high-voltage electric-current electrotherapy techniques is prohibited due to the risk of damage to the implant system. However, low-voltage electrotherapy may be considered only if the electrodes are not placed in areas of the head or neck.

- **Electrosurgery:** Avoid the use of monopolar electrosurgical instruments. These instruments can create radio-frequency fields with electrical voltages that could create coupling between the tip of the instrument and the electrode-array. Induced currents could damage the cochlear tissue or produce permanent damage to the implant. As soon as a cochlear implant is removed from its packaging in the operating room, any monopolar scalpel should be turned off to avoid any damage to the implant. **However, bipolar electric scalpels** may be used as long as it is not near or in direct contact with the cochlear implant.

- **Electroshock therapy and defibrillation:** Sending several thousand volt electrical shocks through the body is not advised in a patient wearing a cochlear implant. Electrical shocks can cause tissue damage in the cochlea or permanently damage the implant.

In a life-threatening situation, only the medical team is in a position to make a decision. If the situation is not life-threatening for the patient, please contact Neurelec.

- **Medical diathermy:** Electromagnetic diathermy cannot be used on patients with implants containing metal. This may cause irreversible damage to tissue (burning of the cochlea) or the implant. **However, ultrasonic diathermy may be considered** on areas of the body that do not touch the head or neck.

- **Diagnostic tests and treatment with ultrasound:** The implant should not be exposed to therapeutic levels of ultrasonic energy. The device may inadvertently concentrate the ultrasound and be damaged.

- **Radiation therapy with ionizing radiation:** During radiation therapy sessions, we strongly advise against direct irradiation of the implant area. Direct, massive exposure of the implant to radiation could lead to partial or total loss of implant functions. This damage is not necessarily immediately noticeable.

- Non-ionizing electromagnetic radiation:

- **MRI (Magnetic Resonance Imaging):** The Digisonic® SP implant contains a permanent magnet. An MRI exam or the application of an intense magnetic field in the area of the implant could lead to damage to the implant and/or the patient. It is however possible to perform an MRI exam at 1.5 Tesla by following the recommendations. In this case, the radiologist should fill out an exam form available from the website www.neurelec.com.

Refer to the processor's user manual for more information.

10. INSTRUCTIONS FOR USE

Please note!

Electrosurgery: Once the implant is put into place, the rest of the surgical procedure should not include the use of any monopolar electrosurgical instruments. These instruments can create radio-frequency fields with electrical voltages that could damage the implant. **Bipolar electric scalpels** can be used as long as they are not used near or in contact with the cochlear implant.

Preoperative management

- Before performing the first Digisonic® SP implantation, the operator should become familiar with the particular technical specifications of the Digisonic® SP cochlear implant and the associated surgical technique.
- Before implantation, patients should be informed of the benefits of a cochlear implant and its possible secondary risks (see §7).


Intraoperative management

- Placing the incision and positioning the receiver: Before making the incision for the skin flap, determining the optimal site of the implanted system is recommended. In doing so, the location of the incision should account for the fact that the implant should not be placed below the auricle but at a large enough distance for the patient to wear the processor behind the ear without conflicting with the position of external antenna. In addition, the incision line should be far enough away from the implant to prevent the risk of postoperative extrusion or infection. Using a sterile skin pencil to draw the retroauricular incision line offset about 1cm from the retroauricular crease and the receiver's position about 2 cm away from the auricle of the ear is therefore recommended. These locations can be determined by using the implant template and the processor template (supplied separately upon request) placed in position over the ear.

After making the incision, use the implant template to prepare the location necessary for the proper positioning of the receiver. The surface of the bone should be checked to make sure it is flat enough to subsequently attach the receiver in place easily with screws.

- Handling the implant: Remove the implant from the sales packaging only after completing the standard surgical procedures up to starting the cochleostomy. Open the outer packaging making sure the implant is not damaged. After carefully cutting the suture threads holding the implant in place, gently remove the implant from its sterile holder. Avoid touching and/or bending the electrode array. Do not use sharp surgical instruments that could damage the electrode-array.

- Orientating the implant: The engraved surface of the Digisonic® SP implant (NEURELEC, BOTTOM) should not be visible in configuration for use and should therefore be turned toward the skull. This titanium metal-plate contains, as follows, important information with which to identify the implant:

- Attaching the implant: No milling of the bone is needed for the implant bed. Simply slide the receiver under the temporal muscle into the location prepared using the silicone template. Confirm that the implant is positioned properly and cannot rock with application of finger pressure. Only access to the implant's mounting tabs is required for attaching the system. The implanted receiver has to be stabilized to prevent any possible displacement, which could create tension and lead to possible damage to the electrode wires. It is therefore recommended to secure it in place with the two self-tapping screws provided in the packaging, and pre-positioned in the implant holder. Follow the steps below to remove the screws from the holder:

- Insert the screwdriver (provided with the initial surgical kit) into the screw using firm axial pressure.
- Unscrew the screw while slowly withdrawing it from the holder.
- The screw is now attached to the screwdriver and can be used.

Gently position the first screw into one of the fixation system's titanium inserts. It is recommended to hold the screwdriver vertical to the implant's axis to make attachment easier. Tighten the screw until feeling resistance, then progressively free the screwdriver by applying a slight circular motion to it. Check that it is secure, then repeat the same procedure for the second screw.

Note: Neither the surgeon nor any person not authorized by Neurelec may make changes to the implant (such as to the fixation system).

- Inserting the electrode-array: Orient the electrode-array so as to insert following the cochlear spiral. Guide the end of the electrode-array towards the base of the scala tympani using the insertion claw or small forceps, then progressively insert the electrode-array forcing as little as possible. Finish the insertion by pressing down on the silicone extracochlear push rings. Once insertion is complete, the rings should block the cochleostomy. The electrode array may be attached to prevent the risk of migration. The attachment method and attachment points will depend on the surgical access and the surgeon's preferences.

- Positioning the reference electrode: Place the extracochlear reference ball electrode against the bone under the temporal muscle.

- Confirming functioning of the implant: Impedance measurements, before or after closing the incision, can confirm the proper functioning of the implanted device.

11. EXPLANTATION

If malfunction of the Digisonic® SP implant is suspected, an expert review of the system must be done with the help of Neurelec clinical support. If malfunction of the implant is confirmed and the medical team decides to explant, it is important to contact Neurelec, who will inform you of the explantation procedure to follow. It is especially important to request an explantation kit from Neurelec so that the explanted system may be sent back, and thus the expert review of the device may be performed.

12. TECHNICAL SPECIFICATIONS OF THE DIGISONIC® SP IMPLANT	
Primary function	Cochlear implant
Stimulation mode	Balanced biphasic stimulation
Stimulation rate	24,000 pps maximum (pps: pulses per second)
Other available functions	- Impedance measurement - Measurement of the implant's power - Integrity test* - EABR (Evoked Auditory Brainstem Response)* - Stapedius reflex Psychoacoustic tests (gap test, etc.)* * With associated equipment
Weight	10.5 g
Dimensions	Diameter 30.2 mm – Thickness: ranging from 4.9 (edge) to 5.75 mm (center)
Volume	4 cm ³
Material in direct contact with human tissue	- LSR 40 shore A silicone - Silicone HCR 35 shore A and HCR 50 shore A - Silicone adhesive - Platinum iridium 10% - Titanium grade 2 - Titanium grade 5
13. PERFORMANCE CHARACTERISTICS	
Characteristics of the output signal (on resistance of 1 k Ω)	Max 2V -255 μ s V: 0 to 4V I: 10 μ A to 2 mA Δ t: 10 μ s to 255 μ s
Impedance measurement	Yes: Normal values: 500 Ω – 5k Ω
MRI safety level	Compatible with 1.5 Tesla. Refer to the recommendations
Recommended methods for determining the proper functioning of the system	Yes. Impedance measurement and integrity test (with collection equipment)
14. SPECIFICATIONS AND CHARACTERISTICS OF THE ELECTRODE-ARRAY	
Monopolar electrical configurations	Common ground
Number of independent active electrodes	20
Reduced cochleostomy size	1 mm diameter
General shape	Straight with shape memory
Shape at apex	Soft tip
Shape at the base	Two 1.5 mm diameter push rings
Constituent materials	- Connecting wire: Platinum iridium 10% - Stimulation electrodes: Platinum iridium 10%
Insulation	PEI (polyesterimide)
Dimensions	- Length: 26mm - Diameter at the base and at the apex: 1.07 mm and 0.5 mm - Minimum and maximum surface area of the stimulation electrodes: 0.39 mm ² and 0.77 mm ² - Distance between electrodes: 0.7 mm - Maximum distance between proximal and distal electrodes: 22.3 mm

- In case of doubt about the specified performance, please contact Neurelec customer service for your safety.
- Should you have any comments or if the information provided is incomplete, please contact your manufacturer or your local distributor.

Digisonic® SP is a registered trademark of Neurelec-France

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