



ampliSUITE

Intuitive audiometry
software application



For everyone



ampliSUITE

INTUITIVE AUDIOMETRY SOFTWARE APPLICATION

ampliSuite is designed for audiologists, ENT, hearing aid dispensers, private practitioners and other hearing health care personal.

Consisting of an audiometry and tympanometry module, you can read and review your patient's results with ease, compare thresholds, explain audiograms through built-in counselling tools and process test results quickly.



KEY FEATURES

- Audiometry and tympanometry modules in one suite
- Customisable PDF printouts with your own clinic/logo
- Simple and intuitive interface
- Large graphs for detailed review
- Pass and refer evaluation for each test result (tymp and reflexes)
- 17 available languages
- EMR connectivity including Noah, Auditbase, GDT and OtoAccess®
- Counselling tools for tone audiometrical test results, including adjustable PTA and baseline comparison.

ampliSUITE

AUDIOMETRY MODULE

The audiometry module has both pure tone audiometry and speech audiometry options which can be uploaded directly from the audiometer to ampliSuite. In order to help the patient understand their hearing situation (e.g. conductive, sensorineural or a combined hearing loss) and audiogram, counselling tools are available to be displayed within the audiogram in order to provide the best possible support and outcomes.

TYMPANOMETRY MODULE

The tympanometry module includes both tympanometry and acoustic reflex graphs, to make an effortless review of the measurements available. Normative boxes and Jerger classifications can also be added to the tympanometry graphs. Within an additional

CONNECTIVITY

ampliSuite is the perfect companion for the following devices:
Screening audiometers: Model 116, Model 170, PC850 and Otosure
Diagnostic audiometers: Model 240, Model 260, Model 270 and Model 270+
Tympanometers: all Otowave tympanometers can be used with ampliSuite

acoustic reflex sub-module, each reflex test can be selected and reviewed for middle-ear measurement diagnosis.

SOFTWARE INTEGRATION

ampliSuite flexibility stretches from stand-alone modules to the integration into third-party databases. ampliSuite is available as a free stand-alone software option as well as providing Noah and Auditbase integration, empowering hearing care personal all over the world to review audiological test data and support their patients in the best possible way.

Amplivox also offers the integration of ampliSuite into the patient database system OtoAccess® (version ampliSuite Pro), which is available to purchase.



The audiometry module has both pure tone audiometry and speech audiometry options



The tympanometry module includes both tympanometry and acoustic reflex graphs

START USING ampliSUITE TODAY

All three versions of ampliSuite are available to download from our website

DOWNLOAD NOW

CONTACT US

EMAIL US

Click on the relevant email address to message Amplivox directly.

Sales: sales@amplivox.com

General enquiries: hello@amplivox.com

UNITED KINGDOM

Tel: +44 (0) 1865 880 846

Amplivox Limited, 3800 Parkside, Solihull Parkway, Birmingham Business Park, Birmingham, West Midlands, B37 7YG, United Kingdom.

UNITED STATES

Tel: +1 888 941 4208

Fax: 952 903 4100

Amplivox, 10393 West 70th Street Eden Prairie, MN 55344, United States.

[Visit website](#)



VIDEO

Amplivox is trusted by many industry professionals globally, watch our video to learn more.



For everyone