

COMPARISON ► SCREENING AUDIOMETERS



COMPARE USER FEATURES	116	170	Otosure	PC850
Portable	◆	◆	◆	◆
PC controlled	-	-	◆	◆
Data transfer to PC	◆	◆	-	-
Thermal printer output	◆	◆	-	-
Manual test	◆	◆	◆	◆
Automatic test	-	◆	◆	◆
Self recording/Bekesy type test	-	-	-	◆
Automatic audiogram categorisation	-	-	◆	◆
Internal memory	◆	◆	-	◆
Talk forward	◆	◆	-	◆
Compatible with optional USB charger	◆	◆	-	◆
Carry case supplied	◆	◆	◆	◆
Patient response switch supplied	◆	◆	◆	◆
Mains power adaptor supplied	◆	◆	N/A	◆

COMPARE TECHNICAL FEATURES

	116	170	Otosure	PC850
Frequency range (Hz)	125-8K	125-8K	125-8K	250-8K
Attenuation range	-10 to 100 max	-10 to 100 max	-10 to 100 max	-10 to 100 max
Attenuation step size (dBHL)	5	5	5	5 / 2.5
Pure tone signal	◆	◆	◆	◆
Warble tone signal	◆	-	-	-
Continuous signal	◆	-	-	-
Pulsed signal	◆	-	-	-
Dimensions (L x W x D mm)	270 x 165 x 60	270 x 165 x 60	120 x 86 x 28	270 x 165 x 60
Weight (g)	720	710	165	710
External PSU	◆	◆	N/A	◆
Optional battery power	◆	-	-	-

For more information, please contact

e: sales@amplivox.com | t: +44 (0) 1865 880 846 | w: www.amplivox.com

Amplivox Ltd, 3800 Parkside, Solihull Parkway, Birmingham Business Park, Birmingham, West Midlands, B37 7YG, United Kingdom