

Spirometers

Technical features comparison



Features	Spirodoc	Spirolab	ML3500-Mk8	SpiroUSB
Portable	✓	✓	✓	✓
Automatic interpretation	✓	✓	✓	✓
Automatic quality control messages	✓	✓	✓	✓
On screen predictive values	✓	✓	✓	✓
Equivalent 'lung age' assessment	✓	✓	✓	✓
Internal database	✓	✓	✓	✓
Integral printer	-	✓	✓	-
Stand alone and PC operation (in real time)	✓	✓	✓	-
Touch screen display	✓	✓	✓	-
Disposable turbine/mouthpiece for maximum hygiene control	✓	✓	-	-
Disposable bacterial viral mouthpieces for hygiene control	✓	✓	✓	✓
Carry case supplied as standard	✓	✓	✓	✓
PC database software supplied as standard	✓	✓	-	✓
Weight	130g	1.6kg	630g	129g