

# COMPARISON ► SPIROMETERS



COMPARE FEATURES AND BENEFITS	Spirodoc	Spirolab	ML3500 Mk8 (microlab)	SpiroUSB
Portable	◆	◆	◆	◆
Automatic interpretation	◆	◆	◆	◆
Automatic quality control messages	◆	◆	◆	◆
On screen predictive values	◆	◆	◆	◆
Equivalent 'lung age' assessment	◆	◆	◆	◆
Internal database	◆	◆	◆	◆
Integral printer	-	◆	◆	-
Stand alone and PC operation (in real time)	◆	◆	◆	-
Touch screen display	◆	◆	◆	-
Disposable turbine/mouthpiece for maximum hygiene control	◆	◆	-	-
Disposable bacterial viral mouthpieces for hygiene control	◆	◆	◆	◆
Carry case supplied as standard	◆	◆	◆	◆
PC database software supplied as standard	◆	◆	-	◆
Weight	130g	1.6kg	630g	129g

For more information, please contact

e: sales@amplivox.com | t: +44 (0) 1865 880 846 | w: www.amplivox.com

Amplivox Ltd, 3800 Parkside, Solihull Parkway, Birmingham Business Park, Birmingham, West Midlands, B37 7YG, United Kingdom