

COMPARISON ► TYMPANOMETERS



COMPARE TECHNICAL FEATURES	102-1	102-4	102-C	202	202-H	302	302+
226Hz probe tone $\pm 0.1\%$	◆		◆		◆	◆	◆
1000Hz probe tone $\pm 0.1\%$	-		-	-	◆	-	◆
User selectable sweep speeds	-		-		◆	◆	◆
+200daPa to -400daPa pressure range ± 10 daPa	◆		◆		◆	◆	◆
0.2 to 5ml volumetric range ± 0.01 ml (226 Hz probe tone)	◆		◆		◆	◆	-
0.1 to 5ml volumetric range ± 0.01 ml (1kHz probe tone)	-		-	-	◆	-	◆
1kHz ipsi reflex function only	◆	-	-	-	-	-	-
500, 1k, 2k and 4k ipsi reflex function	-	◆	◆		◆	◆	◆
500, 1k 2k and 4k contralateral reflex function	-		-		◆	◆	◆
70 to 100dB reflex range	◆		◆		◆	◆	◆
0.01 to 0.5ml reflex measurement range	◆		◆		◆	◆	◆
Dimensions (W x H x D mm)	80 x 210 x 40		85 x 190 x 40		85 x 190 x 40	270 x 175 x 70	270 x 175 x 70
Weight (g)	380		380 (650 with cradle)		430	760	760

For more information, please contact

e: sales@amplivox.com | t: +44 (0) 1865 880 846 | w: www.amplivox.com

Amplivox Ltd, 3800 Parkside, Solihull Parkway, Birmingham Business Park, Birmingham, West Midlands, B37 7YG, United Kingdom

COMPARISON ► TYMPANOMETERS



COMPARE USER FEATURES	102-1	102-4	102-C	202	202-H	302	302+
Portable	◆		◆		◆	◆	◆
Six programmable testing profiles	-		-		-	-	◆
PC connectivity (USB)	-		◆		◆	◆	◆
PC connectivity (IrDA)	◆		-		-	-	-
Noah compatible	◆		◆		◆	◆	◆
Thermal printer output	◆		◆		◆	◆	◆
Internal memory	◆		◆		◆	◆	◆
Programmable 1 frequency (1kHz) ipsilateral reflex	◆	-	-		-	-	-
Programmable 4 frequency ipsilateral reflex function	-	◆	◆		◆	◆	◆
Programmable 4 frequency contralateral reflex function	-		-		◆	◆	◆
Cradle option	-		◆		-	-	-
Battery power function	◆		◆ (built-in)		◆	-	-
Mains power function (adaptor supplied)	-		◆ (via cradle)		◆	◆	◆
Compatible with optional USB Portable Charger	-		◆		◆	◆	◆

For more information, please contact

e: sales@amplivox.com | t: +44 (0) 1865 880 846 | w: www.amplivox.com

Amplivox Ltd, 3800 Parkside, Solihull Parkway, Birmingham Business Park, Birmingham, West Midlands, B37 7YG, United Kingdom