

# COMPARISON ► VISION TESTING



COMPARE FEATURES AND BENEFITS	Optec 5000P	Optec 5500P	Titmus V4	Titmus V2	Visiolite Master	Visiolite Master GT
Suitable for industrial vision screening	◆	◆	◆	◆	◆	◆
Dual lens system for accurate testing bifocal patients	◆	◆	-	-	◆	◆
6/4.5 to 6/20 acuity range	◆	◆	◆	◆	-	-
6/4.8 to 6/30	-	-	-	-	◆	◆
Meets recent DVLA 6/7.5 requirements	◆	◆	◆	◆	◆	◆
Remote control panel	-	◆	◆	-	Optional	Optional
Computerised Device with PDF output	-	-	-	-	◆	◆
Maintenance free lighting system	◆	◆	◆	◆	◆	◆
Portable	◆	◆	◆	◆	◆	◆
Far visual acuity	◆	◆	◆	◆	◆	◆
Near visual acuity	◆	◆	◆	◆	◆	◆

For more information, please contact

e: [sales@amplivox.com](mailto:sales@amplivox.com) | t: +44 (0) 1865 880 846 | w: [www.amplivox.com](http://www.amplivox.com)

Amplivox Ltd, 3800 Parkside, Solihull Parkway, Birmingham Business Park, Birmingham, West Midlands, B37 7YG, United Kingdom

# COMPARISON ► VISION TESTING



COMPARE FEATURES AND BENEFITS	Optec 5000P	Optec 5500P	Titmus V4	Titmus V2	Visiolite Master	Visiolite Master GT
Intermediate visual acuity	◆	◆	◆	◆	◆	◆
Depth perception test	◆	◆	◆	◆	◆	◆
Colour vision test	◆	◆	◆	◆	◆	◆
Lateral phoria test	◆	◆	◆	◆	◆	◆
Vertical phoria test	◆	◆	◆	◆	◆	◆
Peripheral vision test	◆	◆	◆	◆	◆	◆
Special tests (ARMD, Hyperopia, Astigmatism and Contrast)	-	-	-	-	◆	◆
Optional carry case	◆	◆	◆	◆	Included	Included
Optional intermediate lens	◆	◆	◆	◆	Included	Included
Weight	6.8 kg	6.8 kg	6 kg	6 kg	4.8 kg	4.8 kg
Dimensions (mm) width x length x height	280 x 450 x 490	280 x 450 x 490	280 x 420 x 190 (closed position)	280 x 420 x 190 (closed position)	470 x 430 x 350 (closed position)	470 x 430 x 350 (closed position)

For more information, please contact

e: sales@amplivox.com | t: +44 (0) 1865 880 846 | w: www.amplivox.com

Amplivox Ltd, 3800 Parkside, Solihull Parkway, Birmingham Business Park, Birmingham, West Midlands, B37 7YG, United Kingdom