

AMPLISUITE 1.2.3

We have released a new free to use version of ampliSuite to enhance functionality and correct minor bugs found in previous versions. The software will be included in all Amplivox audiometers and tympanometers manufactured from the 23rd April 2020 and can be downloaded from the Amplivox [website](#). Key updates include:

Run Otosure and PC850 from ampliSuite

It is now possible to start an audiometry test with the PC850 and Otosure audiometer directly from ampliSuite. This solution is ideal for those who wish to run screening programs and export results via printer or PDF into an external database.

Extra fields

Extra fields have been added to assist the operator through the test. New fields that have been added include the examiner details as well as a notes section which can be printed and are automatically integrated into the processed medical record.

Editable test date and results management

When test records are processed after the test date or where incorrect information has been extracted from the machine, information can easily be copied or edited.

South African Percentage Hearing Loss

South African Percentage Hearing Loss (PHL) calculation has been added to the software which can be turned on or off depending on requirement. Suitable for all South Africa occupational health screening professionals, especially those linking to third part Electronic Medical Record (EMR) systems or requiring printed results.

Connection Issues

Amplivox medical equipment could not always connect to computers where other fitting systems such as HiPro's were connected. This has now been fixed and means Amplivox medical equipment can run and transfer results alongside fitting systems such as the HiPro.

Printing Issues

Tympanometry results were not always displayed correctly depending on print settings. These have been rectified along with the addition of being able to print/PDF the speech banana and WHO hearing level tools.