

## TITMUS V2 AND V4

### TITMUS V SERIES MODEL COMPARISON

	TITMUS V4	TITMUS V2
Membrane control panel	●	
Built-in membrane switches		●
Ergonomic	●	●
New Fluorescent light source	●	●
Portability	●	●
Soft foam headrest	●	●
Head positioning sensor	●	
Peripheral vision testing	●	●
Test slide	●	●
Advance/reverse	Control panel	Manual
Night vision testing	●	
Standard accessories	●	●

### TECHNICAL SPECIFICATIONS

Power supply:	110-240 VAC, 0.4A-0.2A, 50/60 Hz
Cables:	Power cable, length 3m V4 Control panel cable; FireWire connector, length 1m
Illumination:	As per ISO 8596 and ANSI Z80.21 Standards

### PHYSICAL DATA

Dimensions (L x W x H):	280 x 430 x 190mm units (closed condition)
Weight:	6kg
Temperature:	+50°F to +104°F (+10°C to +40°C)
Humidity:	30% to 75% (non-condensing)
Pressure:	80-106 kPa

### TEST SLIDES AVAILABLE

Occupational:	Used by Occupational Health Professionals for screening vision of adults for various job functions
Professional:	Used by Family Physicians, General Practitioners for general screening of vision for both adults and children
Aeromedical:	Used by Aeromedical Examiners for screening vision of pilots
Paediatric:	Used by child healthcare professionals for general screening of vision for younger children
School:	Used by Family Physicians, Paediatricians for general screening of vision for all school children
Customised:	There is also an extensive range of test slides available should you require a customised solution for your vision screening requirements



### ELECTRICAL STANDARDS

Meets the following CE, UL, and CSA Standards.

**Europe:** EN 60601-1, IEC 60601-1, 3rd ed Class I, Annex IX of MDD 93/42/EEC, Class I Device

**CB Scheme:** IEC 60601-1

### EQUIPMENT

#### STANDARD

- 8 test slides made from high grade photographic film sealed between optical-quality glass
- Training manual
- Record forms – pad of 50 forms to record test results
- 10 lens cleaning wipes
- Package of 3 fog eliminator cloths
- Accessory case for storing lens cleaning wipes, fog eliminator cloths, and lenses
- Dust cover for instrument storage
- Carrying case.

For more information, please contact

e: sales@amplivox.com | t: +44 (0) 1865 880 846 | w: www.amplivox.com

The Amplivox policy is one of continuous development and consequently the equipment may vary in detail from the description and specification in this publication. Product specification may differ by country.