



Hygiene Consumables

Product catalogue
2022



For everyone



Hygiene for health screening

With the emergence of coronavirus (COVID-19) it's important for employers to take reasonable steps to protect their workers and others from infection. It's therefore critical to ensure that health screening programmes reduce and eliminate the potential for possible infection.

We offer a range of disposable hygiene solutions to make sure you can work hygienically and take care of your patient's health during coronavirus (COVID-19).

Why choose Amplivox?

We supply healthcare instruments, accessories and software for those who provide audiometry, spirometry and general health assessment programmes.

We support a wide variety of companies such as (and not limited to) the military, police forces, the NHS, private medical practices, manufacturers and the aviation industry. We are trusted by these industries to deliver the support they need to ensure they deliver the best service for everyone and we are here to provide you with reliable products to support you on a daily basis.

Audiometry

Audiometry: Infection control

Ear cushions can serve as potential carriers of micro-organisms (from the air, skin contact, draining ears, etc.) (Bankaitis, 2012).

Ear cushion covers serve as a physical barrier for headphone users where there is direct contact of the pinnae (outer ears) with the headset ear cushions.

Amplivox ear cushion covers

Made of acoustically transparent material, the covers are designed to fit over the earphones.

Ordering information:
Covers for standard headsets:
Amplivox stock code: 8107419
£63.63 per box (500)

Covers for Audiocups headsets:
Amplivox stock code: 8513404
£21.21 per box (100)





Audiometry

Ensure all equipment is cleaned between patients. Use of Clinell Universal wipes is recommended for:

- Audiometer front panel or pc keyboard
- Headset headband and cables
- Patient response switch (would also consider test subject wearing a disposable glove)
- Use appropriate ear cushion covers on the headset
- Consider testing outside the booth from a contamination perspective.

Spirometry

Ensure all equipment is cleaned between patients. Use of Clinell Universal wipes recommended for:

- Actual spirometer keys and screens
- Transducer housing
- Cables to transducer
- Consider test subject using a disposable glove to hold the transducer.

Amplivox would recommend use of PeraSafe@cleaning solution on a daily basis.

Ordering information:
PeraSafe@cleaning solution
Amplivox stock code: 1012009,
£27.03 per bottle (81g)

Clinell Universal wipes
Amplivox stock code: 1100000
£10.71 per pack (200)



Spirometry

We would like to advise the following if you are planning to complete spirometry sessions, with regard to safety critical medicals.

Amplivox recommend the use of single use bacterial viral filters (BVF) mouthpieces in place of the standard one-way filter or no-filter design mouthpieces. These BVF mouthpieces will provide:

- Improved hygiene
- Reassurance for the test subject undergoing the spirometry test
- The BVF mouthpieces will fit MicroMedical, MIR & Vitalograph spirometers.

Ordering information:
BVF mouthpieces for spirometers
Amplivox stock code: 8531510
£48.20 per box of 50
(individually sealed)

BVF mouthpieces for spirometers
Amplivox stock code: 1060011
£96.39 per box of 100
(individually sealed)

Spirometry: Additional sources of information

We would also recommend clients consider the environment in which spirometry is completed with regard to hygiene/cleaning and the patient's health.

Three helpful sources of information in this respect can be accessed by the links on the right:



Health and safety executive

Health surveillance and safety critical medicals during the coronavirus (COVID-19) pandemic

ARTP

ARTP COVID-19 update

British Thoracic Society

COVID-19: information for the respiratory community.





Blood pressure monitoring

Single use BP cuff covers

BP measurements are an essential part of medical and health surveillance testing. We would like to advise you of our new disposable single-use BP cuff covers that can provide safe and effective BP measurements.

The cuffs ensure a transmission barrier is created between the patient/employee when BP tests are to be completed.

- Can be used with any blood pressure cuff and offers additional protection for patients

- Does not affect blood pressure readings
- Comprises of a soft fabricel tissue fused to a leak-proof polyethylene film that offers 99.99% protection from bacteria pass-through
- Provide reduced cleaning requirements for cuffs
- Cuffs are latex-free.

Ordering information:

Single use BP Cuff Barriers
(Adult, 50qty, 46x17cm)

Amplivox stock code: 1011022 £16.07

Single use BP Cuff Barriers
(Adult Large, 50qty, 57x20cm)

Amplivox stock code: 1011023 £21.42

OMRON M6 Comfort Digital BP Monitor

Amplivox stock code: 1011020 £113.53

Blood pressure monitor

OMRON are recognised as a leading manufacturer of BP monitors. Amplivox are pleased to advise availability of the M6 monitor. A product that combines ease of use, performance and a competitive price.

- High performance, quality digital blood pressure monitor designed to withstand daily use within healthcare settings
- Omron's patented Intellisense inflation and deflation technology for patient comfort
- The clear LCD display shows heart rate and systolic and diastolic pressures
- Automatic indication of normal and high blood pressure.



General hygiene items

Nitrile examination gloves

8515551	Nitrile examination gloves small / box 100pcs	GBP	£19.40
8515552	Nitrile examination gloves medium / box 100pcs	GBP	£19.40
8515553	Nitrile examination gloves large / box 100pcs	GBP	£19.40
8515554	Nitrile examination gloves extra large / box 100pcs	GBP	£19.40

Disposable aprons

8519348	Disposable apron box 100pcs	GBP	£29.18
---------	-----------------------------	-----	--------

Disposable shoe covers

8519349	Disposable shoe covers box 100pcs	GBP	£18.24
---------	-----------------------------------	-----	--------

Disposable head covers

8519350	Disposable head covers box 100pcs	GBP	£21.28
---------	-----------------------------------	-----	--------

General hygiene items

General

8522417	Safety goggles 1pc	GBP	£6.74
1100003	Clinical 30L waste sacks - yellow 10x50	GBP	£27.25
1100007	Torch couch roll 48cm - blue 56m	GBP	£4.65
1100001	Clinell hand sanitiser 250ml pump action	GBP	£5.67
1100002	Clinell hand sanitiser 500ml pump action	GBP	£6.80



For everyone



Amplivox Ltd, 3800 Parkside, Solihull Parkway, Birmingham Business Park,
Birmingham, West Midlands, B37 7YG, United Kingdom

www.amplivox.com | +44 (0)1865 880846 | hello@amplivox.com

The Amplivox policy is one of continuous development and consequently the equipment may vary in detail from the description and specification in this publication.