



For everyone



Amplisuite

Audiological data management software

Amplisuite is an audiological software application designed to process test results that have been completed with Amplivox screening and diagnostic audiometry and tympanometry equipment.

Designed for use by ENT experts, hearing aid dispensers, private practitioners and other hearing health care personnel.

Audiometry module

The audiometry module has both pure tone audiometry and speech audiometry options where patient results can be uploaded from the audiometer's internal memory directly to Amplisuite.

To help the patient understand their hearing situation (e.g. conductive, sensorineural or a combined hearing loss) and audiogram, counselling tools are available to be displayed within the audiogram to provide the best possible support and outcomes.

Tympanometry module

The tympanometry module is designed to be used with our range of Otowave tympanometers to provide effortless data transfer of test records for an improved patient management workflow.

The tympanometry module includes both tympanometry and acoustic reflex graphs to review the measurements available.

Normative boxes and Jerger classifications can also be added to the tympanometry graphs. Within an additional acoustic reflex sub-module, each reflex test can be selected and reviewed for middle-ear measurement diagnosis.

EMR integration

Amplisuite has the flexibility to be used as a stand-alone product or part of a larger patient management record system.

Amplisuite can interface directly with NOAH and other NOAH compatible patient management systems such as Auditbase and Practice Navigator, to provide a streamlined and efficient audiology service.

Amplivox also offers the integration of Amplisuite into OtoAccess® (Amplisuite Pro), which is available to purchase.

Key features

- Audiometry and tympanometry modules in one suite
- Customisable PDF printouts with your own clinic/logo
- Simple and intuitive interface
- Large graphs for detailed review
- Pass and refer evaluation for each test result (tymp and reflexes)
- 17 available languages
- EMR connectivity including NOAH, Auditbase, GDT and OtoAccess®
- Counselling tools for tone audiometric test results, including adjustable PTA and baseline comparison

Visit www.amplivox.com/software/amplisuite to learn more about the features and benefits of Amplisuite

www.amplivox.com | +44 (0)1865 880846 | hello@amplivox.com

Connectivity

Screening audiometers

- Model 116
- Model 170
- PC850
- Otosure

Diagnostic audiometers

- Model 240
- Model 260
- Model 270
- Model 270+

Screening tympanometers

- Otowave 102-1
- Otowave 102-4
- Otowave 102-C

Diagnostic tympanometers

- Otowave 202
- Otowave 202-H
- Otowave 302
- Otowave 302+

Software updates

Amplisuite version 1.2.3 is now available to download. This latest version includes the following features:

- Otosure and PC850 compatibility
- Session comment section
- Option to edit and copy uploaded tests
- Option to add examiner name
- Percent hearing loss calculation used in South Africa



For everyone

Amplivox Ltd, 3800 Parkside, Solihull Parkway, Birmingham Business Park,
Birmingham, West Midlands, B37 7YG, United Kingdom

www.amplivox.com | +44 (0)1865 880846 | hello@amplivox.com

The Amplivox policy is one of continuous development and consequently the equipment may vary in detail from the description and specification in this publication.

EN-ASS40-VI