














Product portfolio overview















Alpha 9

 Hybrid Sound Processing™ Frequency bandwidth 10 kHz 24 fitting bands	
Frequency Composition™	
 Hybrid Balancing™ Speech Balancer (3 options)	
Hybrid Noise Management™ Smart Directionality Dynamic States (3 options) Omni States (2 options)	
Directionality controls <ul style="list-style-type: none">DynamicAdaptive Full DirectionalityFixed DirectionalityFixed Omni<ul style="list-style-type: none">Omni DirectionalTrue Directionality Plus	
Transition (4 options)	
Low Frequency Enhancer	
 Hybrid Balancing™ Noise Balancer (4 options)	
Hybrid Noise Management™ Smart Noise Reduction (4 options)	
Hybrid Feedback Canceller™	
Dynamic Range Extender	
Binaural Noise Manager	
Transient Noise Reduction (4 options)	
Wind Noise Manager	
Soft Noise Manager	
 Binaural coordination	
RC-A	
SoundClip-A	
EasyControl-A app	
iPhone Mic	
TV-A adapter	
Direct streaming from iOS and Android™*	
Made for all phones*	
 Personalization	
Tinnitus SoundSupport	
Adaptation Manager	
Comfort in Airplane	
Music Experience program	
Sound optimizer	












Alpha 7






 Hybrid Sound Processing™ Frequency bandwidth 8 kHz 20 fitting bands	
Frequency Composition™	
 Hybrid Balancing™ Speech Balancer (2 options)	
Hybrid Noise Management™ Smart Directionality Dynamic States (2 options) Omni States (2 options)	
Directionality controls <ul style="list-style-type: none">DynamicAdaptive Full DirectionalityFixed DirectionalityFixed Omni<ul style="list-style-type: none">Omni DirectionalTrue Directionality Plus	
Transition (3 options)	
Low Frequency Enhancer	
 Hybrid Balancing™ Noise Balancer (2 options)	
Hybrid Noise Management™ Smart Noise Reduction (4 options)	
Hybrid Feedback Canceller™	
Dynamic Range Extender	
Binaural Noise Manager	
Transient Noise Reduction (3 options)	
Wind Noise Manager	
Soft Noise Manager	
 Binaural coordination	
RC-A	
SoundClip-A	
EasyControl-A app	
iPhone Mic	
TV-A adapter	
Direct streaming from iOS and Android*	
Made for all phones*	
 Personalization	
Tinnitus SoundSupport	
Adaptation Manager	
Music Experience program	
Sound optimizer	








Alpha 5

 Hybrid Sound Processing™ Frequency bandwidth 8 kHz 18 fitting bands	
Frequency Composition™	
 Hybrid Noise Management™ Smart Directionality	
Directionality controls <ul style="list-style-type: none">DynamicAdaptive Full DirectionalityFixed DirectionalityFixed Omni<ul style="list-style-type: none">Omni Directional	
Transition (2 options)	
Low Frequency Enhancer	
 Hybrid Noise Management™ Smart Noise Reduction (3 options)	
Hybrid Feedback Canceller™	
Transient Noise Reduction (3 options)	
Wind Noise Manager	
Soft Noise Manager	
 Binaural coordination	
RC-A	
SoundClip-A	
EasyControl-A app	
iPhone Mic	
TV-A adapter	
Direct streaming from iOS and Android*	
Made for all phones*	
 Personalization	
Tinnitus SoundSupport	
Adaptation Manager	
Music Experience program	
Sound optimizer	

Viron 3

 Continuous Environment Detection	
ChannelFree™	
Speech Cue Priority™	
Frequency Composition™	
Frequency bandwidth 8 kHz	
10 fitting bands	
 Continuous Environment Detection	
Dynamic Noise Management™ (DNM™) Dynamic Directionality (Low focus)	
Directionality controls (Dynamic Directionality / Fixed Dir / Fixed Omni)	
Low Frequency Enhancer	
 Continuous Environment Detection	
Dynamic Noise Management™ (DNM™) Dynamic Noise Reduction	
Dynamic Feedback Canceller™ (DFC™)	
Transient Noise Reduction	
Wind Noise Manager	
Soft Noise Manager	
 Binaural coordination	
RC-A	
SoundClip-A	
EasyControl-A app	
iPhone Mic	
TV-A adapter	
Made for all phones*	
 Tinnitus SoundSupport	
Adaptation Manager	
Sound optimizer	

Viron 1

 Continuous Environment Detection	
ChannelFree™	
Speech Cue Priority™	
Frequency bandwidth 8 kHz	
8 fitting bands	
 Continuous Environment Detection	
Dynamic Noise Management™ (DNM™) Dynamic Directionality (Low focus)	
Directionality controls (Dynamic Directionality / Fixed Dir / Fixed Omni)	
Low Frequency Enhancer	
 Continuous Environment Detection	
Dynamic Noise Management™ (DNM™) Dynamic Noise Reduction	
Dynamic Feedback Canceller™ (DFC™)	
Wind Noise Manager	
Soft Noise Manager	
 Binaural coordination	
RC-A	
SoundClip-A	
EasyControl-A app	
iPhone Mic	
TV-A adapter	
Made for all phones*	
 Tinnitus SoundSupport	
Adaptation Manager	
Sound optimizer	

User benefits



Hybrid Technology™
Bernafon's Hybrid Technology™ is an industry first. It joins two processing approaches into one hybrid system. Prioritizing optimally in countless listening situations, it augments both speech understanding and listening comfort for best sound without compromise.



Sound
Strong signal processing provides users with better speech understanding in dynamic listening situations and a natural sound experience across all frequencies.



Speech
The right amount of directionality and speech boost offer an optimized binaural signal-to-noise ratio so that users can concentrate on what's important to them.



Hearing comfort
Different noise and feedback management systems reduce disturbing noise without adversely affecting speech providing users with more hearing comfort.



Operating comfort
Various accessories and connectivity options enable direct sound streaming into hearing aids for more control and operating comfort.



Individualization
Various individualization options are available for a more personalized fit.



*For information on compatibility, please visit www.bernafon.com/hearing-aid-users/hearing-aids/connectivity. Please note that availability of some features may be dependent on technical constraints of hearing aid style. Please check Product Information of desired style for more details.

Model overview



Technical overview	miniRITE	miniRITE T	miniRITE T R	BTE I05	IIC	CIC	ITC ¹⁾	ITE HS ¹⁾	ITE FS ¹⁾
Alpha 9 7 5		•	•						
Viron 9 7 5 3 1	•	•	•	•					
Zerena 9 7 5 3 1	•	•		•	•	•	•	•	•
Battery size	312	312	312+	13	10	10	312	312, 13	312, 13
Fitting level	Speaker 60, 85, 100, 105	Speaker 60, 85, 100, 105	Speaker 60, 85, 100, 105	Speaker 105	Speaker 75, 85	Speaker 75, 85	Speaker 75, 85, 90, 100	Speaker 75, 85, 90, 100	Speaker 75, 85, 90, 100
miniFit thin tube				•					
Earhook				•					
Li-ion rechargeable battery			•						
Dual-radio 2.4 GHz technology	•	•	•	•			○	○	○
NFMI technology	•	•	•	•		○	•	•	•
Single push button	•					○	○	○	○
Double push button		•	•	•					
Volume wheel							○	○	○
Telecoil		•	•	•			○	○	○
Autophone ²⁾	•	•		•			○	○	○
Hydrophobic coating	•	•	•	•	•	•	•	•	•
IP rating	IP68	IP68	IP68	IP68	IP68	IP68	IP68	IP68	IP68

Connectivity ³⁾ and accessories overview									
Made for iPhone®	•	•	•	•			•	•	•
Works with Android ⁴⁾		•	•						
Direct Audio Input (DAI) and FM adapter				•					
EasyControl-A app	•	•	•	•			•	•	•
RC-A (remote control)	•	•	•	•			•	•	•
TV-A (TV adapter)	•	•	•	•			•	•	•
SoundClip-A	•	•	•	•			•	•	•
Portable Charger Plus ⁴⁾			•						
Stable desk charger			•						
Wireless CROS/BiCROS ⁴⁾			•						

• Included ○ Optional
¹⁾ ITE hearing aids with battery size 312 can accommodate either 2.4 GHz wireless or telecoil. Choose battery size 13 if both options are required.
²⁾ Only available in Zerena hearing aids.
³⁾ For hearing aids with dual-radio 2.4 GHz technology.
⁴⁾ For Alpha hearing aids only.

BTE color options

Top shell

Base shell

midnight blue¹⁾

jet black²⁾

cocoa brown²⁾

metallic anthracite

metallic silver

sand beige

antique bronze

metallic anthracite

metallic silver²⁾

sand beige²⁾

midnight blue

metallic silver

metallic anthracite

¹⁾ Only available for Alpha 9 miniRITE T R
²⁾ Not available for Viron 3|1 and Zerena 3|1
Apple, the Apple logo, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

Accessories



SoundClip-A



RC-A (remote control)



TV-A (TV adapter)



EasyControl-A app

bernafon

Your hearing · Our passion

ITE color options

beige

light brown

medium brown

dark brown

black*

* For IIC only