

easyCart

The easyCart Facilitates Baby-Friendly, Room-to-Room Hearing Screening with your easyScreen ABR/OAE Device

Mobile Screening Station

The easyCart is our dedicated mobility solution for newborn hearing screening. It provides a convenient place to store the easyScreen, printer and additional accessories right where you need them.

Specifically developed for our easyScreen, the easyCart is the best solution to create an organized workplace ensuring that you can focus on the needs of your patients instead of the device itself.

Organized Setup

The easyScreen and the printer can be attached to the easyCart top shelf, ensuring secure connection during transport. Other accessories and disposables can be stored in the dedicated compartments which are also securely mounted on the easyCart.

The easyCart is engineered to perfectly fit the MAICO easyScreen, providing a practical tool to help you perform reliable hearing screening in any location.

Keep the Device Ready-to-Go

With your power strip mounted safely to the cart, you can always charge your devices easily with one plug-in. Just be ready for that next hearing screening without much fuss.

Optional Tablet Accessory Kit

A custom gooseneck pole holds a tablet PC in a locked enclosure. Thus, you may keep everything on one compact cart. You do not need to move the system to a different location, in order to transfer your screening data. Just complete the entire screening process including data transfer at crib-side.

Features at a Glance

- Helps to get organized, every part has its place
- Structured setup enables an easy workflow
- Top Shelf to attach device and accessories
- Compartment for disposables
- Device and cradle fixation to make sure they remain in place on the go
- 5 large wheels create a secure and mobile cart to keep pace in a busy hospital or clinic
- Dedicated additional parts available like tablet holder



Technical Data easyCart

Compartments

Utility bin front (H x W x D)	12.7 x 30.5 x 30.5 cm
Chart bin back (H x W x D)	8.3 x 36.2 x 22.9 cm
Fixation	Clamps for quick release

Top Shelf

Dimension (W x D)	35.5 x 28 cm
Fixations with screw for:	easyScreen cradle device in cradle optional label printer
Cut outs for	power supply cables preamplifier insert phones OAE probe

Handle to push the cart

Wheelbase

Width	53.3 cm
Wheels	5 (2 with locks)
Wheel diameter	7.6 cm
Wheel material	Rubber



Total Dimensions

Height / Weight	115.5 cm / 11.4 kg
-----------------	--------------------

Standard Components

Wheelbase, stand, power strip mount, two bins with clamps, top shelf, cradle and easyScreen fixations, printer mounting adapter

Optional Accessories

- BERAPhone holder
- Gooseneck pole with VESA mounting plate
- Tablet enclosure for Surface Pro 4
- No power strip included

Specifications are subject to change without notice



The devices and disposables illustrated are not included in the scope of delivery of the easyCart.

Distributed by:

MAICO Diagnostics GmbH

Sickingenstr. 70-71 · 10553 Berlin · Germany

Tel.: +49 30 / 70 71 46-50 · Fax: +49 30 / 70 71 46-99

sales@maico.biz · www.maico.biz