


	Oticon Opn 1	Oticon Opn 2	Oticon Opn 3
Taleforståelse			
OpenSound Navigator™	Niveau 1	Niveau 2	Niveau 3
- Effekt af balancering	100%	50%	50%
- Maks. fjernelse af støj	9 dB	5 dB	3 dB
Speech Guard™ LX	Niveau 1	Niveau 2	Niveau 3
Spatial Sound™ LX	4 estimatorer	2 estimatorer	2 estimatorer
Soft Speech Booster LX	•	•	•
Speech Rescue™ LX	•	•	•
Lydkvalitet			
Clear Dynamics	•	•	-
Spatial Noise Management	•	•	-
Frekvensbåndbredde	10 kHz	8 kHz	8 kHz
Processeringskanaler	64	48	48
Bass Boost (streaming)	•	•	•
Lyttekomfort			
Transient Noise Management	4 konfigurationer	Til/Fra	Til/Fra
Feedback shield LX	•	•	•
Wind Noise Management	•	•	•
Individualisering & optimeret tilpasning			
YouMatic™ LX	3 konfigurationer	2 konfigurationer	1 konfiguration
Tilpasningskontroller	16	14	12
Lytteprogrammer	•	•	•
Direktionalitetsmuligheder	•	•	•
Tilvænningsstrin	•	•	•
Rationaler	VAC+, NAL-NL1+2, DSL v5.0	VAC+, NAL-NL1+2, DSL v5.0	VAC+, NAL-NL1+2, DSL v5.0
Konnektivitet			
Stereo-streaming (2.4 GHz)	•	•	•
Made for iPhone®	•	•	•
Oticon ON App	•	•	•
ConnectClip	•	•	•
Remote Control 3.0	•	•	•
TV Adapter 3.0	•	•	•
Særlige behov			
Tinnitus SoundSupport™	•	•	•

OTICON | Opn

Produktoversigt



TwinLink
NFMI + 2.4 GHz

oticon
PEOPLE FIRST