

**Acto**  
**All together better**  
for everyone



Introducing **Oticon Acto**, the new high performance mid-market hearing instrument. Available in the miniRITE style and built on the RISE 2 platform, Oticon Acto empowers people with hearing loss to **communicate freely, interact naturally and participate actively** in greater comfort and convenience. This is what our People First promise is all about.

Oticon Acto is the first mid-market Oticon hearing instrument to offer **binaural synchronisation** as just one best in class feature in a full technology package. We believe that in terms of sound quality, connectivity and flexibility, Oticon Acto is the **all together better** option for your clients and your business.





All together better  
**comfort, clarity and  
understanding**

Oticon | **Acto**



All together better  
**communication  
and  
entertainment**



All together better  
**to work with**



Available in  
**miniRITE**



**The all together better  
mid-market hearing solution**

Oticon Acto features technologies and benefits previously only associated with premium hearing instruments. It delivers fully on Oticon's commitment to providing better sound quality while also raising the bar in terms of flexibility and convenience.



Oticon Acto's exceptional performance is driven by the RISE 2 processing platform. RISE 2 is designed for superior sound and processing power.

In addition, Oticon Acto's embedded wireless connectivity capabilities allows the user to switch seamlessly between entertainment and communication devices, all intuitively managed with just one remote control.

All together better  
**comfort, clarity and understanding**

**Oticon Acto presents the ultimate opportunity to enjoy life in greater detail and comfort. Your clients can truly wear it and forget it.**

With this hearing instrument Oticon has striven to create the full package for people with mild to severe hearing loss.

Oticon Acto integrates a variety of best in class sound technologies such as Multi-band Adaptive Directionality, TriState Noise Management and DecisionMaker2. The enhanced processing power and 8 kHz bandwidth of the RISE 2 platform also contribute to a fuller, richer and more balanced sound picture.

Meanwhile the open fit RITE (Receiver-in-the-Ear) design provides a richer frequency response and a more comfortable mix of natural and amplified sound.





## All together better **communication and entertainment**

ConnectLine restores the full benefits of modern communication and entertainment technology to people with hearing loss.

Oticon Acto and ConnectLine are an ideal option for people who wish to enjoy communication and entertainment without impediment or inconvenience. They have full control over telephone calls, music and television at the touch of a button.

ConnectLine channels sound directly from the TV, telephone or music player to the hearing instrument. Combining ease of use and significantly improved audibility, ConnectLine enables the user to stay in touch with friends and family.


### TV

The TV sound goes directly into the user's ears at the touch of a button.



### Landline phone

Direct hands-free connectivity for the landline phone. Phone calls can even be answered remotely at the touch of a button.



### MP3

The hearing instrument can be used as a wireless headset for MP3 players, PC etc. (easily connected via Bluetooth or cable).

### Mobile phone

The hearing instrument can be used as a wireless hands-free solution for mobile phones with excellent clarity and audibility.

## ConnectLine

enables users to enjoy all the entertainment and convenience offered by modern media and communication technology. They can, for example, listen to music or watch TV at a volume that also suits the other people in the room.

When the telephone rings, ConnectLine discreetly notifies the user who can answer the call directly with the touch of a button. When the call ends, the TV sound returns.



All together better  
to work with


With Oticon Acto you can offer your clients a mid-market product giving access to state of the art technology and features - with a style to match everyone.

Oticon Acto is available in a variety of styles that allow you to cater to clients' individual needs and cosmetic preferences.

**miniRITE** is exceptionally small and lightweight and fits comfortably behind the ear. Despite its size, the organically shaped miniRITE offers the industry's most advanced wireless capacities and fits up to 110 dB HL with the flexible RITE speaker system offering 3 power levels.

For you, Oticon Acto represents a far more patient-centric and flexible solution. For example, it enables you to use in-situ audiometry for flexible threshold verification and fitting. Oticon Acto also features configurable controls in Genie and individual sound source characteristics adjustments.

To ensure maximum client satisfaction, Oticon also supplies our unique eCAPS counselling software tool. It is just one aspect of our ongoing business support for dispensers and our ongoing commitment to our People First promise.





## All together better coverage of mid-market needs

### All together better comfort, clarity and understanding

- RISE 2 signal processing platform
- Advanced binaural automatics
- 8 kHz bandwidth
- TriState Noise Management
- Multiband Adaptive Directionality
- Styles to match any personal preference

### All together better

### All together better communication and entertainment

- One Streamer that gives access to all communication; TV, phone, mobile etc.
- Easy and intuitive user interface developed with real users
- Small and lightweight

### All together better to work with

- Broad range of styles, including the stylish miniRITE design
- RITE speaker system with 3 power levels allows for flexibility and future proof fitting
- Flexible in-situ audiometry
- Wireless fitting option



miniRITE and RITE are available with three power levels for even greater flexibility. Just change the speaker if more power is needed.

Built on Oticon's  
most advanced  
processing  
platform, RISE 2



# Meeting all individual needs

## miniRITE/RITE




miniRITE




RITE

### STANDARD / MEDIUM



Open Dome



Plus Dome



Micro Mould

### POWER




Power Dome



Power Micro Mould



Power Mould



STANDARD

MEDIUM

POWER

■ Custom Mould - Power Dome

□ Open Dome - Plus Dome

## CUSTOM



ITE\*



ITC Power (Omni or Dir)\*\*



ITC\*




CIC/MIC Power\*\*



CIC/MIC\*\*

\* Also available as a non-wireless style  
\*\* Not wireless



■ ITC Power

□ CIC/MIC


## BTE 13



Hook



Corda² thin tube




■ Earmould

□ Corda²

## BTE POWER



Hook



### Main Features

	Acto Pro	Acto
Binaural Synchronisation	•	•
Binaural Coordination	•	•
Binaural Dynamic Feedback Cancellation 2 (DFC2)	•	•
DecisionMaker2™	•	•
Multiband Adaptive Directionality	•	•
TriState Noise Management	•	•



### Fitting options (Genie 2010.2)

	Acto Pro	Acto
nEARcom enabled	•	•
In-situ audiometry	•	•
NAL-NL1 & DSL v5.0a m[i/o]	•	•
Identities	3	3
Fitting Bands	8	6
Memory	•	•

### More

	Acto Pro	Acto
Bandwidth	8 khz	8 khz
3 RITE receivers (S, M and P)	•	•
RISE 2 platform	•	•


### Connectivity

	Acto Pro	Acto
ConnectLine enabled	•	•
DAI and FM	•	•




44

Silver




91

Silver Grey




92

Steel Grey




63

Diamond Black




93

Chestnut Brown



94

Terracotta



90

Chroma Beige

## People First

People First is our promise  
to empower people  
to communicate freely,  
interact naturally and  
participate actively