

Introducing the most powerful, sophisticated and smallest Super Power ever





**Designed
with you in mind**

Oticon Chili is developed in collaboration with people who know the everyday challenges in communicating and participating. Not just people with a hearing loss, but also their families and friends, to find the best solution for whatever the day brings.

That is how Oticon Chili has become the next generation of super power instruments introducing groundbreaking new approaches in audiology, connectivity and reliability. It is much more than super-powered. It's a new way forward in meeting the particular needs of those with severe-to-profound hearing loss. It's designed specifically to empower you, every waking moment of the day, every day of the week.

It's small and cosmetic and the modern design even incorporates fast control features desired for instant muting, volume adjustment and program selection.

Never before has super power been so understated!



The concepts that make the difference

Speech Guard - clear undistorted speech

Oticon Chili introduces Speech Guard, a unique amplification technique that helps you achieve maximal audibility as well as speech recognition without the jumble and distortion you could encounter in even the most advanced hearing aids today. Speech Guard has been well documented to provide a clear and natural signal reducing the clatter and increasing the ability to follow conversations, even in more complex environments where noise or more voices are competing for your attention.

Effortless access to phones, TV and mobiles

Oticon ConnectLine offers unmatched ability to stay in touch, keep informed and enjoy entertainment. ConnectLine lets you watch TV at the preferred loudness level with sound directly in the hearing instruments while the rest of the family enjoys the TV sound comfortably from the TV-speakers. If the phone rings you will be discreetly notified by a tone in your ears, and you can pick up the call just by pushing one button without getting out of the couch. When the call ends, the system automatically switches back to the TV sound.

And most importantly everything is easily controlled by you with one small and very accessible Streamer control unit.



Build to last - quality to rely on

Chili is designed with double protection - inside out - to protect the instrument from moisture, sweat, dirt and shocks. For example Chili is:

- Shock absorbent. The unique shock-absorbing fixation of the most important electronic components prevent them from breaking should the instrument be dropped or fall off the ear.
- Resistant to water. Moisture and dirt are kept out by the high-tech water repellent coating and the battery door is designed to reject water and sweat that might trickle down the instrument.
- Deep coated. Inside a permanent coat effectively seals the electronic parts resulting in a fully water, moisture and dirt repellent hearing instrument.

Chili is not only ideal for those who live active lives but simply for all super power users who demand 24/7 reliability!



Inconspicuous in every way

Have you ever asked yourself what it would take to make you change from your current hearing instruments to try something new? We fully acknowledge that it appears to be a big step, and we understand any hesitation.

On the other hand, we sincerely believe that now is the time to do so. The progress achieved in sound quality and speech clarity is very tangible and we have experienced both surprise and delight from users who have already tried Oticon Chili. Surprise because they very fast and with no or very few fine tunings reached a better result than what they were used to. Delight because the clarity and connectivity offered by Chili reach a level never before obtained.



People First

People First is our promise to empower people to communicate freely, interact naturally and participate actively