

Oticon CROS compatibility overview



	Oticon CROS miniRITE T	Oticon CROS PX miniRITE R
Oticon More 1* : miniRITE R, miniRITE T	•	•
Oticon More 2* : miniRITE R, miniRITE T	•	•
Oticon More 3* : miniRITE R, miniRITE T	•	•
Oticon Opn S 1** : miniRITE, miniRITE T, miniRITE R, BTE PP	•	•
Oticon Opn S 2** : miniRITE, miniRITE T, miniRITE R, BTE PP	•	•
Oticon Opn S 3** : miniRITE, miniRITE T, miniRITE R, BTE PP	-	-
Oticon Xceed 1** : BTE SP, BTE UP	•	•
Oticon Xceed 2** : BTE SP, BTE UP	•	•
Oticon Xceed 3** : BTE SP, BTE UP	•	•
Oticon Ruby 1*** : miniRITE, miniRITE T, miniRITE R, BTE, BTE PP	•	•
Oticon Ruby 2*** : miniRITE, miniRITE T, miniRITE R, BTE, BTE PP	•	•
Oticon Opn Play 1** : miniRITE, miniRITE T, miniRITE R, BTE PP	•	•
Oticon Opn Play 2** : miniRITE, miniRITE T, miniRITE R, BTE PP	-	-
Oticon Xceed Play 1** : BTE SP, BTE UP	•	•
Oticon Xceed Play 2** : BTE SP, BTE UP	•	•

* Firmware 1.1.0 and onwards
 ** Firmware 8.0 and onwards
 *** Firmware 1.0 and onwards