

Introducing EduMic

Stay active and engaged



The perfect partner
for your Oticon
hearing aids

oticon
life-changing technology

The perfect partner throughout your day



A hearing aid accessory with many functions

EduMic helps when you need direct and clear sound, by transmitting audio or a speaker's voice into your hearing aids.

- ✓ **Remote microphone** - Transmits a speaker's voice directly to your hearing aids
- ✓ **Telecoil receiver** - Captures audio from public hearing loop systems and sends it to your hearing aids
- ✓ **Media streamer** - Sends audio from your computer, tablet, TV and more, directly to your hearing aids

Exceptional sound quality

EduMic gives you the cutting-edge technology known from Oticon hearing aids. The combination of EduMic and Oticon hearing aids* gives you a clear signal with exceptional sound quality.

Very easy to use

With just a single push-button, EduMic is simple and easy to operate. It only needs to be paired with your Oticon hearing aids* once via Bluetooth® wireless technology, and there's no bulky receiver to attach to your hearing aids.

EduMic easily clips onto clothing or is worn on a lanyard, where it looks sleek and discreet.



* EduMic is compatible with Oticon 2.4 GHz Bluetooth low energy technology hearing aids: Oticon Opn™, Oticon Opn S™, Oticon Xceed, Oticon Opn Play™, Oticon Xceed Play, Oticon Ruby, and Oticon Siya.



Hear the person speaking with incredible clarity

Transmit a speaker's voice directly to overcome distance and noise

EduMic gives you an extra ear that's closer to what you want to hear. It is ideal when a person is speaking from far away, or if their voice is struggling to cut through the noise. EduMic clips onto their clothes, to send their voice to your hearing aids directly and clearly.

For learning, EduMic helps you when teachers are hard to follow from a distance. With one press of the button you'll be hearing their voice better and getting much more from the lesson. You'll be focusing on what counts: learning, and keeping up with the class.

EduMic can send speech to many Oticon hearing aids at once - so your partner, friends, or relatives can benefit from exceptional sound quality as well.

For meetings at work and elsewhere, EduMic can help in noisy environments, so you can participate actively.

For sports on the playing field or in the gym, EduMic lets you hear the trainer more easily, so you can focus more closely on what you're doing.



Direct, wireless audio for those special moments



Get clear sound from hearing loop systems - without a telecoil in your hearing aids

EduMic lets you connect to hearing loop systems in public buildings such as theatres, churches, museums, lecture halls, or concert venues - even if you don't have telecoil hearing aids.

EduMic has a built-in telecoil, so you can connect to hearing loop systems anywhere you see the symbol.

Then you can enjoy the beauty of clear sound even in noisy public venues.

EduMic is also compatible with other remote microphone systems via its standard Europin interface.

Stream quality audio directly from your devices

EduMic plugs into devices with a standard 3.5 mm audio jack, so you can stream audio from your computer, tablet, TV, music player, and more.

Once plugged in, EduMic sends the sound directly to your hearing aids.

Control the sound with our mobile app

The Oticon ON App lets you adjust the sound transmitted from EduMic to your hearing aids using your smartphone.



Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

Contact your local hearing care professional:

oticon.global/edu-mic

oticon
life-changing technology