


	Oticon Opn 1	Oticon Opn 2	Oticon Opn 3
Speech Understanding			
OpenSound Navigator™	Niveau 1	Niveau 2	Niveau 3
- Balancing power effect	100%	50%	50%
- Max. noise removal	9 dB	5 dB	3 dB
Speech Guard™ LX	Niveau 1	Niveau 2	Niveau 3
Spatial Sound™ LX	4 estimators	2 estimators	2 estimators
Soft Speech Booster LX	•	•	•
Speech Rescue™ LX	•	•	•
Sound Quality			
Clear Dynamics	•	•	-
Spatial Noise Management	•	•	-
Aanpas bandbreedte	10 kHz	8 kHz	8 kHz
Verwerkingskanalen	64	48	48
Bass Boost (streaming)	•	•	•
Luistercomfort			
Transient Noise Management	4 configuraties	Aan/Uit	Aan/Uit
Feedback shield LX	•	•	•
Windruismanagement	•	•	•
Personalisatie & Optimaliseren aanpassing			
YouMatic™ LX	3 configuraties	2 configuraties	1 configuratie
Aanpasbanden	16	14	12
Luitsterprogramma's	•	•	•
Meerdere directionaliteitsopties	•	•	•
Adaptatiemanagement	•	•	•
Aanpasregels	VAC+, NAL-NL1+2, DSL v5.0	VAC+, NAL-NL1+2, DSL v5.0	VAC+, NAL-NL1+2, DSL v5.0
Connecting to the World			
Stereo streaming (2.4 GHz)	•	•	•
Made for iPhone®	•	•	•
Oticon ON App	•	•	•
ConnectClip	•	•	•
Afstandsbediening 3.0	•	•	•
TV Adapter 3.0	•	•	•
Speciale behoeften			
Tinnitus SoundSupport™	•	•	•

OTICON | Opn

Product- vergelijking



TwinLink
NFMI + 2.4 GHz

oticon
PEOPLE FIRST