

Release Note for Users

# HearSIM™ 3.3.5


Title:	HearSIM™ 3.3.5 Release Note for Users
Project:	HearSIM™ 3.3.5
Date issued:	October 4 <sup>th</sup> 2023

## Table of Contents

Introduction .....	2
Installation .....	2
Windows® requirements .....	2
New HearSIM™ features .....	3

## Introduction

We're pleased to announce that a new version of the HearSIM™ PC Software is now available. This software can be used with several newborn hearing screening devices.

If you already have the HearSIM™ application on your PC, the software version number can be checked at the bottom of the main menu  accessible in the HearSIM™ start screen. This document provides a summary of features related to the software changes from **HearSIM™ version 3.3.2 to the new version named above.**

The new HearSIM™ software and instructions for use will be available as of the effective date **(only available for authorized service partners).**

If you upgrade HearSIM™, it is necessary to check the version of the device firmware and upgrade to the latest version (Version 1.24.21).

## Installation

This version can be installed as a new installation or as an update. Just run the setup file from the installation package and follow the installation guide. For details about the installation, consult the Instructions for Use documents for OtoAccess® Database and for HearSIM™:

There are some preconditions that have to be fulfilled prior to installation of HearSIM™:

- OtoAccess® Database must be installed first. OtoAccess® Database software and instructions for use are available for download at [otoaccess.com](http://otoaccess.com).
- OtoAccess® Database requires entry of a license code for each PC installation.

The HearSIM 3.3.5 is available for download at the download section of the OtoAccess webpage:

<https://www.otoaccess.com/downloads>

### Windows® requirements

	HearSIM™
Current released version	3.3.5
Windows® compatibility:	
- 10 Professional – 32 bit	Yes
- 10 Professional – 64 bit	Yes
- 11 Professional – 32 bit	Yes
- 11 Professional – 64 bit	Yes

## New HearSIM™ features

- Additional server upload option:
  - Direct upload option of tracking export files to a WebDAV server for easier tracking file transfer to the tracking center.
- Support of multiple instances of SQL databases on an OtoAccess Server with different names
- Bug fixes
  - Minor bug fixes

Please find further details in the HearSIM Instruction for Use, which is available under Manuals at <https://www.otoaccess.com/downloads>