

Bone Anchored Hearing System


Ponto 4

Open Sound. Open Life



Made for
 iPhone | iPad | iPod

Always read the instructions for use. Follow the instructions for use.



oticon
MEDICAL

Join the open sound revolution and a future of connected hearing



Velox S™ platform
– speed and power




OpenSound Navigator™
– 360 degree sound

A paradigm shift in bone anchored hearing

Hearing aid users around the world have already had their lives changed by the groundbreaking technology in Oticon Opn™ and Oticon Opn S™. For the first time, this is now available to bone anchored users. This unique technology gives constant access to 360-degree sound and lets your brain decide what to focus on. The result is a full and open soundscape.

Always read the instructions for use. Follow the instructions for use.



The world's smallest bone anchored sound processor*
– 27% smaller than Ponto 3

2.4 GHz wireless

First internet connected bone anchored sound processor**

Unbelievably small

Ponto 4 is the smallest bone anchored device in the market* for discreet, reliable performance.

And fully wireless

Ponto 4 can be connected to all your favourite compatible electronic devices; and with the Oticon ON App and the IFTTT network, it can also connect to any internet enabled device, such as home appliances, email programmes and more***.

**Data on file at Oticon Medical*

***Using Oticon ON App and IFTTT network*

****IFTTT compatible devices required*




Experience open sound with Ponto 4


Introducing the revolutionary OpenSound Navigator™

With Ponto 4, traditional sound processing strategies is a thing of the past. The unique, groundbreaking chip in Ponto 4 processes sound 50 times faster* and gives constant access to 360-degree sound.

In this more complete and open soundscape, you are in charge of the sounds you want to focus on – even in noisy environments.


**Compared to the Inium Sense platform*






OpenSound Navigator

– 3 steps in a split second




1: Analyse

Scans the 360° sound environment more than 100 times per second to identify and separate noise from speech.



2: Balance

Rapidly reduces the levels of loud noise coming from specific directions, while preserving speech.



3: Noise removal

Rapidly attenuates remaining diffuse noise, even between individual words.

*Always read the instructions for use.
Follow the instructions for use.*

The world's first truly connected bone anchored sound processor



Truly connected

Ponto 4 is more than a sound processor – it is a piece of wearable technology that can connect you with your favourite compatible devices.

It uses 2.4 GHz technology and is also the world's first internet connected bone anchored sound processor.*

However, there is more to it. The chip in Ponto 4 is prepared for the rapid changes in technology in hearing healthcare so you can continue to stay connected in the future.

**Using the Oticon ON App and the IFTTT network.*



Connect to IFTTT – the Internet of Things

Through a unique Oticon cloud solution, Ponto 4 and the Oticon ON App can connect to any other internet enabled devices and services**.

For example:

- Get a voice alert when someone rings your doorbell
- Get a text message when your battery is low
- Switch to home programme when entering the front door

Visit oticon.global to explore the world of opportunities.

***IFTTT compatible devices required.*

Always read the instructions for use. Follow the instructions for use.

Wireless connection to your compatible phone and favourite devices



Made for iPhone

Ponto 4 connects directly to your iPhone® so you can talk on the phone and hear the conversation in your sound processor. You can also stream music directly from your compatible phone without needing an extra device.

Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

Control your Ponto 4 with the Oticon ON App

The Oticon ON App lets you discreetly adjust the volume, switch between programmes, use “Find my hearing aids”, access the Internet of Things and so much more.

Download the free app on the App Store® and on Google Play™.



For information on compatibility, please visit:
www.oticonmedical.com/wireless-compatibility

Wireless accessories

Use your Ponto 4 in even more situations:

- **ConnectClip** – Turn your Ponto 4 into a hands-free mobile headset with any compatible smartphone
- **TV Adapter 3.0** – Listen to compatible TV or music through your Ponto 4
- **Remote Control 3.0** – Adjust the volume and change programmes




ConnectClip



Remote Control 3.0

TV Adapter 3.0



Additional accessories are available for purchase and not included as part of the Ponto Bone Anchored Hearing System package.

Always read the instructions for use. Follow the instructions for use.

The world's smallest bone anchored sound processor*

Tried, tested and more than ready for daily life

The Ponto 4 delivers the renowned standards of the Ponto quality in the smallest size to date. It has been thoroughly tested to meet high standards of durability and reliability, including the demanding IP 57 rating for water and dust resistance.



1:1 scaling



Featuring:

- LED notification light
- Tamper resistant battery drawer
- Robust metal spring coupling
- Skins and stickers for personalisation



Chroma Beige
C090



Terracotta
C094



Chestnut Brown
C093



To find out what Ponto 4 can do for you,
talk to your hearing care professional today or go to:
www.oticonmedical.com



Steel Grey
C092



Diamond Black
C063



Silver
C044

Additional accessories are available for purchase and not included as part of the Ponto Bone Anchored Hearing System package.

Always read the instructions for use. Follow the instructions for use.

Ponto 4 personalised for you


Stickers and skins

Choose a colour that matches your hair or make your personality come through by personalising your Ponto 4 with a range of skins and stickers.



Skin colours

1:1 scaling



Stickers



Additional accessories are available for purchase and not included as part of the Ponto Bone Anchored Hearing System package.

Because sound matters

Oticon Medical is a global company in implantable hearing solutions, dedicated to bringing the magical world of sound to people at every stage of life. As part of the Demant group, a global leader in hearing healthcare with 14,500 people in over 130 countries, we have access to one of the world's strongest research and development teams, the latest technological advances and insights into hearing care.

Our competencies span more than a century of innovations in sound processing and decades of pioneering experience in hearing implant technology. We work collaboratively with patients, physicians and hearing care professionals to ensure that every solution we create is designed with users' needs in mind. We share an unwavering commitment to provide innovative solutions and support that enhance quality of life for people wherever life may take them. Because we know how much sound matters.



 **Oticon Medical AB**
Datavägen 37B
SE-436 32 Askim
Sweden
Tel: +46 31 748 61 00

Always read the instructions for use. Follow the instructions for use.

oticon
MEDICAL | Because
sound matters

www.oticonmedical.com