


# Ponto comparison chart

Get an overview of different Ponto models and see the benefits of upgrading to the new Ponto 5 Mini

Product information	Ponto 3	Ponto 4	Ponto 5 Mini	Why choose Ponto 5 Mini?
Fitting range	Up to 45 dB HL	Up to 45 dB HL	Up to 45 dB HL	
Size (L x W x H)	34 x 21 x 11 mm	26 x 19 x 11 mm <i>World's smallest</i>	26 x 19 x 11 mm <i>Still world's smallest</i>	
Weight (without battery)	14 g	13.2 g	13.2 g	
Hardware certification (including battery compartment)	IP57	IP57	IP57	
Battery type	13	312	312	
Frequency range	200 Hz – 9.5 kHz	200 Hz – 9.5 kHz	200 Hz – 9.5 kHz	
LED indicator light	-	✓	✓	
Performance				
BrainHearing technology	✓	✓	✓	All Ponto models support the brain to make sense of sound for better speech understanding with less listening effort. <sup>1</sup>
Open sound experience	-	OpenSound Navigator <i>Enables a 360° soundscape</i>	OpenSound Navigator <i>Enables a 360° soundscape</i>	For a truly open and balanced sound experience that gives users access to 360° sound – proven to improve speech understanding by more than 20%. <sup>1</sup>
	-	-	OpenSound Optimizer <i>Detects and prevents audible feedback from occurring</i>	For 6 dB more stable gain – more sound, not less, all day. <sup>2*</sup>
Noise reduction	TriState Noise reduction	OpenSound Navigator <i>Precise analysis of the sound environment for ultra-fast noise removal</i>	OpenSound Navigator <i>Precise analysis of the sound environment for ultra-fast noise removal</i>	Allows users to move seamlessly between different environments, providing better speech understanding. <sup>1</sup>
Feedback management	Inium Sense feedback shield	Feedback shield LX	Feedback shield LX	Feedback shield LX removes feedback with traditional feedback handling technologies.
	-	-	OpenSound Optimizer <i>Detects and prevents audible feedback from occurring</i>	Prevents audible feedback without reducing gain – up to 6 dB more gain without risk of feedback. <sup>2*</sup>
Directionality	FreeFocus <i>Traditional directionality</i>	-	-	Ponto 5 Mini uses OpenSound Navigator instead of traditional directionality.
Technology platform	Inium Sense <i>15 frequency channels</i>	Velox S™ <i>50 times faster data processing than Inium Sense. Outstanding resolution with 64 frequency channels</i>	Velox S™ <i>50 times faster data processing than Inium Sense. Outstanding resolution with 64 frequency channels. New acoustic measurements 56,000 per sec. enabling OpenSound Optimizer</i>	Velox S is Oticon Medical's fastest and most advanced technology platform with the speed and precision to follow rapidly changing conversations with multiple speakers, even in noise.
RemoteCare	-	-	✓	Enables remote appointments, adjustments and upgrades.
Special purpose programme	-	-	3 special purpose programmes <i>Speech in noise / Comfort / Music</i>	Different programmes for different listening situations.
Connectivity and streaming				
Stereo streaming	Via Oticon Medical Streamer	2.4 GHz Bluetooth Low Energy	2.4 GHz Bluetooth Low Energy	Connects with the Oticon ON App and ConnectClip and lets users control the sound processor, take hands-free calls, stream music and more.
Made for iPhone®	-	✓	✓	Direct connection to your iPhone, so you can stream phone calls, stereo music and video sounds directly to your sound processor.
ConnectClip	-	✓	✓	Allows direct streaming from any modern smartphone, connects to laptops via the USB dongle and can be used as a remote microphone.
EduMic	-	✓	✓	Streams conversation up to 20 metres and gives access to FM, telecoil or jack modes.
Remote Control App	✓ <i>Via Oticon Medical Streamer</i>	✓ <i>Via Oticon ON App or Remote Control 3.0</i>	✓ <i>Via Oticon ON App or Remote Control 3.0</i>	Discreet control of volume, programmes and connectivity sources.
TV Adapter	✓ <i>TV Adapter 2.0</i>	✓ <i>TV Adapter 3.0</i>	✓ <i>TV Adapter 3.0</i>	Streams stereo sound directly from the TV.
Apps				
Control app	ConnectLine App	Oticon ON App	Oticon ON App	Features in Oticon ON App only: Discreet volume and program control. Includes a "find my Ponto" location feature and access to Hearing fitness, IFTT and the Hearing Guide.
Support and guidance app	Ponto Care App	Ponto Care App	Ponto Care App	Supports users from trial to aftercare to enable a good experience.
Accessories				
Softband and SoundConnector	✓	✓	✓	Provides access to Ponto for those too young or unable to undergo surgery. Lets users try Ponto before they decide to have the implant.
Safety Line	✓	✓	✓	Prevents loss of sound processor and protects against damage.
Tamper-resistant battery drawer	✓	✓	✓	Prevents accidental opening of the battery door.
Skins and stickers	✓	✓	✓	Enables personalisation.
Fitting				
Wireless fitting	-	✓	✓	Fit Ponto with Noahlink Wireless.
Remote fitting	-	-	✓ <i>Oticon RemoteCare App</i>	Users can attend appointments online and have programming and adjustments carried out remotely.
NAL-NL 1 BC	✓	✓	✓	Choose the most relevant fitting rationale for the user.
DSL BC	✓	-	✓	

<sup>1</sup> Manuscript in preparation, Data on File – Clinical study BC102

<sup>2</sup> Data on file at Oticon Medical

\* When fitted according to recommendation

Apple, the Apple logo, iPhone®, iPad®, and iPod touch® are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.