

The upgrade you’ve been waiting for

Ponto 5 SuperPower is packed with the latest open sound benefits. For an experience beyond compare.

Product information		Ponto 3 SuperPower	Ponto 5 SuperPower	Why upgrade to Ponto 5 SuperPower?
Fitting range		Up to 65 dB HL BC	Up to 65 dB HL BC	
Size (L x W x H)		34 x 21 x 14 mm	31 x 21 x 15 mm <i>World's smallest SuperPower device</i>	
Weight (without battery)		17 g	20 g	
Hardware certification (including battery compartment)		IP57	IP57	
Battery type		675P	675P	
Frequency range		260 Hz – 9.5 kHz	200 Hz – 9.5 kHz	
LED indicator light		–	✓	
Push button		Volume wheel and push button	Double push button	
Performance				
BrainHearing™ technology		✓	✓	All Ponto models support the brain to make sense of sound for significantly better speech understanding.* ¹
Open sound experience		–	OpenSound Navigator™ <i>Enables a 360° soundscape</i>	For a truly open and balanced sound experience that gives users access to 360° sound – proven to significantly improve speech understanding.* ¹
		–	OpenSound Optimizer™ <i>Detects and prevents audible feedback from occurring**²</i>	For 5 dB more stable gain – more sound, not less, all day.** ²
Noise reduction		TriState Noise reduction	OpenSound Navigator <i>Precise analysis of the sound environment for ultra-fast noise removal</i>	Allows users to move seamlessly between different environments, providing significantly better speech understanding.* ¹
Feedback management		Inium Sense feedback shield	Feedback shield LX	Feedback shield LX removes feedback with traditional feedback handling technologies.
		–	OpenSound Optimizer <i>Detects and prevents audible feedback from occurring**²</i>	Prevents audible feedback without reducing gain – up to 5 dB more stable gain without the risk of feedback.** ²
Directionality		FreeFocus <i>Traditional directionality</i>	–	Ponto 5 SuperPower uses OpenSound Navigator instead of traditional directionality.
Technology platform		Inium Sense <i>15 frequency channels</i>	Velox S™ <i>50 times faster data processing than Inium Sense. Outstanding resolution with 64 frequency channels. New acoustic measurements 56,000 per sec. enabling OpenSound Optimizer</i>	Velox S is Oticon Medical's fastest and most advanced technology platform with the speed and precision to follow rapidly changing conversations with multiple speakers, even in noise.
Special purpose programme		–	3 special purpose programmes <i>Speech in noise / Comfort / Music</i>	Different programmes for different listening situations.
Connectivity and streaming				
Stereo streaming		Via Oticon Medical Streamer	2.4 GHz Bluetooth® Low Energy	Connects with the Oticon ON app and ConnectClip and lets users control the sound processor, take hands-free calls, stream music and more.
Made for iPhone®		–	✓	Direct connection to iPhone, so users can stream phone calls, music and video sounds directly to the user's sound processor.
ConnectClip		–	✓	Allows direct streaming from any compatible smartphones, connects to laptops via the USB dongle and can be used as a remote microphone.
EduMic		–	✓	Streams conversation up to 20 metres and gives access to FM, telecoil or jack modes.
Remote Control app		✓ <i>Oticon Medical Streamer and ConnectLine app</i>	✓ <i>Via Oticon ON app or Remote Control 3.0</i>	Discreet control of volume, programmes and connectivity sources.
TV Adapter		✓ <i>Via Oticon Medical Streamer and TV Adapter 2.0</i>	✓ <i>TV Adapter 3.0</i>	Streams stereo sound directly from a TV.
Apps				
Oticon RemoteCare app		–	✓	Enables remote appointments, adjustments and upgrades in the comfort of the users' own homes – or in other sound environments that they would like optimized.
Control app		ConnectLine app via Oticon Medical Streamer	Oticon ON app	Features in Oticon ON app only: Discreet volume and programme control. Includes a “find my Ponto” location feature and access to HearingFitness™ and the Hearing Guide.
Support and guidance app		Ponto Care app	Ponto Care app	Supports users from trial to aftercare to enable a good experience.
Accessories				
Softband and SoundConnector		✓	✓	Provides access to Ponto for those too young or unable to undergo surgery. Lets users try Ponto before they decide to have the implant.
Safety Line		✓	✓	Prevents loss of sound processor and protects against damage.
Tamper-resistant battery drawer		✓	✓	Prevents accidental opening of the battery door.
Skins and stickers		✓	✓	Enables personalisation.
Fitting				
Wireless fitting		–	✓	Fit Ponto with Noahlink Wireless.
Remote fitting		–	✓ <i>Oticon RemoteCare app</i>	Users can attend appointments online and have programming and adjustments carried out remotely.
NAL-NL 1 BC		✓	✓	Choose the most relevant fitting rationale for the user.
DSL BC		✓	✓	

¹ Data on file at Oticon Medical – Clinical study BC109.

² Data on file at Oticon Medical.

* Increased speech understanding when comparing performance with OpenSound Navigator ON relative to OpenSound Navigator OFF.

** For prescribed fittings, according to best practice and during normal use.

Apple, the Apple logo, iPhone®, iPad®, and iPod touch® are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.