

# Set-up Guide MONO procedure

## Implant components and disposable instruments

		
<p><b>220701</b> MONO Surgery Kit A: Cannula B: MONO drill C: Soft healing cap D: Insertion indicator</p>	<p>Ø 4 mm or Ø 5 mm <b>Biopsy punch</b> Not supplied by Oticon Medical</p>	<p><b>M52169</b> Ponto BHX Implant, 4 mm, with abutment, 9 mm</p>

## Non-disposable instruments

				
<p><b>M50428</b> Sound processor indicator</p>	<p><b>M52143</b> Ruler</p>	<p><b>M52196</b> Double-ended dissector</p>	<p><b>M51695</b> Abutment inserter, machine</p>	<p><b>M51690</b> Counter torque wrench</p>

## Back-up components

	
<p><b>M52170</b> Ponto BHX implant, 4 mm, with abutment, 12 mm</p>	<p><b>M52171</b> Ponto BHX implant, 4 mm, with abutment, 14 mm</p>

Set-up should always include planning for back-up implant components with different abutment lengths to correspond to the soft tissue thickness. Consider whether there is a need for back-up components and instruments for an alternative technique.\*

**Caution**  
The MONO drill must always be used together with the cannula. The cannula provides the stop that prevents drilling deeper than intended.

The Set-up Guide does not replace the Surgical Manual for the Ponto System or the Addendum including the MONO procedure. It is important to read the Surgical Manual and Addendum for a description of patient indications, contraindications and recommended procedures including warnings and precautions. \*For complete product range, refer to Ponto™ System Product Catalogue. Product availability may vary between markets.