

Ponto Streamer

New wireless communication possibilities



Ponto™
– The Bone Anchored
Hearing System



oticon
MEDICAL

Your ideal companion

New communication possibilities

With Ponto Streamer you get access to Oticon ConnectLine and new communication possibilities for your Ponto Plus sound processor. ConnectLine can help you experience better connections – with people, with information and with entertainment.



Cell Phone



TV



Tele Loop

What is ConnectLine?

ConnectLine is a system that helps you connect to a wide variety of communication and entertainment applications – transforming your Ponto Plus sound processor into a personal wireless headset.



Phone



Microphone





Music

Different applications

With ConnectLine you can improve listening and communication possibilities related to phone conversations, one-to-one communication in noisy environments, watching TV, listening to music, video chatting, and more. Start with the application you find most useful for your needs. You can always add more applications at any time.

Talk to your hearing centre today about trying ConnectLine.



Computer



Office Phone



Remote



Ponto Streamer – the heart of ConnectLine

Ponto Streamer is your gateway to a wide range of different electronic devices. Most audio sources can be transmitted through Ponto Streamer wirelessly to your Ponto Plus sound processor. Ponto Streamer is worn around the neck and features three buttons for the core applications: Phone, TV and microphone.

Ponto Streamer also acts as a remote control for your Ponto Plus sound processor and is an easy way to change program or to adjust the volume.

Ponto Streamer:

- Is a single device used to control a wide range of communication possibilities
- Is designed for ease of use
- Turns your sound processor into a headset
- Has a built-in telecoil
- Is an integrated system. Advanced features ensure that, for example, a phone call can be heard and answered whilst watching TV
- Acts as a remote control for your Ponto Plus sound processor





Easy communication

Landline Phone

No more running for the phone and no missed calls.

The ConnectLine phone adapter works with your existing phone (analog/PSTN) and turns your Ponto Plus sound processor into a wireless headset, so you don't have to get up to answer the phone. To take the call simply press a button on the Ponto Streamer.

Video chatting

Video chat* across town or across the world with family and friends!

ConnectLine offers simple and intuitive ways to connect the Ponto Streamer directly to your computer so the sound is sent directly to your sound processor for easy online conversation and direct streaming of sound.

Connect the Ponto Streamer directly to the computer with a cable or go for a wireless connection and you are ready to chat!

**Like Skype and FaceTime on a computer*



Enjoy **entertainment**

ConnectLine opens up entertainment possibilities and connects you to TV, radio, music and more – either alone or with your family.

Watching TV

With ConnectLine you can enjoy TV sound with new freedom. So relax – hear the score, get the news and never miss a punch line.

You can hear the television sound directly in the Ponto Plus sound processor while leaving the volume at a level the whole family can enjoy. Your sound processor acts as your personal headset and the Ponto Streamer as your personal volume control!





The ConnectLine TV adapter transmits sound from the TV, through the Ponto Streamer and into your sound processor.

Listening to music and audio

Ponto Streamer creates a connection to a variety of popular entertainment devices from MP3 players to other audio sources such as computers and audio books. Music and audio are then streamed directly into your sound processor.



Enjoy going out

Hearing the voices you want to hear while in public can be a challenge – on the street, in a restaurant or participating in a service or lecture. More noise often demands “extra help”. ConnectLine offers a range of solutions to help in these situations.

ConnectLine Microphone

The microphone is a small, discreet wireless device that can be worn by the person you’re talking to. It helps you in situations where your companion might be far away, your surroundings might be noisy, or you’re faced with a combination of both.





Cell Phone

The Ponto Streamer helps you make and take calls on your mobile phone – safely and hassle-free. The mobile phone sound is transferred directly to your Ponto Plus sound processor, so it works as a wireless headset.

The Ponto Streamer can connect to a broad range of cell phones with built-in Bluetooth.

Tele loop systems

Enjoy easy access to loop systems in theatres, lecture halls and more. With the built-in telecoil of the Ponto Streamer, a simple push of a button connects you to looped locations.



Improve communication while **driving**

When driving, communicating with a person or talking on the phone can sometimes be a challenge. With ConnectLine we have made it easier for you.

Mobile Phone

The Ponto Streamer transmits cell phone sound wirelessly into the Ponto Plus sound processor, so you will be able to use the cell phone more easily.





Chatting in the car

With the ConnectLine microphone placed on your companion, you can focus on the journey. The microphone transmits the voice of your companion wirelessly to your sound processor, so you can concentrate on the road ahead.



More possibilities **at work**

ConnectLine offers different solutions for you at work. Use the office phone or video conference with colleagues or business partners. In meetings you can also use a discreet clip-on microphone which can send the sound directly to your Ponto Plus sound processor.

Office Phone

The ConnectLine office phone kit enables you to use your sound processor as a wired or wireless headset for office phones.

Computer communication

ConnectLine offers simple and intuitive ways to connect the Ponto Streamer directly to your computer. Audio is transmitted directly to your sound processor for easy online conversation and listening.

About Oticon Medical

Oticon Medical creates innovative hearing implant systems that provide everyday benefits and lifelong positive outcomes for users. Through a commitment to putting user needs first and ongoing investment in product development, research and knowledge sharing, Oticon Medical aims to drive advances in hearing implant systems that benefit users.



Distributor

Oticon Medical AB
Datavägen 37B
SE-436 32 Askim
Sweden
www.oticonmedical.com

Manufacturer

Oticon A/S
Kongebakken 9
2765 Smørum
Denmark
www.oticon.com

United Kingdom

Oticon Medical
Cadzow Industrial Estate
Low Waters Road
Hamilton, ML3 7QE
Phone: +44 (0)1698208234
info@oticonmedical.co.uk

Israel

Medton Ltd.
51 Ben Yehuda Street
Tel Aviv 63431
Israel
Phone: +972 3 620 3992
Info@oticonmedical.co.il

Australia

Oticon Medical/Oticon Australia Pty Ltd
Suite 4, Level 4, Building B
11 Talavera Road
North Ryde NSW 2113
Phone: + 61 (0)2 9857 8288
info@oticonmedical.com.au

South Africa

Oticon Medical
Warich Office Close,
39 Van Vuuren Street
Constantia Kloof 1709, Gauteng
Phone: +27 11 675 6104
info@oticonmedical.co.za

New Zealand

Oticon Medical/Oticon New Zealand Ltd
142 Lambton Quay
Wellington 6021
Phone: +64 (0)4 473 3330
info@oticonmedical.co.nz