

## easyTransport cart

Facilitates baby-friendly, room-to-room hearing screening with your MAICO easyScreen ABR/OAE system.

### Features at a Glance

- The easyTransport cart – a place for everything!
- Options with or without tablet PC
- Secure attachment of the easyScreen cradle and label printer
- "Block-it" bracket keeps easyScreen in the cradle during transport while still permitting easy access for hand-held use when needed
- Screen without removing easyScreen from the cart. Roll up to crib-side; cables are long enough to connect to the baby
- Custom cut-outs for storage of cables
- Storage bins for supplies and documents
- Built-in handle for easy maneuverability
- Charge easyScreen and printer with one hospital grade power strip



Pictured with gooseneck and tablet PC

### Tablet accessory kit for the easyTransport cart

A custom gooseneck holds a tablet PC in a locked enclosure. Keep everything on one compact cart, there is no need to move the system to a different location to transfer your screening data. Complete the entire screening process at crib-side.



Field engineered with direct input from our customers.



## easyTransport cart dimensions

- Measurements include gooseneck and tablet PC
- Weight including power strip = 25 lbs + 4 lbs for the optional tablet accessory kit
- Height = 45.5" from floor to top of shelf + 18" for the optional tablet accessory kit
- Diameter (wheel base) = 21"

## Major components



5-leg wheel base: 21 in diameter



Utility bin: 3-5 in x 12 in x 12 in



Chart bin: 3.25 in x 14.25 in x 9 in



Mounting bracket with hospital-grade coiled cord power strip with medical 4 outlets

If you rejected a hand-held newborn hearing screening instrument before because you preferred the convenience of a cart-based system, we invite you to reconsider. Take a look at MAICO easyScreen with the easyTransport cart.

Contact us for a demonstration or for more information.



### MAICO Diagnostics

10393 West 70th Street  
Eden Prairie, MN 55344, USA

Tel.: (888) 941-4201  
Fax: (952) 903-4100

info@maico-diagnostics.com  
www.maico-diagnostics.com/us

**ONE YEAR LIMITED WARRANTY:** This warranty is extended to the original purchaser of the instrument, by MAICO, through the Distributor from whom it was purchased. The warranty covers defects in material and workmanship for a period of one year from date of invoice of the instrument to the original purchaser. Accessories which are purchased from MAICO at the same time as the instrument are also under warranty for one year from the date of invoice. For additional information contact your MAICO Distributor.