

# MA 1

## Air Conduction Audiometer

Smallest handheld portable air-conduction audiometer.

### Features At a Glance

- Portable air conduction screener (500Hz, 1K, 2k, 4KHz)
- Battery operated (2AA alkaline)
- Easy one hand push button operation
- Exceptionally lightweight, only 1 pound
- Carrying case included

### Sleek, Small Design

With its super-lightweight design, the MA 1 is small enough to hold in the palm of your hand. Tucked into its carrying case (included) the MA 1 is our smallest, lightest portable audiometer. Although it is small, it offers a full range of screening.

### Accurate, Dependable Engineering

Because it's portable, the MA 1 is engineered to be tough and reliable. Take the MA 1 anywhere with confidence. Its user-friendly push button design makes adjusting frequency and intensity easy.



## Technical Data

Test Frequencies	500, 1000, 2000, 4000 Hz
Rise/Fall Time	35 msec typical
Attenuator Range	15 to 50 db <sub>HL</sub> 5 db steps
Distortion	0.5% THD typical, 3% THD maximum
Stimulus	Continuous
Dimensions	2.5 in x 6 in x 7/8 in
Weight	2.1 lbs
Power Supply	Battery operated (2 AA alkaline)

### STANDARD ACCESSORIES

- DD45 headphones with RE-7 headband
- Carrying Case
- Operation Manual

Specifications subject to change without notice.

Standard  
Components







Carrying Case

DD45 headset with  
RE-7 headband

MA 1 Handheld  
Device